

GB	Continuous peripheral nerve block catheter set	5
DE	Katheterbesteck für die kontinuierliche Nervenblockade	6
BG	Контиплекс Туохи – Катетърен комплект за продължителна периферна нервна блокада	7
CZ	Souprava katétrů pro nepřetržitou blokádu periferních nervů	8
DK	Katetersæt til kontinuerlig ledningsanalgesi	9
EE	Püsiva perifeerse närviblokaadi kateetri komplekt	10
ES	Set de catéter para el bloqueo nervioso continuo	11
FI	Katetrisarja ääreishermoston jatkuvaan salpaukseen	12
FR	Set de cathétérisme pour bloc nerveux périphérique continu	13
GR	Σετ καθετήρα για το συνεχή αποκλεισμό νεύρων	14
HR	Kateter set za kontinuirane periferne živčane blokove	15
HU	Katéterkészlet folyamatos perifériás idegblokádhoz	16
IT	Set catetere per blocco nervoso periferico continuo	17
LT	Kateterio rinkinys nuolatinei periferinių nervų blokada	18
LV	Ilgstoša perifērā nervu bloka katetra komplekts	19
NL	Kateterset voor continue perifere zenuwblokkering	20
PL	Zestaw do ciągłej blokady nerwu obwodowego z cewnikiem	21
PT	Cânula para bloqueio de plexo contínuo	22
RO	Set cu cateter pentru anestezie plexală continuă	23
SE	Kateterset för kontinuerlig perifer nervblockad	24
SI	Katetrski set za neprekinjeno blokado živcev	25
SK	Súprava katétra na kontinuálnu blokádu periférneho nervu	26







GB Line drawings of the individual components
DE Strichzeichnungen der Komponenten
BG Скици на отделните съставни части
CZ Perokresby jednotlivých částí
DK Stregtegninger af de enkelte dele
EE Üksikute komponentide skeemid
ES Esquemas de los componentes individuales
FI Osien merkit
FR Croquis des différents constituants
GR Γραμμικά σχέδια των ξεχωριστών στοιχείων

HR Prikaz sastavnih dijelova
HU Az egyes alkotórészek vonalas rajzai
IT Schemi dei singoli componenti
LT Kiekviena detalė pažymėta brūkšniu kodu
LV Atsevišķo materiāla sastāvdaļu svītrkodsī
NL Tekeningen van de verschillende onderdelen
PL Rysunki pojedynczych części składowych
PT Traços da marca de cada componente
RO Desene liniare ale componentelor
SE Bilder på de enskilda komponenterna

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




GB Instructions for use

Contents

1. Insulated Contiplex® S needle with integrated wire (Fig. A)
2. Hemostasis Valve/Sideport assembly with 20 in. (51 cm) tubing (Fig. B)
3. 20 Ga Perifix® catheter with threading assist guide (Fig. C)
4. Perifix® catheter connector (Fig. D)
5. Perifix® filter (bacteria filter) (Fig. E)
6. Perifix® PinPad (Fig. F)
7. Omnifix® 5ml syringe (Fig. G)
8. Catheter label (Fig. H)

Materials used

ABS, beryllium, brass, copper, IR, PA, PE, PET, PMMA, POM, PP, PVC, SAN, stainless steel, TPE, medical quality foam with synthetic polyisoprene rubber, tin

Important!

Contiplex® S Sets with NRFit® components are not compatible with Luer devices. They are only suitable for products with NRFit® connector!

Indications

Plexus and peripheral nerve blocks using electrical nerve stimulation as block anesthesia for surgical procedures involving the upper and lower extremities with various block techniques, as well as pain treatment, analgesia for early mobilization and physiotherapy or acute and long-term analgesia (e.g. cancer diseases).

Contraindications

- Pre-existing neurological diseases
- Infection/inflammation at the side of insertion
- Excessive preoperative anxiety (not responsive to interventions) and inability to tolerate positioning
- Pre-existing coagulopathies, either endogenous (e.g. thrombocytopenia) or iatrogenic (e.g. warfarin therapy) due to increased risk of hematoma
- Anatomic anomalies that make identification of physical landmarks difficult
- Hepatic disease that may interfere with clearance
- Failure to obtain patient's consent
- Known hypersensitivity to one or more materials used
- For information on other contraindications against providing continuous regional anesthesia, please refer to the specialist literature

NOTE:

The product should only be used by physicians that have been adequately trained in this technique.

Risks

- Nervous lesions
- Pneumothorax
- Pain and haematoma at the site of insertion of the needle
- Incomplete block
- Local anaesthetics related toxicity: cardiac, neurologic, allergic
- Increased risk of thrombosis and risk of thromboembolism
- Paravertebral technique:
 - epidural/intrathecal injection

Warning

- Products with NRFit® connector are not compatible with Luer devices. They are only suitable for products with NRFit® connector!
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to

contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

- In order to achieve the optimum connection between the stimulation cannula and the nerve stimulator, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.
- Do not resterilize.
- Use only if the packaging is intact.
- Use the properly stored product by the date stated on the packaging.

DIRECTIONS: Use Aseptic Technique.

1. Attach sideport to hub of insulated Contiplex® S needle.
2. Attach syringe filled with Sodium Chloride or anesthetic solution to proximal hub of sideport tubing. Prime tubing and needle with solution.
3. If a peripheral nerve stimulator is used, attach electrode pin to the conducting wire.
4. Disinfect thoroughly twice (allow time to act!) the relevant region and cover with a sterile non-fenestrated and/or fenestrated drape.
5. After raising skin wheal with anesthetic agent, introduce needle through the puncture site towards targeted neurovascular bundle at an angle according to the used insertion technique.
6. Advance the needle in the direction of the nerve until visible muscle contractions occur in the innervated area.
 - **NOTE:** Contiplex® S cannulas have length markings (every 1 cm). Additional markings: 5 cm – 1 broad ring mark / 10 cm – 2 broad ring marks
7. Reduce the current and optimize needle position until muscle contractions occur at lower current levels. The tip of the needle has reached an optimal position when noticeable contractions occur at a current of approximately 0.2-0.5 mA (higher current levels may be required for certain lumbar blocks). Aspirate for possible intravascular placement.
 - **CAUTION:** If paresthesia is caused inadvertently as a result of direct contact with the nerve, the stimulation cannula must not, under no circumstance, be inserted any further.
8. Following negative aspiration, a test dose of local anesthetic can be injected through the extension tubing. Muscle contractions should cease within 5 to 10 seconds. Prior to indwelling catheter placement, the desired dose of anesthetic agent can be injected through the sideport tubing.
 - **CAUTION:** Once the local anesthetic has been injected, the needle cannot be repositioned as it is usually no longer possible to get the nerves stimulated and the patient no longer experiences pain in the event of the cannula's inadvertent intraneural placement.
9. The plastic threading assist guide should be attached to the hemostasis valve of the sideport. The catheter should then be inserted through the threading assist guide. It will continue through the valve and through the insulated cannula.
10. Insert catheter to desired depth. Hold the catheter in place and slowly remove the insulated needle/sideport device.
 - **NOTE:** Do not use any adhesive spray which contains organic solvents to fix the catheter.

11. Check that the connector hinged cap is in the open position. Insert proximal end of catheter into opening on connector. (Figure 1)
 - **NOTE:** Inserting the catheter when hinged cap is not in the open position may cause the catheter to not seat properly.
12. Ensure that the catheter is fully inserted and seated to the hinged area of the connector. (see dotted line on Figure 2)
 - **NOTE:** Catheter may occlude or leak if not seated all the way to the end. If this should happen, please verify the catheter is seated all the way to the end of the connector.
13. Close the top hinged cap of the connector body by depressing it until an audible "click" is heard. If not securely attached, re-open connector and repeat steps 11 and 12. (Figure 3)
14. After removal of the protective cap attach the connector to filter or anesthetic delivery device with a twist action (NRFit® connector). (Figure 4)
 - **CAUTION:** After initial connection to the catheter, do not attach any anesthetic fluid administration devices to the connector if the replacement sterility cap is missing. You must assume that the connector has been contaminated and it should be replaced.
 - **CAUTION:** During connection to a Perifix® catheter connector, the catheter may not be fully inserted into the connector resulting in the potential for leakage of anesthetic agent intended for delivery to the patient via the catheter. Loosen the top hinged piece of the connector and completely re-insert the catheter until it cannot be advanced any further. Then engage the connector by depressing the top hinged piece.
 - **CAUTION:** During manipulation of the connector to attach syringes, pump or IV tubing, filters, or extension set tubing, the connector may be inadvertently unlocked. This action could permit the catheter to come loose from the connector risking contamination of the proximal end of the catheter, and/or interruption of anesthetic delivery.
15. Aspirate and administer additional anesthetic as needed.
 - **CAUTION:** If a Perifix® filter (bacteria filter) is being used, do not use syringes smaller than 10 ml as they may cause the filter housing to be damaged as a result of possibly excessive pressure from the syringe.
16. To re-open the connector, firmly hold the connector with two fingers on both ends and press on both ends firmly flexing the connector as if you are trying to snap the connector in half. The lid will spring open with an audible click. (Figure 5)
 - **CAUTION:** Do not withdraw catheter through needle due to possible danger of shearing.

CAUTION:

Do not use any adhesive spray which contains organic solvents to fix the catheter.

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