Contiplex® Tuohy Ultra



GB	Continuous peripheral nerve block catheter set	5
DE	Katheterbesteck für die kontinuierliche Nervenblockade	6
BG	Контиплекс Туохи – Катетърен комплект за продължителна периферна нервна блокада	7
CN	连续外周神经阻滞导管套件	8
CZ	Souprava katétrů pro nepřetržitou blokádu periferních nervů	9
DK	Katetersæt til kontinuerlig ledningsanalgesi	10
EE	Püsiva perifeerse närviblokaadi kateetri komplekt	11
ES	Set de catéter para el bloqueo nervioso continuo	12
FI	Katetrisarja ääreishermoston jatkuvaan salpaukseen	13
FR	Set de cathétérisme pour bloc nerveux périphérique continu	14
GR	Σετ καθετήρα για το συνεχή αποκλεισμό νεύρων	15
HR	Kateter set za kontinuirane periferne živčane blokove	16
HU	Katéterkészlet folyamatos perifériás idegblokádhoz	17
ID	Set kateter untuk blok saraf perifer secara terus-menerus	18
IT	Set catetere per blocco nervoso periferico continuo	19
LT	Kateterio rinkinys nuolatinei periferinių nervų blokada	20
LV	Ilgstoša perifērā nervu bloka katetra komplekts	21
NL	Kateterset voor continue perifere zenuwblokkering	22
NO	Katetersett til kontinuerlig perifer nerveblokkade	23
PL	Zestaw do ciągłej blokady nerwu obwodowego z cewnikiem	24
PT	Cânula para bloqueio de plexo contínuo	25
RO	Set cu cateter pentru anestezie plexală continuă	26
RS	Set katetera za kontinuiranu blokadu perifernog nerva	27
RU	Контиплекс Туохи – Набор для продленной блокады нервных стволов и сплетений	28
SE	Kateterset för kontinuerlig perifer nervblockad	29
SI	Katetrski set za neprekinjeno blokado živcev	30
SK	Súprava katétra na kontinuálnu blokádu periférneho nervu	31
TR	Kontinü periferik sinir blokaj kateter seti	32

B | BRAUN

GB Line drawings of the individual components
DE Strichzeichnungen der Komponenten
BG Скици на отделните съставни части
CN 单独组件的图示
CZ Perokresby jednotlivých částí
DK Stregtegninger af de enkelte dele
EE Üksikute komponentide skeemid
ES Esquemas de los componentes individuales
FI Osien merkit
FR Croquis des différents constituants







GR Γραμμικά σχέδια των ξεχωριστών στοιχείων
HR Prikaz sastavnih dijelova
HU Az egyes alkotórészek vonalas rajzai
ID Ilustrasi setiap komponen
IT Schemi dei singoli componenti
LT Kiekviena detalė pažymėta brūkšniu kodu
LV Atsevišķo materiāla sastāvdaļu svītrkodi
NL Tekeningen van de verschillende onderdelen
NO Strektegninger av de enkelte bestanddelene
PL Rysunki pojedynczych części składowych

PT Traços da marca de cada componente
RO Desene liniare ale componentelor
RS Skice pojedinačnih komponenti
RU Рисунки отдельных компонентов
SE Bilder på de enskilda komponenterna
SI Shematski prikazi posameznih sestavnih delov
SK Perokresby jednotlivých částí
TR Komponentlerin ayrı çizimleri

Kanülenlänge Needle length	Abstand Spitze -> Stop Block Distance Catheter tip -> Stop Block
Contiplex Tuohy	Measure*
40 mm	94 mm
50 mm	106 mm
100 mm	156 mm
150 mm	206 mm






Catheter-	Measure**
0,84 mm	29,5 ± 0,5 mm

GB Directions
DE Anwendungshinweise
BG Указания за употреба
CN 说明
CZ Pokyny pro použití
DK Betjeningsvejledning
EE Kasutusjuhendid
ES Aplicaciones
FI Käyttöohje
FR Conseils d'utilisation

GR Υποδείξεις εφαρμογής
HR Upute
HU Használati útmutatás
ID Petunjuk
IT Modalità d'uso
LT Įspėjimai
LV Norādījumi par lietošanu
NL Gebruiksaanwijzingen
NO Bruk
PL Przepisy stosowania

PT Instruções de uso
RO Instrucţiuni de folosire
RS Uputstva
RU Указания по применению
SE Användningsanvisningar
SI Napotki za uporabo
SK Aplikac né pokyny
TR Kullanım talimatları

GB	Do not re-use
DE	Nicht wiederverwenden
BG	Да не се използва повторно
CN	禁止重复使用
CZ	Nepoužívat opětovně
DK	Må ikke genbruges
EE	Mitte kordukasutada
ES	No reutilizable
FI	Ei saa käyttää uudelleen
FR	Ne pas réutiliser
GR	Να μην επαναχρησιμοποιείται
HR	Nemojte ponovno potrijebiti
HU	Ne használja fel újra
ID	Jangan gunakan kembali
IT	Non riutilizzare
LT	Negalima naudoti kartotinai
LV	Neizmantot atkārtoti
NL	Niet hergebruiken
NO	Kun engangsbruk
PL	Nie używać ponownie
PT	Não reutilizar
RO	A nu se reutiliza
RS	Ne koristiti ponovo
RU	Запрет на повторное применение
SE	Får ej återanvändas
SI	Samo za enkratno uporabo
SK	Opätovne nepoužívať
TR	Yeniden kullanmayınız

Consult instruction for use
Gebrauchsanweisung beachten
Вижте инструкциите за употреба
阅读使用说明
Čtěte návod k použití
Læs brugsanvisningen
Vt kasutusjuhendit
Consulte las instrucciones de uso
Katso käyttöohje
Consulter le mode d'emploi
Συμβουλευτείτε τις οδηγίες χρήσης
Pogledajte upute za uporabu
Nézze meg a használati utasítást
Baca petunjuk penggunaan
Consultare le istruzioni per l'uso
Žr. naudojimo instrukcijas
Lasīt lietošanas instrukciju
Raadpleeg gebruikersinformatie
Se i bruksanvisningen
Przed użyciem należy zapoznać się z instrukcją
Consultar instruções para uso
Consultați instrucțiunile de utilizare
Pogledajte uputstvo za upotrebu
Обратитесь к инструкции по применению
Se bruksanvisningen
Glejte navodila za uporabo
Pozri návod na použitie
Kullanma talimatına bakınız

Batch number
Chargennummer
Партиден номер
批号
Kód dávky
Lot-nr.
Parti number
Código de lote
Eränumero
Numéro de lot
Αριθμός παρτίδας
Broj serije
Gyártási sorozat száma
Nomor batch
Numero di lotto
Partijos numeris
Sērijas numurs
Chargennummer
Batch/LOT-nummer
Numer serii
Número do lote
Lot nr.
Broj partije
Номер серии
Batchnummer
Številka serije
Číslo šarže
Parti kodu

Caution
Achtung
ВНИМАНИЕ
警告
Pozor (výstraha)
Forsigtig!
Ettevaatus
Precaución
Huomio
Attention
Προσοχή
Oprez
Figyelem!
Perhatian
Attenzione
Atsargiai!
Uzmanību
LET OP!
Forsiktig
Uwaga
Cuidado
Avertizare
Oprez
Осторожно! Обратитесь к инструкции по применению
Varning!
Previdno
Varovanie
Dikkat

STERILE

GB	Sterile
DE	Steril
BG	Стерилен
CN	无菌
CZ	Sterilní
DK	Steril
EE	Steriilne
ES	Estéril
FI	Steriili
FR	Stérile
GR	Αποστειρωμένο
HR	Sterilno
HU	Steril
ID	Steril
IT	Sterile
LT	Sterilu
LV	Sterils
NL	Steriel
NO	Steril
PL	Sterylny
PT	Esterilizado
RO	Steril
RS	Sterilno
RU	Стерильно
SE	Steril
SI	Sterilno
SK	Sterilné
TR	Steril

Use-by date
Verwendbar bis
Срок на годност
有效期
Použit do data
Anvendes inden
Kasutada kuni:
Fecha de caducidad
Viimeinen käyttöpäivä
Date limite d'utilisation
Ανόλωση έως
Upotrijebiti do
Szavatossági idő
Gunakan sebelum
Da utilizarsi entro
Tinka iki datos
Derīguma termiņš
Te gebruiken tot
Holdbarhetsdato
Data przydatności do użycia
Prazo de validade
Data de expirare
Rok trajanja
Использовать до
Används före
Rok uporabnosti
Použitelné do
Son kullanım tarihi

Date of manufacture
Herstellungsdatum
Дата на производство
生产日期
Datum výroby
Fremstillingsdato
Tootmiskuupäev
Fecha de fabricación
Valmistuspäivä
Date de fabrication
Ημερομηνία κατασκευής
Datum proizvodnje
Gyártási dátum
Tanggal produksi
Data di produzione
Pagaminimo data
Izgatavošanas datums
Productiedatum
Produksjonsdato
Data produkcji
Data de fabrico
Data fabricației
Datum proizvodnje
Дата изготовления
Tillverkningsdatum
Datum izdelave
Dátum výroby
İmalat tarihi

DEHP-free
DEHP-frei
Не съдържа DEHP
不含 DEHP
Neobsahuje DEHP
DEHP-fri
DEHP-vaba
No contiene DEHP
Ei sisällä DEHP:a
Sans DEHP
Δεν περιέχει DEHP
Bez DEHP-a
DEHP-mentes
Bebas DEHP
Privo di DEHP
Be DEHP
Nesatur DEHP
DEHP-vrij
DEHP-fri
Nie zawiera DEHP
Sem DEHP
fără DEHP
Ne sadrži DEHP
Не содержит фталаты
DEHP-fri
Brez DEHP-ja
Neobsahuje DEHP
DEHP iğermez

GB Instructions for use

Contents

1. Insulated Tuohy needle with integrated wire (Fig. A)
2. Hemostasis Valve/Sideport assembly with 20 in. (51 cm) tubing (Fig. B)
3. 20 Ga Perifix® ONE catheter with threading assist guide (Fig. C)
4. Perifix® catheter connector (Fig. D)
5. Perifix® filter (bacteria filter) (Fig. E)
6. Perifix® PinPad (Fig. F)
7. Omnifix® 5ml syringe (Fig. G)
8. Catheter label (Fig. H)

Materials used

ABS, beryllium, brass, copper, HDPE, IR, PA, PE, PET, POM, PMMA, PP, PUR, PVC, SAN, stainless steel, TPE, medical quality foam with synthetic polyisoprene rubber

Indications

Plexus and peripheral nerve blocks using percutaneous electrical nerve stimulation (PENS) as block anesthesia for surgical procedures involving the upper and lower extremities with various block techniques, as well as pain treatment, analgesia for early mobilization and physiotherapy or acute and long-term analgesia (e.g. cancer diseases).

Contraindications

- Pre-existing neurological diseases
- Infection/inflammation at the side of insertion
- Excessive preoperative anxiety (not responsive to interventions) and inability to tolerate positioning
- Pre-existing coagulopathies, either endogenous (e.g. thrombocytopenia) or iatrogenic (e.g. warfarin therapy) due to increased risk of hematoma
- Anatomic anomalies that make identification of physical landmarks difficult
- Hepatic disease that may interfere with clearance
- Failure to obtain patient's consent
- Known hypersensitivity to one or more materials used
- For information on other contraindications against providing continuous regional anesthesia, please refer to the specialist literature

NOTE:

The product should only be used by physicians that have been adequately trained in this technique.

Risks

- Nervous lesions
- Pneumothorax
- Pain and haematoma at the site of insertion of the needle
- Incomplete block
- Local anaesthetics related toxicity: cardiac, neurologic, allergic
- Increased risk of thrombosis and risk of thromboembolism
- Paravertebral technique: - epidural/intrathecal injection

Warning

- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- In order to achieve the optimum connection between the stimulation cannula and the nerve stimulator, thereby ensuring patient

safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.

- Do not resterilize.
- Use only if the packaging is intact.
- Use the properly stored product by the date stated on the packaging.

DIRECTIONS: Use Aseptic Technique.

1. Attach sideport to hub of insulated Tuohy needle.
2. Attach syringe filled with Sodium Chloride or anesthetic solution to proximal hub of sideport tubing. Prime tubing and needle with solution.
3. If a peripheral nerve stimulator is used, attach electrode pin to the conducting wire.
4. Disinfect thoroughly twice (allow time to act!) the relevant region and cover with a sterile non-fenestrated and/or fenestrated drape.
5. After raising skin wheal with anesthetic agent, introduce needle through the puncture site towards targeted neurovascular bundle at an angle according to the used insertion technique.
6. Advance the needle in the direction of the nerve until visible muscle contractions occur in the innervated area.
NOTE: Contiplex® Tuohy cannulas have length markings (every 1 cm). Additional markings: 0,5 cm – 1 broad ring mark / 10 cm – 2 broad ring marks
7. Reduce the current and optimize needle position until muscle contractions occur at lower current levels. The tip of the needle has reached an optimal position when noticeable contractions occur at a current of approximately 0.2–0.5 mA (higher current levels may be required for certain lumbar blocks). Aspirate for possible intravascular placement.
CAUTION: If paresthesia is caused inadvertently as a result of direct contact with the nerve, the stimulation cannula must not, under no circumstance, be inserted any further.
8. Following negative aspiration, a test dose of local anesthetic can be injected through the extension tubing. Muscle contractions should cease within 5 to 10 seconds. Prior to indwelling catheter placement, the desired dose of anesthetic agent can be injected through the sideport tubing.
CAUTION: Once the local anesthetic has been injected, the needle cannot be repositioned as it is usually no longer possible to get the nerves stimulated and the patient no longer experiences pain in the event of the cannula's inadvertent intraneural placement.
9. The plastic threading assist guide should be attached to the hemostasis valve of the sideport. The catheter should then be inserted through the threading assist guide. It will continue through the valve and through the insulated cannula.
10. Insert catheter to desired depth. Hold the catheter in place and slowly remove the insulated needle/sideport device.
NOTE: Do not use any adhesive spray which contains organic solvents to fix the catheter.
11. Check that the connector hinged cap is in the open position. Insert proximal end of catheter into opening on connector. (Figure 1)
NOTE: Inserting the catheter when hinged cap is not in the open position may cause the catheter to not seat properly.
12. Ensure that the catheter is fully inserted and seated to the hinged area of the connector. (see dotted line on Figure 2)

NOTE: Catheter may occlude or leak if not seated all the way to the end. If this should happen, please verify the catheter is seated all the way to the end of the connector.

13. Close the top hinged cap of the connector body by depressing it until an audible "click" is heard. If not securely attached, re-open connector and repeat steps 11 and 12. (Figure 3)
14. After removal of the protective cap attach the connector to filter or anesthetic delivery device with a twist action (Luer Lock). (Figure 4)

CAUTION: After initial connection to the catheter, do not attach any anesthetic fluid administration devices to the connector if the replacement sterility cap is missing. You must assume that the connector has been contaminated and it should be replaced.

CAUTION: During connection to a Perifix® catheter connector, the catheter may not be fully inserted into the connector resulting in the potential for leakage of anesthetic agent intended for delivery to the patient via the catheter. Loosen the top hinged piece of the connector and completely re-insert the catheter until it cannot be advanced any further. Then engage the connector by depressing the top hinged piece.

CAUTION: During manipulation of the connector to attach syringes, pump or IV tubing, filters, or extension set tubing, the connector may be inadvertently unlocked. This action could permit the catheter to come loose from the connector risking contamination of the proximal end of the catheter, and/or interruption of anesthetic delivery.

15. Aspirate and administer additional anesthetic as needed.

CAUTION: If a Perifix® filter (bacteria filter) is being used, do not use syringes smaller than 10 ml as they may cause the filter housing to be damaged as a result of possibly excessive pressure from the syringe.

16. To re-open the connector, firmly hold the connector with two fingers on both ends and press on both ends firmly flexing the connector as if you are trying to snap the connector in half. The lid will spring open with an audible click. (Figure 5)

CAUTION: Do not withdraw catheter through needle due to possible danger of shearing.

CAUTION:

Do not use any adhesive spray which contains organic solvents to fix the catheter.

When using ultrasound, reflections of the cannula tip geometry and the characteristic shape of the distal cannula shaft should be visible and should indicate the cannula position relative to the nerve. Please be aware that the visibility is influenced by a number of factors, e.g. the correct setting of the ultrasound device, patient-related tissue disorders, artefacts as well as the correct alignment of the ultrasound probe and the cannula (recommended angle: 45°).

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B.Braun Melsungen AG, Carl-Braun Str. 1,
Distributed by: B. Braun Medical Inc.