Espocan® CSE Set

NRFit®

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









For further technical characteristics such as needle and catheter size, see product packaging.

Epidural needles, spinal needles, catheters and other medical devices or accessories can be used with devices which are compatible or in conformance with NRRFit® connectors (e.g. ISO 80369-6).

Espocan® Sets with NRRFit® components are not compatible with Luer devices. They are only suitable for products with NRRFit® connector!

External devices can be attached to the female connector of the Perifix® Filter, which is the point of access to the applied Perifix® Catheter. Main applications will be syringes for aspirating and infusing anesthetics or connecting elastomeric infusion pumps respectively extension lines of syringes intended for infusion pumps for continuous infusion of anesthetics or analgesics.

Assessment of Magnetic Resonance (MR) compatibility

Espocan® epidural needle NRRFit® / Espocan® Sleeved spinal needle NRRFit® / Perifix® LOR syringe NRRFit®: The device has not been evaluated for safety in the MR environment, since it is not expected to enter the MR environment in a typical application. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

Perifix® Catheter Connector NRRFit® / Perifix® Filter 0.2 µm NRRFit® / Perifix® PinPad / Docking System NRRFit®: MR safe.

Perifix® Catheter NRRFit® / Perifix® SoftTip Catheter NRRFit®: MR conditional (please check Operating Instructions section for additional conditions).

Sterility

- Sterilized using ethylene oxide.
- Do not use if package is damaged.
- Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Set for combined spinal and epidural anesthesia.

Indication

Used for combination of spinal and epidural anesthesia during surgery and epidural analgesia for regional pain treatment. If necessary, anesthesia can be prolonged during surgery. Analgesia is possible at any time after surgery (pain treatment).

Patient population

Adults, pregnant women

Intended User

The products should only be used by healthcare workers that have been adequately trained in this technique.

Contraindications

Known contraindications are the same as those for spinal and epidural anesthesia and include, in particular, the following absolute and relative contraindications:

Absolute: patient refusal, sepsis, local infections at the injection site, immune deficiency, severe decompensated hypovolemia, shock, cardiovascular diseases of myocardial, ischemic, or valvular origin, acute cerebral or spinal cord diseases, increased intracranial pressure, indeterminate neurological disease, hypersensitivity to local anesthetic agents or any of the materials employed.

Relative: coagulopathy (coagulation disorders or anticoagulant therapy), spinal pathology (including severe spinal deformities, arthritis, osteoporosis, intervertebral disc prolapse, spinal canal stenosis, post-spinal surgery and spinal metastases), unknown duration of surgery, mild vascular disease (aortic stenosis), preexisting neurological deficits (radiculopathies, peripheral neuropathies, multiple sclerosis).

Taken from: Jankovic (Jankovic D. Combined spinal-epidural anaesthesia. In: Jankovic D, Peng PWH. Regional nerve blocks in anesthesia and pain therapy. Cham [u.a.]: Springer. 2015; 575-583).

Please also consult local medical textbooks for epidural and spinal anesthesia.

Residual Risks/Side effects

Known risks of spinal and epidural anesthesia techniques include, in particular, hypotension, dyspnea / apnea, cardiac arrest, bradycardia, urinary retention, toxic effects of the local anesthetic / adjuvant (please refer to the Summary of Product Characteristics (SPC) of the drug), temporary or permanent neurological complications and local or systemic infections (post-dural puncture headache (PDPH), paresthesia, meningitis / cauda equine syndrome, transverse myelitis or anterior spinal artery syndrome with permanent paralysis, epidural hematoma / epidural space infections, abscesses, skin and subcutaneous tissue infections), inadvertent puncture of blood vessels.

Technical complications with epidural catheters:

- Catheter dislodgement, migration, bent or knotted catheters, which may lead to block failures or occlusion
- Shearing off of catheter fragments due to trauma during removal / insertion

Precautions

- Trained staff and sufficient emergency equipment, in particular for artificial ventilation and hemodynamic support, must be nearby
- Monitor all connections of the catheter set-up for tightness and disconnection during whole anesthesia and analgesia procedure (e.g. pain treatment)
- Use strictly aseptic techniques
- Do not use if package is damaged
- Do not use if product is incomplete or damaged
- Test dose to check the position of the epidural catheter

- Visually inspect the sterile barrier system for breaches and integrity prior to use

Warning

- Single use only. Do not re-use. Do not re-sterilize.
- Espocan® NRRFit® components have standardized connectors acc. to ISO 80369-6 and are not compatible with Luer devices. They are only suitable for products with NRRFit® connector!
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Never pull the catheter through the needle as it may shear off otherwise!
- The tip of the needle becomes blunt or can even split after bone contact. The needle might bend and eventually break when applying too strong force during advancement. Needles with damaged tip increase the risk of intrathecal and intravascular placement.
- Do not wet Perifix® Filter NRRFit® 0.2 µm with disinfectants containing alcohol as this may result in cracks in the housing.
- Do not use syringes smaller than 10 ml as they may cause the filter housing to be damaged as a result of excessive pressure from the syringe.
- The catheter must be inserted into the connector as far as possible. Injection is impossible otherwise and the catheter can become disconnected.
- Do not use an adhesive spray which contains organic solvents to fix the catheter.
- Inspect the insertion site of the catheter daily. Remove the catheter in the case of any signs of local or systemic inflammation of unknown origin and according to hospital protocols.
- If paresthesia occurs, do not insert the spinal needle further. If necessary, pull it back until there is no paresthesia.
- Only slightly fasten the docking system as otherwise the spinal needle cannot be inserted and moved.
- Always use Espocan® Sets with Docking System together with the docking system (Fig. J) as otherwise the spinal needle is too long (tip protrudes) and the dura on the ventral side may be punctured unintentionally.
- In order to achieve the optimum connection between epidural needle, spinal needle, catheter and other medical devices or accessories, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.
- Perifix® Catheter NRRFit® / Perifix® SoftTip Catheter NRRFit®: MR conditional (please check Operating Instructions section for additional conditions).

Perifix® LOR Syringe NRRFit®:

Special syringe for the Loss of Resistance (LOR) technique for regional anesthesia. Syringe, only suitable to a limited extent for aspiration of liquids. Do not aspirate beyond the mechanical stop. For hygiene reasons, do not fill the syringe until shortly before use.

Operating Instructions

- Illustrations on page 2/3.

► Directions for using Espocan® sets NRFit®

1. The Espocan® epidural needle NRFit® (Fig. C) and the Perifix® LOR syringe NRFit® (Fig. B) are used to puncture the epidural space using the usual Loss of Resistance (LOR) technique. The chosen intervertebral space must also be suitable for subsequent spinal anesthesia.
2. The epidural space is identified, and the Espocan® Sleeved spinal needle NRFit® (Fig. A) is inserted into the epidural needle. A mark on the spinal needle (Fig. 2) indicates where the spinal needle protrudes from the back-eye of the epidural needle. The subsequent puncture of the dura can usually be felt as a significant change in resistance.
3. Remove the stylet of the spinal needle slowly and evenly. When the spinal needle is in the correct position, there will be a free flow of cerebrospinal fluid (CSF). In this case spinal anesthesia can be performed. If there is not a free flow of CSF, the spinal needle should be pulled back under constant observation. If pulling the needle back does not produce a free flow of CSF, remove the spinal needle completely and carefully correct the position of the epidural needle or make a new puncture.
4. After spinal anesthesia remove the spinal needle. Use the threading assist guide (Fig. E) (Fig. 1) and insert the Perifix® Catheter NRFit® (Fig. D) through the epidural needle into the epidural space and up to the required position. In this case, as in the case of paresthesia, observe whether there is a free flow of CSF (intrathecal position) or a backflow of blood (intravascular position). The catheter is marked on the tip and starting at 55 mm (59 mm for SoftTip) in distances of 10 mm. The 100, 150, and 200 mm positions are marked with 2, 3, and 4 bands respectively. The position when the catheter tip protrudes from the epidural needle is marked with a broad band (Fig. D*). The catheter openings (holes) are located within the first 40 mm of the catheter tip.
5. Pull the epidural needle back over the catheter completely, then insert the catheter as far as possible into the Perifix® Catheter Connector NRFit® (Fig. F) and close (Fig. 4–7). The Catheter Connector works according to the Click & Ready principle.
 - Open the catheter connector lid as far as possible (Fig. 4).
 - Insert the catheter into the funnel-shaped opening on the connector tip (Fig. 5).
 - Push the catheter well into the connector until it clearly reaches a stop (Fig. 5).
 - Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 6).
 - Remove air from the Perifix® Filter NRFit® 0.2 µm (Fig. G, if provided) using physiological saline solution.
 - Remove closing cone from connector, attach the connector to the Perifix® Filter NRFit® with a twist action (NRFit® Lock) (Fig. 7).
6. Rinse the system with 1–2 ml physiological saline solution. If necessary: Re-opening of the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 8).
7. Secure the filter to the patient using the B. Braun Perifix® PinPad (Fig. H).

Epidural catheter:

A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
 - Maximum spatial field gradient of 2,500 G/cm (25 T/m)
 - Maximum force product of 50,000,000 G²/cm (50 T²/m)
 - Maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) without time limit in continuous scanning.

Test dose to check the position of the catheter

A test dose cannot be administered immediately if spinal anesthesia has been performed beforehand. If the effects of epidural anesthesia need to be extended or prolonged during surgery, small quantities of the local anesthetic can be administered at intervals of 3 to 5 minutes until the desired level of anesthesia is reached. During the action of spinal anesthesia, the effect of local anesthetics administered through the epidural needle is greater and the time to act faster than with epidural anesthesia alone. If surgery was finished without a subsequent epidural injection, allow the effect of spinal anesthesia to decline before administering the usual test dose.

► Directions for use of a set with the Espocan® Docking System NRFit®

1. The Espocan® epidural needle NRFit® (Fig. C) and the Perifix® LOR syringe NRFit® (Fig. B) are used to puncture the epidural space using the usual Loss of Resistance (LOR) technique. The chosen intervertebral space must also be suitable for subsequent spinal anesthesia.
2. Identify the epidural space. Place the Docking System (Fig. J) on the end of the epidural needle and rotate it clockwise into the end of the epidural needle until a slight resistance can be felt. (Fig. 3)
3. When using the docking system, push the Espocan® Sleeved spinal needle NRFit® (Fig. A) into the epidural needle. A mark on the spinal needle (Fig. 2) indicates where the spinal needle protrudes from the back-eye of the epidural needle.
4. Puncture the dura with the spinal needle. The puncture of the dura can usually be felt as a significant change in resistance.
5. Remove the stylet of the spinal needle slowly and evenly. When the spinal needle is in the correct position, there will be a free flow of CSF. In this case spinal anesthesia can be performed. If there is not a free flow of CSF, the spinal needle should be pulled back under constant observation. If pulling the needle back does not produce a free flow of CSF, remove the spinal needle completely and carefully correct the position of the epidural needle or make a new puncture.
6. Secure the spinal needle in the epidural needle after making a puncture (a free flow of CSF) by rotating the docking system (Fig. J) clockwise as far as possible. This secures the spinal needle to prevent unintentional movement although allowing rotation. If not using the docking system, secure the spinal and epidural needles by appropriate measures.
7. Subsequently, perform local anesthesia.
8. Loosen the docking system (Fig. J) from the end of the epidural needle by turning it anti-

clockwise. Remove the spinal needle and the docking system together. If the docking system is not used, the spinal needle can be removed directly from the epidural needle.

9. Use the threading assist guide (Fig. E) (Fig. 1) and insert the Perifix® Catheter NRFit® (Fig. D) through the epidural needle into the epidural space and up to the required position. Observe whether there is a free flow of CSF (intrathecal position) or a backflow of blood (intravascular position) and whether paresthesia occurs. The catheter is marked on the tip and starting at 55 mm (59 mm for SoftTip) in distances of 10 mm. The 100, 150, and 200 mm positions are marked with 2, 3, and 4 bands respectively. The position when the catheter tip protrudes from the epidural needle is marked with a broad band (Fig. D*). The catheter openings (holes) are located within the first 40 mm of the catheter tip.
10. Pull the epidural needle back over the catheter completely, then insert the catheter as far as possible into the Perifix® Catheter Connector NRFit® (Fig. F) and close (Fig. 4–7). The catheter connector works according to the Click & Ready principle.
 - Open the catheter connector lid as far as possible (Fig. 4).
 - Insert the catheter into the funnel-shaped opening on the connector tip (Fig. 5).
 - Push the catheter well into the connector until it clearly reaches a stop (Fig. 5).
 - Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 6).
 - Remove air from the Perifix® Filter NRFit® 0.2 µm (Fig. G, if provided) using physiological saline solution.
 - Remove closing cone from connector, attach the connector to the Perifix® Filter NRFit® with a twist action (NRFit® Lock) (Fig. 7).
11. Rinse the system with 1–2 ml physiological saline solution. If necessary: Re-opening of the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 8).
12. Secure the filter to the patient using the B. Braun Perifix® PinPad (Fig. H).

Epidural catheter:

A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
 - Maximum spatial field gradient of 2,500 G/cm (25 T/m)
 - Maximum force product of 50,000,000 G²/cm (50 T²/m)
 - Maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) without time limit in continuous scanning.

Test dose to check the position of the catheter

A test dose cannot be administered immediately if spinal anesthesia has been performed beforehand. If the effects of epidural anesthesia need to be extended or prolonged during surgery, small quantities of the local anesthetic can be administered at intervals of 3 to 5 minutes until the desired level of anesthesia is reached. During the

action of spinal anesthesia, the effect of local anesthetics administered through the epidural needle is greater and the time to act faster than with epidural anesthesia alone. If surgery was finished without a subsequent epidural injection, allow the effect of spinal anesthesia to decline before administering the usual test dose.

Duration of use

The Perifix® Catheter NRFit® is intended for short-term use (< 30 days). Maximum intended application duration is 30 days as defined by the biological evaluation. Standard application usually ranges from 2–4 days for postoperative regional anesthesia and may be up to few weeks for Analgesia in pain therapy. However, duration of use is dependent on the intended therapy, healthcare workers discretion and individual patient conditions in accordance with hospital protocols and/or national guidelines.

Inspect the insertion site of the catheter daily. Remove the catheter in case of any signs of local or systemic inflammation of unknown origin and according to hospital protocols.

Disposal

Dispose products according to local regulations and/or clinical protocols.

Storage and handling conditions

Storage conditions:

Keep away from sunlight

Keep dry

Clinical benefits & Performances

The summary of safety and clinical performance is available via EUDAMED at ec.europa.eu/tools/eudamed using the following Basic UDI-DI: 40392390000023862Y and can be requested to the Manufacturer.

Notice to the user

All serious incidents related to this product must be reported to B.Braun and the competent authority in the country where the product is being used.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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AU B. Braun Australia Pty. Ltd.