
PERIFIX® Catheter Connector / PERIFIX® Catheter Connector NRFit®


en	Catheter Connector for Regional Anesthesia	4
de	Katheterkupplung für die Regionalanästhesie	6
bg	Конектор за катетър за регионална анестезия	8
cs	Katétrový konektor na regionální anestezii	10
da	Kateterkonnektor til regional anæstesi	12
el	Σύνδεσμος καθετήρα για περιοχική αναισθησία	14
es	Conector de catéter para anestesia local	16
et	Kateetriiliitmik regionaalanesteesiaks	18
fi	Katetriiliitin regionaaliiseen anestesiaan	20
fr	Raccord de cathéter pour anesthésie loco-régionale	22
hr	Konektor katetera za regionalnu anesteziju	24
hu	Katétersatlakozó regionális anesztéziához	26
it	Connettore per catetere per anestesia regionale	28
lt	Vietinės nejautos kateterio jungtis	31
lv	Katetru savienotājs, lai veiktu vietējo anestēziju	33
nl	Katheterconnector voor regionale anesthesie	34
no	Kateterkobling for regionalanestesi	37
pl	Łącznik cewnika do znieczulenia regionalnego	39
pt	Conector de cateter para anestesia regional	41
ro	Conector de cateter pentru anestezie regională	43
ru	Коннектор катетера для регионарной анестезии	45
sk	Katétrová spojka na regionálnu anestéziu	47
sl	Priključek katetra za regionalno anestezijo	49
sr	Priključak katetera za regionalnu anesteziju	51
sv	Kateterkoppling för regionalanestesi	53
tr	Bölgesel Anestezi İçin Kateter Konnektörü	55

B | BRAUN

LLDorder 6068 – lw 9301498 (vorm qw 230410)

en	Directions
de	Anwendungshinweise
bg	Указания за употреба
es	Pokyny pro použití
da	Betjeningsvejledning
el	Υποδείξεις εφαρμογής
es	Aplicaciones
et	Kasutusjuhendid
fi	Käyttöohje
fr	Conseils d'utilisation
hr	Upute
hu	Használati útmutatás
it	Modalità d'uso
lt	Įspėjimai
lv	Norādījumi par lietošanu
nl	Gebruiksaanwijzingen
no	Veiledning
pl	Przepisy stosowania
pt	Instruções de uso
ro	Instructiuni de folosire
ru	Инструкция по применению
sk	Aplikac'né pokyny
sl	Napotki za uporabo
sr	Uputstvo
sv	Användningsanvisningar
tr	Kullanım talimatları

en Instructions for use

Device description

Perifix® Catheter Connector /

Perifix® Catheter Connector NRFit®

Content (depending on version):

- Perifix® Catheter Connector 20/24G (Luer) or
- Perifix® Catheter Connector 19G (Luer) or
- Perifix® Catheter Connector 20/24G (NRFit®) or
- Perifix® Catheter Connector 19G (NRFit®)

For further technical characteristics, see product packaging.

An epidural catheter or a peripheral nerve block catheter is inserted in the catheter connector on one end. On the other end, the catheter connector can be connected to a filter, a syringe or an infusion line (anesthesia, analgesia).

Luer version: The Perifix® Catheter Connector can be used with devices which are compatible or in conformance with Luer connectors (e.g., ISO 80369-7). As risk mitigation measures for neuraxial applications with Luer connectors, the protective caps and the epidural catheter label are yellow-colored.

NRFit® version: The Perifix® Catheter Connector NRFit® can be used with devices which are compatible or in conformance with NRFit® connectors (e.g., ISO 80369-6).

Assessment of Magnetic Resonance (MR) compatibility

Perifix® Catheter Connector / Perifix® Catheter Connector NRFit®: MR safe.

NOTE: Components of a set, once connected, take the highest level of

MR compatibility assessment (e.g., the Perifix® Catheter Connector - Perifix® Catheter system is considered to be MR conditional).

Sterility

Sterilized using ethylene oxide.

Do not use if package is damaged. Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Catheter Connector for Regional Anesthesia.

Indication

Regional Anesthesia procedures (e.g., Epidural Anesthesia, Combined Spinal & Epidural Anesthesia, Continuous Peripheral Nerve Blocks, etc.).

Patient population

Patient can be every person that requires Regional Anesthesia. Specific patient population of the selected Regional Anesthesia method must be taken into account.

Intended User

The product should only be used by healthcare workers that have been adequately trained in the technique of the Regional Anesthesia applied.

Contraindications

There are no device-specific contraindications for the Perifix® Catheter Connector. General contraindications of the particular Regional Anesthesia technique apply:

Epidural and CSE technique:

- **Absolute contraindications** include patient refusal, severe uncorrected hypovolemia, increased intracranial pressure, infection at the site of injection and known hypersensitivity to local anesthetics.

- **Relative contraindications** include coagulation disorders, anticoagulant therapy, anatomical abnormalities of the vertebral column, chronic low back pain, neurological diseases (e.g., nerve lesions, polyneuropathy, multiple sclerosis), cardiovascular diseases (inability to increase cardiac output) and sepsis.
- **Additional contraindications** for cervical epidural anesthesia are sympathetic blockade in patients with congestive heart failure and severe aortic or mitral stenosis and congestive lung disease.

PNB technique:

- Known pre-existing neurological diseases
- Infection / inflammation at the site of insertion
- Excessive preoperative anxiety (not responsive to interventions) and inability to tolerate positioning
- Pre-existing coagulopathies, either endogenous (e.g., thrombocytopenia) or iatrogenic (e.g., warfarin therapy) due to increased risk of hematoma
- Anatomic anomalies that make identification of physical landmarks difficult
- Hepatic disease that may interfere with clearance
- Failure to obtain patient's consent
- Stimulation is not to be carried out during defibrillation or cardioversion
- Contraindications for axillary blocks: lymphangitis and lymphadenectomies

Please consult medical textbooks for standard contraindications to performing continuous Regional Anesthesia procedures.

Residual Risks/Side effects

No specific side effects for the device itself.

The general side effects and known risks of peripheral and epidural anesthesia techniques apply and need to be considered.

Precautions

- Use strictly aseptic techniques.
- Do not use if package is damaged.
- Do not use if product is incomplete or damaged.
- Visually inspect the sterile barrier system for breaches and integrity prior to use.
- Monitor all connections of the catheter set-up for tightness and disconnection during whole anesthesia and analgesia procedure (e.g., pain treatment).
- General precautions of the particular Regional Anesthesia technique apply.
- Additional for Epidural Anesthesia:
 - Trained staff and sufficient emergency equipment, in particular for artificial ventilation and hemodynamic support, must be nearby.
 - During and after cervical epidural anesthesia, cardiovascular and respiratory monitoring is important to exclude cardio-pulmonary insufficiency.

Warning

- Single use only. Do not re-use. Do not sterilize.
- Products with NRRFit® connector are not compatible with Luer devices. They are only suitable for products with NRRFit® connector!
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination

and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

- The catheter must be inserted into the connector as far as possible as injection is impossible otherwise and the catheter can become disconnected.
- Use the properly stored product by the date stated on the packaging.
- In order to achieve the optimum connection between the catheter connector and other medical devices or accessories, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.

Directions for using Perifix® Catheter Connector (incl. NRRFit®-versions)

The Catheter Connector works according to the Click & Ready principle.

- a) Open the catheter connector lid as far as possible (Fig. 1).
- b) Insert the catheter into the funnel-shaped opening on the connector tip (Fig. 2).
- c) Push the catheter well into the connector until it clearly reaches a stop (Fig. 2).
- d) Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 3).
- e) After removal of the protective cap, attach the connector to the Perifix® Filter (if provided) with a twist action (Fig. 4).

If necessary: Re-opening the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 5).

Duration of use

Duration of contact: Perifix® Catheter Connector (incl. NRRFit®) - 30 days maximum.

Disposal

Dispose products according to local regulations and/or clinical protocols.

Storage and handling conditions

Keep away from sunlight. Keep dry.

Notice to the user

All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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Perifix® Catheter Connector / Perifix®

Catheter Connector NRFit®; MR-sicher.

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