Perifix® ONE Paed Set

NRFit®

GB	Set for continuous epidural anaesthesia	6
DE	Set zur kontinuierlichen Epiduralanästhesie	7
BG	Комплект за непрекъсната епидурална анестезия	8
CZ	Set pro podávání kontinuální epidurální anestézie	9
DK	Sæt til kontinuerlig epiduralanæstesi	10
EE	Pideva epiduraalanesteesia komplekt	11
ES	Kit para anestesia epidural continua	12
FI	Jatkuvan epiduraalanestesian antosarja	13
FR	Set pour anesthésie péridurale continue	14
GR	Σετ για συνεχή επισκληρίδιο αναισθησία	15
HR	Set za kontinuiranu epiduralnu anesteziju	16
HU	Készlet folyamatos epidurális érzéstelenítéshez	17
IT	Set per anestesia epidurale continua	18
LT	Nuolatinės epidurinės nejautos rinkinys	19
LV	Komplekts, lai veiktu nepārtrauktu epidurālo anestēziju	20
NL	Set voor continue epidurale anesthesie	21
PL	Zestaw do ciągłego znieczulenia zewnątrzoponowego	22
PT	Definido para anestesia epidural continua	23
RO	Set pentru anestezie epidurală continuă	24
SE	Sats för kontinuerlig epidural anesthesi	25
SI	Set za neprekinjeno epiduralno anestezijo	26
SK	Súprava na podávanie kontinuálnej epidurálnej anestézie	27

B | BRAUN

LOT

STERILE EO

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











GB Instructions for use

Content

1. Perican® with NRFit® connector, Epidural needle with Tuohy bevel (Fig. A)
 - For older children: with plastic stylet, Ø 1.3 x 50 mm (18 G x 2") with 5 mm markings
 - For younger children: with metal stylet, Ø 0.9 x 50 mm (20 G x 2") with 5 mm markings
2. Perifix® ONE Catheter with threading assist guide (Fig. B+C)
 - 20 G (for insertion into the 18 G needle) Ø 0.85 x 1000 mm
 - 24 G (for insertion into the 20 G needle) Ø 0.6 x 720 mm
3. Perifix® Catheter Connector with NRFit® connector (Fig. E)
4. Perifix® Filter 0.2 µm with NRFit® connector and PinPad (Fig. F+G)
5. L.O.R. Syringe with NRFit® connector, latex-free (Fig. D)
6. Single-use syringe with NRFit® connector and single-use filter needle with NRFit® connector (Fig. H)
7. Catheter label (Fig. I)

Important!

Perifix® ONE Paed Sets with NRFit® components are not compatible with Luer devices. They are only suitable for products with NRFit® connector!

Materials used

Stainless steel, PC, PE, PP, PA, PUR, SAN, ABS, POM, PE foam with polyacrylate adhesive, PS, Polyisoprene rubber, TPE, EPDM, epoxid resin, silica, PDMS

Indications

Epidural anaesthesia is a central neuraxial block procedure with use in surgery, obstetrics and pain control. Continuous epidural anaesthesia / caudal anaesthesia / analgesia in children. Epidural anaesthesia (lumbar, thoracic) for fast and reliable anaesthesia during surgery. If necessary, anaesthesia can be extended and prolonged during surgery. Analgesia is possible at any time after surgery.

Contraindications and restrictions on use

The product should only be used by physicians that have been adequately trained in this technique. Absolute contraindications include patient refusal, severe uncorrected hypovolemia, increased intracranial pressure, infection at the site of injection and known hypersensitivity to local anaesthetics. Relative contraindications include coagulation disorders, anticoagulant therapy, anatomical abnormalities of the vertebral column, chronic low back pain, neurological diseases (e.g. nerve lesions, polyneuropathy, multiple sclerosis), cardiovascular diseases (inability to increase cardiac output) and sepsis. Please consult medical textbooks for standard contraindications to performing continuous epidural anaesthesia procedures.

Risks

Known risks of epidural anaesthesia techniques include, in particular, hypotension, dyspnoea / apnoea, cardiac arrest, bradycardia, hypothermia, urinary retention, toxic effects of the local anaesthetic / adjuvant (please refer to the Pharma-SPCs), inadequate block, temporary or permanent neurological complications

and local or systemic infections (postdural puncture headache (PDPH), paraesthesia, meningitis / cauda equine syndrome, transverse myelitis or anterior spinal artery syndrome with permanent paralysis, epidural haematoma/epidural space infections, abscesses, skin and subcutaneous tissue infections), inadvertent puncture of blood vessels. Additional risks for caudal anaesthesia (technique of choice for epidural catheter placement in babies and children) are back pain, intraosseous injections with subsequent pelvic damage and venous air embolism, and in very rare cases the development of epidermoid tumors. Trained staff and sufficient emergency equipment, in particular for artificial ventilation and haemodynamic support, must be nearby.

Warnings

Products with NRFit® connector are not compatible with Luer devices. They are only suitable for products with NRFit® connector! Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/ or limited functionality of the device may lead to injury, illness or death of the patient.

Should not be used in patients with known hypersensitivity to any of the materials employed.

Do not use if packaging is damaged.

Do not resterilise. Do not wet Perifix® Filter 0.2 µm with disinfectants containing alcohol as this may result in cracks in the filter housing.

Do not use syringes smaller than 10-ml as they may cause the filter housing to be damaged as a result of excessive pressure from the syringe.

L.O.R. Syringe:

Special syringe for the Loss of Resistance technique for regional anaesthesia. Syringe, only suitable to a limited extent for aspiration of liquids. Do not aspirate via the mechanical stop. For hygiene reasons, do not fill the syringe until shortly before use.

Duration of application

Inspect the insertion site of the catheter daily. Remove the catheter in the case of any signs of local or systemic inflammation of unknown origin.

Directions

CAUTION:

- Use aseptic techniques.
- Never pull the catheter through the needle as it may otherwise shear off.

1. To make the puncture use the Loss of Resistance technique (technique of choice), or use the hanging drop technique or use the Macintosh balloon.
2. After making the puncture and identifying the epidural space, use the threading assist guide (Fig. 1) and insert the Perifix® ONE epidural catheter (B) through the epidural needle into the epidural space up to the required position. (The catheter has length markings starting at 2 cm and then every 1-cm. Two rings mark 10-cm, three rings 15 cm and four rings 20 cm.)
3. Pull the epidural needle back over the catheter, insert the catheter as far as possible into the Perifix® catheter connector (E) and up to the black mark on the catheter. Close connector (Fig. 2-5). The Catheter Connec-

tor works according to the Click & Ready principle.

- a) Open the catheter connector lid as far as possible (Fig. 2).
 - b) Insert the catheter into the funnelshaped opening on the connector tip (Fig. 2).
 - c) Push the catheter well into the connector up to the black mark (Fig. 3).
 - d) Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 4).
 - e) After removal of the protective cap attach the connector to the Perifix® Filter® (if provided) with a twist action (Fig. 5) If necessary: Re-opening the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 6).
4. Remove air from the Perifix® Filter 0.2 µm® (F) using physiological saline solution, connect the filter to the catheter and rinse the catheter with 1–2 ml physiological saline solution.

CAUTION:

- The catheter must be inserted into the connector as far as possible and up to the black mark on the catheter. Otherwise injection is impossible and it can become disconnected.

Do not use if the packaging is damaged.

Do not resterilise.

The catheter connector and epidural filter should not be wet with disinfectants containing alcohol as this may result in cracks in the housing.

Test dose

Before any epidural injection, it is recommended that an adequate test dose is administered. This test dose may contain a suitable cardiovascular active substance in order to check whether the catheter is inadvertently in a blood vessel. Monitor heart rate continuously in order to detect complications, i.e. tachycardia, immediately. The catheter allows the continuous or intermittent administration of compounds that have been approved for epidural anaesthesia to the epidural space over a period of hours or days.

CAUTION:

Do not use an adhesive spray which contains organic solvents to fix the catheter.

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B. Braun Melsungen AG