Perifix[®] ONE Paed Set

NRFit[®]








en	Set for continuous epidural anaesthesia	4
de	Set zur kontinuierlichen Epiduralanästhesie	6
bg	Набор Перификс ЕДНО – Набор за непрекъсната епидурална анестезия	8
cs	Souprava pro kontinuální epidurální anestézii	10
da	Sæt til kontinuerlig epiduralanæstesi	12
el	Σετ για συνεχή επισκληρίδια αναισθησία	14
es	Set para una anestesia epidural continuada	16
et	Püsiva epiduraalanesteesia komplekt	18
fi	Jatkuvan epiduraalipuudutuksen sarja	20
fr	Set pour anesthésie épidurale continue	22
hr	Set za kontinuiranu epiduralnu anesteziju	24
hu	Készlet folytonos epidurális érzéstelenítéshez	26
it	Set per anestesia epidurale continua	28
lt	Nenutrūkstamos epiduralinės anestezijos priemonių rinkinys	30
lv	Ilgstošās epidermālās anestēzijas komplekts	32
nl	Set voor continue epidurale anesthesie	34
no	Sett for kontinuerlig epiduralanestesi	36
pl	Zestaw do ciągłego znieczulania zewnątrzoponowego	38
pt	Conjunto para anestesia epidural, continua	40
ro	Set pentru anestezie epidurală continuă	42
sk	Súprava na kontinuálnu epidurálnu anestéziu	44
sl	Set za kontinuirno epiduralno anestezijo	46
sv	Set för kontinuerlig epiduralanestesi	48

B | BRAUN

en Line drawings of the individual components
de Strichzeichnungen der Komponenten
bg Скици на отделните съставни части
cs Perokresby jednotlivých částí
da Stregtegninger af de enkelte dele
el Γραμμικά σχέδια των ξεχωριστών στοιχείων
es Esquemas de los componentes individuales
et Üksikute komponentide skeemid

fi Osien merkit
fr Croquis des différents constituants
hr Prikaz pojedinih komponenti
hu Az egyes alkotórészek vonalas rajzai
it Schemi dei singoli componenti
lt Kiekviena detalė pažymėta brūkšniu kodu
lv Atsevišķo materiāla sastāvdaļu svītrkodi
nl Tekeningen van de verschillende onderdelen







no Strektegninger av de enkelte bestanddelerne
pl Rysunki pojedynczych części składowych
pt Traços da marca de cada componente
ro Desene liniare ale componentelor
sk Perokresby jednotlivých částí
sl Shematski prikazi posameznih sestavnih delov
sv Bilder på de enskilda komponenterna

en Directions
de Anwendungshinweise
bg Указания за употреба
cs Pokyny pro použití
da Betjeningsvejledning
el Υποδείξεις εφ'αρμογής
es Aplicaciones
et Kasutusjuhendid

fi Käyttöohje
fr Conseils d'utilisation
hr Upute
hu Használati útmutatás
it Modalità d'uso
lt Įspėjimai
lv Norādījumi par lietošanu
nl Gebruiksaanwijzingen

no Bruk
pl Przepisy stosowania
pt Instruções de uso
ro Instrucţiuni de folosire
sk Aplikačné pokyny
sl Napotki za uporabo
sv Användningsanvisningar

en Instructions for Use

Content

- **Perican® NRFit®**, epidural needle with Tuohy bevel (Fig. A)
Option 1: For older children: with plastic stylet, Ø 1.3 x 50 mm (18 G x 2") with 5 mm markings
Option 2: For younger children: with metal stylet, Ø 0.9 x 50 mm (20 G x 2") with 5 mm markings
- **Perifix® ONE Catheter NRFit®** (Fig. B – *example only*) with threading assist guide (Fig. C)
Option 1: 20 G (for insertion into the 18 G needle) Ø 0.84 x 1000 mm
Option 2: 24 G (for insertion into the 20 G needle) Ø 0.59 x 720 mm
- **Perifix® Catheter Connector NRFit®** (Fig. E)
- **Perifix® Filter NRFit®** 0.2 µm (Fig. F)
- **Perifix® PinPad** (Fig. G)
- **Perifix® LOR Syringe NRFit®**, latex-free (Fig. D)
- **Omnifix® syringe NRFit®** and **Sterifix® filter needle NRFit®** (Fig. H)
- **Epidural Catheter label** (Fig. I)

For further technical characteristics such as needle and catheter size, see product packaging. Epidural needles, catheters and other medical devices or accessories can be used with devices which are compatible or in conformance with NRFit® connectors (e.g. ISO 80369-6).

Perifix® ONE Sets with NRFit® components are not compatible with Luer devices. They are only suitable for products with NRFit® connector!

External devices can be attached to the female connector of the **Perifix® Filter NRFit®**, which is the point of access to the applied **Perifix® Catheter NRFit®** version. Main applications will be syringes for aspirating and infusing anesthetics or connecting elastomeric infusion pumps respectively extension lines of syringes intended for infusion pumps for continuous infusion of anesthetics or analgesics.

Assessment of Magnetic Resonance (MR) compatibility

Perifix® LOR syringe NRFit® / Omnifix® syringe NRFit® / Perican® NRFit® epidural needle / Sterifix® filter needle NRFit®: The Device has not been evaluated for safety in the MR environment, since it is not expected to enter the MR environment in a typical application. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

Perifix® Catheter Connector NRFit® / Perifix® Filter NRFit® / Perifix® PinPad: MR safe.

Perifix® ONE Catheter NRFit®: MR conditional (Please check Operating Instructions section for additional conditions).

The device contains Cobalt (Co) as indicated on the packaging.

Following components of this set contain cobalt (Co): **Perican® epidural needles 16G, 17G and 18G (if > 80 mm) NRFit®.** Based on the biological safety evaluation and the risk benefit analysis, no risk for the patient/user is expected and no precautionary measures in this regard are needed.

Sterility

Sterilized using ethylene oxide.
Do not use if package is damaged.

Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Set for continuous epidural anesthesia

Indication

Used for epidural anesthesia (lumbar, thoracic) during surgery and epidural analgesia for regional pain treatment.

- Thoracic placement is more difficult and riskier than lumbar placement, particularly in the pediatric population.
- Caudal anesthesia (CA) is the technique of choice for epidural catheter placement in babies and children.

If necessary, anesthesia can be prolonged during surgery. Analgesia is possible at any time after surgery (pain treatment).

Patient population

Pediatric patients

Intended User

The products should only be used by healthcare workers that have been adequately trained in this technique.

Contraindications

Known contraindications include in particular the following:

Absolute contraindications include patient refusal, severe uncorrected hypovolemia, increased intracranial pressure, infection at the site of injection and known hypersensitivity to local anesthetics.

Relative contraindications include coagulation disorders, anticoagulant therapy, anatomical abnormalities of the vertebral column, chronic low back pain, neurological diseases (e.g. nerve lesions, polyneuropathy, multiple sclerosis), cardiovascular diseases (inability to increase cardiac output) and sepsis.

Caution: Do not perform cervical epidural anesthesia for pediatric anesthesia!

Please consult medical textbooks for standard contraindications to performing continuous epidural anesthesia procedures.

Residual Risks/Side effects

Known risks of epidural anesthesia techniques include, in particular, hypotension, dyspnea/apnea, cardiac arrest, bradycardia, hypothermia, urinary retention, toxic effects of the local anesthetic/adjunct (please refer to the Summary of Product Characteristics (SPC) of the drug), inadequate block, temporary or permanent neurological complications and local or systemic infections [postdural puncture headache (PDPH), paresthesia, meningitis/cauda equina syndrome, transverse myelitis or anterior spinal artery syndrome with permanent paralysis, epidural hematoma/epidural space infections, abscesses, skin and subcutaneous tissue infections], inadvertent puncture of blood vessels.

Additional risks for caudal anesthesia (technique of choice for epidural catheter placement in babies and children) are back pain, intraosseous injections with subsequent pelvic damage and venous air embolism, and in very rare cases the development of epidermoid tumors.

Technical complications with epidural catheters:

- Catheter dislodgement, migration, bent or knotted catheters, which may lead to block failures or occlusion
- Shearing off of catheter fragments due to trauma during removal/insertion

Precautions

- Trained staff and sufficient emergency equipment, in particular for artificial ventilation and hemodynamic support, must be nearby.
- Monitor all connections of the catheter set-up for tightness and disconnection during whole anesthesia and analgesia procedure (e.g. pain treatment).
- Use strictly aseptic techniques.
- Do not use if package is damaged.
- Do not use if product is incomplete or damaged.
- Visually inspect the sterile barrier system for breaches and integrity prior to use.

Warning

- Single use only. Do not re-use. Do not re-sterilize.
- Products with NRFit® connector are not compatible with Luer devices. They are only suitable for products with NRFit® connector!
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness, or death of the patient.
- Never pull the catheter through the needle as it may shear off otherwise!
- The tip of the needle becomes blunt or can even split after bone contact. The needle might bend and eventually break when applying too strong force during advancement. Needles with damaged tip increase the risk of intrathecal and intravascular placement.
- Do not wet **Perifix® Filter NRFit®** 0.2 µm with disinfectants containing alcohol as this may result in cracks in the housing.
- Do not use syringes smaller than 10 ml as they may cause the filter housing to be damaged as a result of excessive pressure from the syringe.
- The catheter must be inserted into the connector as far as possible. Injection is impossible otherwise and the catheter can become disconnected.
- Do not use an adhesive spray which contains organic solvents to fix the catheter.
- Inspect the insertion site of the catheter daily. Remove the catheter in the case of any signs of local or systemic inflammation of unknown origin and according to hospital protocols.
- In order to achieve the optimum connection between catheter, epidural needle, and other medical devices or accessories, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.
- **Perican® NRFit®**, epidural needle / **Sterifix® filter needle NRFit®:** The device is evaluated as "MR unsafe". A patient with this device shall not be scanned in an MR system.
- **Perifix® ONE Catheter NRFit®:** MR conditional (Please check Operating Instructions section for additional conditions).

Perifix® LOR Syringe NRFit®: Special syringe for the "Loss of Resistance" (LOR) technique for regional anesthesia. Syringe, only suitable to a limited extent for aspiration of liquids. Do not aspirate beyond the mechanical stop. For hygiene reasons, do not fill the syringe until shortly before use.

Operating Instructions

☞ Illustrations on page 2/3.

Directions for using Perifix® ONE Paed Sets NRFit®

1. To make the puncture use the "Loss of Resistance" (LOR) technique (technique of choice) or use the hanging drop technique or use the Macintosh balloon.
2. After making the puncture and identifying the epidural space, use the threading assist guide (Fig. C) (Fig. 1) and insert the Perifix® ONE catheter NRFit® (Fig. C) through the epidural needle into the epidural space up to the required position. The catheter has a marking at the tip and length markings starting at 20 mm and then every 10 mm. Two rings mark 100 mm, three rings 150 mm and four rings 200 mm. The position when the catheter tip protrudes from the epidural needle is marked with a broad band (Fig B*).
Example for Perican® epidural needle 50 mm: The distance from the catheter tip to the stop block is 85 mm.
The catheter openings are located within the first 40 mm of the catheter tip.
3. Pull the epidural needle back over the catheter completely, then insert the catheter as far as possible into the Perifix® Catheter Connector NRFit® (Fig. E) and up to the black mark on the catheter. Close connector (Fig. 2-5). The Catheter Connector works according to the Click & Ready principle.
 - Open the catheter connector lid as far as possible (Fig. 2).
 - Insert the catheter into the funnel-shaped opening on the connector tip (Fig. 3).
 - Push the catheter well into the connector up to the black mark (Fig. 3).
 - Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 4).
 - Remove air from the Perifix® Filter NRFit® 0.2 µm (Fig. F, if provided) using physiological saline solution.
 - Remove closing cone from connector, attach the connector to the Perifix® Filter NRFit® with a twist action (NRFit® Lock) (Fig. 5).
4. Rinse the system with 1–2 ml physiological saline solution. If necessary: Re-opening the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 6).
5. Secure the filter to the patient using the B. Braun Perifix® PinPad (Fig. G, if provided). In case the set does not contain a PinPad, please secure the filter according to your hospital protocols.

Epidural catheter:

A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
 - Maximum spatial field gradient of 2,500 G/cm (25 T/m)
 - Maximum force product of 50,000,000 G²/cm (50 T²/m)
 - Maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) without time limit in continuous scanning.

Test dose to check the position of the epidural catheter

Before any epidural injection, it is recommended that an adequate test dose is administered. This test dose may contain a suitable cardiovascular active substance in order to check whether the catheter is inadvertently in a blood vessel. Monitor heart rate continuously in order to detect complications, i.e. tachycardia, immediately. The catheter allows the continuous or intermittent administration of compounds that have been approved for epidural anesthesia to the epidural space over a period of hours or days.

Duration of use

The Perifix® Catheter NRFit® is intended for short-term use (< 30 days). Maximum intended application duration is 30 days as defined by the biological evaluation.

Standard application usually ranges from 2–4 days for postoperative regional anesthesia and may be up to few weeks for Analgesia in pain therapy. However, duration of use is dependent on the intended therapy, healthcare workers discretion and individual patient conditions in accordance with hospital protocols and/or national guidelines.

Inspect the insertion site of the catheter daily. Remove the catheter in case of any signs of local or systemic inflammation of unknown origin and according to hospital protocols.

Disposal

Dispose products according to local regulations and/or clinical protocols.

Storage and handling conditions

Storage conditions:

Keep away from sunlight

Keep dry

Clinical benefits & Performances

The summary of safety and clinical performance of the product is available in the European database on medical devices (EUDAMED) at ec.europa.eu/tools/eudamed using the following Basic UDI-DI: 40392390000023842U

Notice to the user

All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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AU B. Braun Australia Pty. Ltd.
Level 5, 7 – 9 Irvine Place
Bella Vista NSW 2153
Australia
For sterilization type see primary packaging

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CA Distributed by: B. Braun Medical Inc.
Bethlehem, PA USA 18018-3524

ID Distributed by:
PT. B. Braun Medical Indonesia
Jakarta, Indonesia.

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MY Authorised representative:
B. Braun Medical Industries Sdn. Bhd.
Bayan Lepas Free Industrial Zone,
11900 Penang, Malaysia.
www.bbraun.com

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US Rx only • Pictorial is for reference only
Distributed by:
B. Braun Medical Inc. Bethlehem,
PA 18018-3524
Phone 1-800-227-2862

Manufacturing site / Fabricado por:
B. Braun Melsungen AG,
Carl-Braun-Str. 1,
34212 Melsungen, Germany.

0824 12410802

UA.TR.001

B. Braun Melsungen AG
34209 Melsungen
Germany
www.bbraun.com