

Stimuplex® D


B | BRAUN

Peripheral

Spinal

Epidural

en	Atraumatic needle for peripheral nerve blocks	4
de	Atraumatische Kanüle für periphere Nervenblockaden	6
bg	Атравматична игла за блокиране на периферните нерви	8
cs	Atraumatická jehla pro periferní nervovou blokádu	10
da	Atraumatisk nål til perifere nerveblokader	12
el	Ατραυματική βελόνα για αποκλεισμούς περιφερικών νεύρων	14
es	Aguja atraumática para bloqueo de plexo	16
et	Atraumaatiline nõel perifeerse närvi blokaadiks	18
fi	Atraumaattinen neula ääreishermon johtopuudutukseen	20
fr	Aiguille atraumatique por blocs périphériques	22
hr	Atraumatska igla za blokade perifernih živaca	24
hu	Atraumatikus tű perifériás idegblokkhoz	26
it	Ago atraumatico per blocchi nervosi	28
lt	Atrauminė periferinių nervų blokadų adata	30
lv	Atraumatiska adata perifēro nervu blokādēm	32

nl	Atraumatische naald voor perifere zenuwblokkades	34
no	Atraumatisk nål for periferiske nerveblokker	36
pl	Bezurazowa igła do blokady nerwów obwodowych	38
pt	Agulha atraumática para bloqueios de nervos periféricos	40
ro	Ac atraumatic pentru blocurile de nervi periferici	42
ru	Игла атравматическая для проводниковой анестезии	44
sk	Atraumatická ihla pre periférne nervové bloky	46
sl	Atravmatska igla za blokade perifernih živcev	48
sr	Atraumatska igla za blokade perifernog nerva	50
sv	Atraumatisk kanyl för perifera nervblockader	52
tr	Periferik sinir bloğu için atravmatik iğne	54
zh	一次性使用外周神经丛刺激穿刺针	56

en Instructions for Use

Device description

Luer-connector

The female luer connector of the Stimuplex® D devices is a standardized connector. It can be connected to syringes with male luer connector. The syringes can be used for flushing the needle, aspirating and infusing anesthetics.

Electrode connector

The Stimuplex® D needles contain a standardized electrode connector. The electrode connector can be connected to a nerve stimulator which is intended to be used for peripheral nerve block procedures. In order to achieve the optimum connection between the stimulation needle and the nerve stimulator, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.

Assessment of Magnetic Resonance (MR) compatibility

Stimuplex® D Needle: The Device has not been evaluated for safety in the MR environment, since it is not expected to enter the MR Environment in a typical application. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

For further technical characteristics such as needle size, see product packaging.

Sterility

Sterilized by EO.

Do not use if package is damaged.

Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Atraumatic needle for peripheral nerve blocks.

Indication

Plexus and peripheral nerve blocks using electrical nerve stimulation and/or ultrasound respectively for precise detection of definite nerve areas. Nerve blocks are used as anesthesia for surgical procedures as well as for pain therapy involving the upper and lower extremities with various block techniques.

Patient population

All patients in need for anesthesia or pain therapy by plexus or peripheral nerve blocks. If, despite adequate explanation, the patient refuses plexus anesthesia, their wishes should be respected. Likewise, mentally handicapped patients and those with psychiatric problems need careful preoperative assessment.

Intended User

The product should only be used by healthcare workers that have been adequately trained in this technique.

Contraindications

- Known pre-existing neurological diseases
- Infection/inflammation at the desired site of insertion/puncturing
- Excessive preoperative anxiety (not responsive to interventions) and inability to tolerate positioning
- Pre-existing coagulopathies, either endogenous (e.g. thrombocytopenia) or iatrogenic (e.g. warfarin therapy) due to increased risk of hematoma
- Anatomic anomalies that make needle placement difficult via nerve stimulation and anatomical landmarks
- Hepatic disease that may interfere with clearance
- Failure to obtain patient's consent
- Known hypersensitivity to one or more materials used
- Contraindications for axillary blocks: lymphangitis and lymphadenectomies
- Stimulation is not to be carried out during defibrillation or cardioversion
- For information on other contraindications against providing regional anesthesia, please refer to the specialist literature.

Residual Risks/Side effects

- Nervous lesions
- Pneumothorax
- Pain and hematoma at the site of insertion of the needle
- Incomplete block
- Local anesthetics related toxicity: cardiac, neurologic, allergic
- Incorrect epidural or intrathecal injection of local anesthetic in the case of peripheral blocks in the area of the spinal column
- Infection (sepsis)

Precautions

- Once the local anesthetic has been injected, the needle cannot be repositioned as it is usually no longer possible to get the nerves stimulated and the patient no longer experiences pain in the event of the needle inadvertent intraneural placement.

WARNING

- Single use only. Do not re-use.
Do not resterilize.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited

functionality of the device may lead to injury, illness or death of the patient.

- If paresthesia is caused inadvertently as a result of direct contact with the nerve, the stimulation needle must not, under no circumstance, be inserted any further.
- An unusually high injection pressure may indicate that an injection has been administered intraneurally. Discontinue the injection, otherwise irreversible nerve damage may occur.
- In order to achieve the optimum connection between the stimulation needle and the nerve stimulator, thereby ensuring patient safety, only use B. Braun products with other B. Braun products. B. Braun Melsungen AG does not assume any liability where its products are combined with equipment from other manufacturers.
- Visually inspect the sterile barrier system for breaches and integrity prior to use.
- Do not use if package is damaged.
- Use the properly stored product by the date stated on the packaging.

Operating Instructions

- Use strictly under aseptic conditions.
- Disinfect block area thoroughly. Repeat this at least two times. Allow time to act between each application and cover with a sterile fenestrated drape.
- Position the patient and inject a skin wheal of local anesthetic at the puncture site and if necessary, infiltrate the puncture track.
- Connect the nerve stimulator's electrode cable to the stimulation needle and the skin electrode.
- Set the initial nerve stimulator parameters considering e.g. the respective block type and patient condition (e.g. diabetes).
- Prime tubing and needle with solution.
- Hold the stimulation needle by the hub and slowly advance towards the target structure.
- If clear, typical muscle contractions are felt in the target nerve's area of distribution the current is reduced until the muscle contractions are hardly visible.
- Reduce the current and optimize needle position until muscle contractions occur at lower current levels. The tip of the needle has reached the target position when noticeable contractions occur at a current of approximately 0.2 - 0.5 mA (the exact current depend on e.g. block type and patient condition (e.g. diabetes)).
- The stimulation needle is now situated an appropriate distance from the nerve.
- The stimulation needle should not be placed any nearer the nerve to avoid nerve lesions.
- Aspirate for possible intravascular placement.

- Following negative aspiration, a test dose of local anesthetic can be injected through the extension tubing.
- Once the local anesthetic has been injected, the needle cannot be repositioned as it is usually no longer possible to get the nerves stimulated and the patient no longer experiences pain in the event of the needle's inadvertent intraneural placement.

Duration of use

Sterile device for single use.

The Stimuplex® D nerve block needle must only be used for detecting the target nerve and for injecting the local anesthetic. It must then be removed.

Disposal

Disposal according to local guidelines and clinical protocols.

Storage and handling conditions

Please refer to information on packaging. Keep away from sunlight. Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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