Adaptors and Connectors

Instructions for use

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Manufactured in Mexico by:
Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093 USA

CE 2797 EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH Rep and importer:** Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK CA 0086 **UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL, UK

Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Distributed in Australia by:
Medline International, Two Australia Pty. Ltd., Sydney, Australia

RA25145 2025-1
HNL02185

REF 1077CE REF 1565CE REF 1648CE REF 1659CE
REF 1664CE REF 1665CE REF 1740CE REF 1751CE

en Adaptors and Connectors

Intended use:

Hudson RCI Aerosol Adaptors and Connectors are intended to be connected in-line with a nebulizer system.
Hudson RCI Adaptors and Connectors are intended for use in both acute and non-acute care settings by qualified medical personnel, home healthcare providers, or patients under the general supervision of qualified medical personnel.
Hudson RCI Adaptors and Connectors are intended for use on general populations, which may include pediatric and adult patients.
These devices are intended to be used for a maximum of 30 days.

Instructions for use:

REF **1077CE, 1565CE, 1648CE, 1664CE, 1665CE, 1740CE**
1. Connect the adaptor or connector in-line. Ensure all connections are secure.

REF **1659CE, 1751CE**
1. Connect the adaptor in-line.
2. Turn the blue arrow in direction medication is to be administered.
3. Insert the canister and dispense medication as prescribed.

Disposal instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warnings:

For devices with arrows, ensure all arrows are pointing in the direction of flow.
For devices with ports, close off port with cap when not in use to prevent leakage.
For devices with drainage bags, empty when full.
These devices are not intended to be used in a pressurized or closed system.
Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
Non-Sterile. Do not sterilize prior to use.
Do not use if device is soiled or damaged.
Do not use in the presence of an ignition source.
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

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