KEYSURGICAL®

Tray Liners and Soaker Sheets

EN - INTENDED USE:
Key Surgical® Tray Liners and Soaker Sheets are intended for placement on the bottom of a sterilisation tray as a means of protection and to absorb wetness and humidity during processing.

Soaker Sheets / Heavy Tray Liners may be used as shelf liners for added protection from rough handling and wet packs.

CONTRAINDICATIONS:
No known contraindications and/or adverse effects.

PRECAUTIONS:
Store away from a source of ignition.
Single use only; reuse may lead to patient cross contamination.

PREPARATIONS FOR USE:
1. Select appropriate size product for tray
 Note: Tray Liner Rolls should be cut to the required length
2. Product may be placed inside and/or under the tray
3. Place product inside the wrapped tray or rigid container system

TLP-xxx-xxx and TL-xxx-xxx products can be used in the following configurations:

Pre-Vacuum Steam/100% Ethylene Oxide (EtO)
- Up to 2 tray liners inside and 2 tray liners outside
- Up to 8 folded layers (4 tray liners)

INSTRUCTIONS FOR USE:
1. Sterilise according to the information below
2. Place product inside or under tray
3. Place instruments in tray and wrap or close container according to manufacturer's instructions

STERILISATION:
Provided non-sterile and are single use.
Not suitable for vaporized hydrogen peroxide/plasma sterilisation.
Validated for sterilisation using the following cycle parameters:

Steam Sterilisation				
Item Number		Cycle	Temperature	Minimum Time – Full Cycle
TLP-012-018 TLP-020-025 TL-101-200 TL-105-250	TL-202-500 TL-203-500 TL-204-500	Pre-Vacuum	132°C (270°F)	4 minutes
TL-101-200 TL-105-250 TL-202-500	TL-203-500 TL-204-500	Gravity	121°C (250°F)	30 minutes

100% Ethylene Oxide (EtO)				
Item Number	Concentration	Temperature	Exposure Time	Humidity
TLP-012-018 TLP-020-025 TL-101-200 TL-105-250 TL-202-500 TL-203-500 TL-204-500	735 mg/L	55° C	60 minutes	50-70%

All other products are compatible for steam sterilisation up to 134°C.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:
Dispose after use according to facility protocol. Do not reuse.

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MANUFACTURER: KEY SURGICAL LLC