ELECTROSURGICAL ELECTRODE

- en INSTRUCTIONS FOR USE
- pt INSTRUÇÕES PARA USO
- it ISTRUZIONI PER L'USO
- es INSTRUCCIONES DE USO
- el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- fr NOTICE D'UTILISATION
- nl GEBRUIKSAANWIJZING
- de GEBRAUCHSANWEISUNG
- fi KÄYTTÖOHJEET
- sv INSTRUKTIONER FÖR ANVÄNDNING
- da INSTRUKTION
- pl INSTRUKCJA OBSŁUGI
- no BRUKSANVISNING
- ar تعليمات الاستخدام
- hr **UPUTE ZA UPORABU**
- tr **KULLANIM TALİMATLARI**
- sl **NAVODILA ZA UPORABO**
- lt **NAUDOJIMO INSTRUKCIJOS**
- ro **INSTRUCTIUNI DE UTILIZARE**
- bg ИНСТРУКЦИИ ЗА УПОТРЕБА
- cs **NÁVOD K POUŽITÍ**
- et **KASUTUSJUHEND**
- hu **HASZNÁLATI ÚTMUTATÓ**
- sk **NÁVOD NA POUŽITIE**

REF SP200-L31S
REF SP200-L35S
REF SP200-L37S
REF SP200-L45S
REF SP200-L36S

2026-02

EN

Instruction for Use

Tungsten Loop LLETZ Electrode

This is a sterile single-use product and it is not designed for re-sterilization.

Electrosurgical electrodes are designed to remove tissue and control bleeding by high frequency current. Intended to be used with electrosurgical pencils.

Precautions

Check any damage on pouch and expiry date before use.

Do not use if the sterile packaging has been opened or damaged. Do not use electrode if any default/nicks are noticed.

Prior to use, check the electrode connection to the accessory you are using to ensure proper fit and compatibility.

Improper electrode installation and/or use may result in injury or burns to the patient or operating room personnel. Refer to generator's manual of instruction for power setting recommendations for various electrosurgical electrode configurations and specifications on voltage output characteristics to ensure that all safety precautions are followed.

Confirm proper electrosurgical generator setting before proceeding with surgery.Power setting guidelines may vary due to differences in surgical technique, patient, electrodes and surgical set-up.

Electrosurgical accessories that are activated or hot from use can cause fire.Do not place them or in contact with flammable materials (such as gauze or surgical drapes). Use a holster to hold electrosurgical pencils or similar accessories safely away from patients, personnel and surgical drapes. When not using active accessories, place the pencils in a clean, dry, nonconductive and highly visible area with no contact with the patient. Inadvertent contact with the patient may result in burns. Restricted to sale by the order of a physician. Safety classification BF type class I & class II.

Assembly instruction

1.Open package and place content on sterile field
2.Securely attach the electrode to the electrosurgical pencil
3.Look electrode in correct position via the hex lock facility
4.Refer to the generator instruction manual for proper use and select appropriate mode for electrode excision procedure

Notice
Loop and wire electrodes are not designed for use in coagulation mode. Do not activate the generator in coagulation while using a loop or wire electrode. Loop or wire breakage may occur.

Before inserting or changing an electrode be sure that the active accessory is not connected to an electrosurgical generator or make sure that the generator is off (standby).

Grasp the electrode by using the insulating sleeve and insert the round shank into the electrosurgical accessory until the insulating sleeve is fully inserted. If the shank and/or insulating sleeve does not fit use is not recommended.

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