REF ORISISO
REF ORISLL
REF ORIS15
REF ORISPUSH

en INSTRUCTIONS FOR USE
de GEBRAUCHSANWEISUNG
fr NOTICE D'UTILISATION
nl GEBRUIKSAANWIJZING
it ISTRUZIONI PER L'USO
es INSTRUCCIONES DE USO
pt INSTRUÇÕES DE UTILIZAÇÃO

MD STERILE EO

STERILIZE

15°C 25°C

LATEX

Latex Free - Latexfrei - Sans latex - Latexvrij - Senza lattice - Sin Látex - Sem látex

Advanced Medical Design Co., Ltd.,
4F~5F, No. 29, Wuquan 5th Rd., Wugu Dist.,
24888 New Taipei City, Taiwan.

CE 2460 EC REP OBELIS S.A,
Bd. Général Wahis, 53,
1030, Brussels, Belgium

CH REP Arazy Group Swiss GmbH,
Bruderholzallee 53, 4059 Basel, Schweiz

UK Representative: Lotus Global Co., Ltd.,
23 Maine Street, Reading, RG2 6AG, England, United Kingdom

Importer Information

Medline International B.V.,
Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Medline Industries Ltd, WA4 6HL, UK

Medline International Switzerland Sàrl,
c/o MN & Associés SA, 1 Place de Longemalle,
1204 Genève, Switzerland

RJ22AMD 2022-10

en INSUFFLATION TUBING SET


Indication for use:

Insufflation are intended to filter and insufflate CO_2 in order to establish pneumoperitoneum prior to laparoscopy procedures.

Description:

Insufflation tubing set comprises a flexible tubing and a filtering device.

Medline's insufflation tubing set is available with universal Luer-Lock connector and multiple connector onto insufflator outlet:

ORISISO, with ISO 15M/22F connector

ORIS15, with 15 mm male connector

ORISLL, with Luer-Lock connector

ORISPUSH, with flexible female connector

Instructions for use:

Select the insufflation tubing set according to the insufflation machine used.

Inspect the insufflation tubing set packaging to ensure it is not broken or damaged.

Open packaging using aseptic technique, remove the insufflation tubing set from inner bag, and inspect overall conditions and physical integrity.

Screw Luer-lock universal connector to the entry port (veress needles or trocar cannula).

Connect the insufflation tubing set to the insufflation machine through the connector closest to the filter.

Proceed with insufflation.

Warnings:

For single patient use.

Do NOT use the device if sterile conditions have been compromised or the packaging is damaged.

The product must be used exclusively by qualified medical staff. Thorough understanding of principles, instrumentation and techniques involving minimal invasive procedures is required to avoid possible injuries to both patient and medical staff.

Do NOT reuse and/or re-sterilise. Reuse of single-use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function, which may lead to injury, illness and/or death of the patient.

Any serious incident that has occurred in relation to the device should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal instructions:

Discard after use. Precautions should be taken when discarding this device. Disposal of the device should be done in accordance with applicable hospital or national regulations for biologically hazardous waste.

Storage:

Store under clean conditions between 15°C~25°C. Avoid exposure to sunlight and humidity.

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