Steriset Containers (p.2-5)

REF TVAxx series: Container with Aluminium Lid
REF TVExx series: Container with Stainless Steel Lid
REF KVAxx series: Mini TASKIT Container with Aluminium Lid

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Instruction for use

en The Steriset System

Steriset containers are reusable, rigid sterilization containers. They are intended to be used to enclose another medical device that is to be sterilized by a healthcare provider and also to maintain sterility until used. The Steriset container is designed for use with vacuum steam sterilization systems.

⚠

Hot-air sterilization, gravity or circulation procedures, formaldehyde, ethylene oxide sterilization or other substitute procedures for the sterilization of thermolabile products such as plasma sterilization or peroxide sterilization may not be used.

⚠

Steriset containers must be used only in combination with components or spare parts (e.g. seals) manufactured by Medline for the Steriset system.

Valve System:

Steriset valve containers are containers with a closed (unperforated) base or a thermostatic condensate drain with a permanent stainless-steel pressuresensitive valve in the inner lid. These valves open and close in response to the changes in pressure during vacuum steam sterilization cycle.

- During the vacuum phase, the valve opens upwards, and the air/steam mixture flows out of the container.
- During the pressurization phase, the valve opens downwards and thus allows steam to flow in. The system is automatically flushed and sterilized by the steam rushing through the valve during each sterilization cycle.
- After sterilization, the valve system is closed and a barrier to microorganisms.

After 5000 cycles of accelerated aging in accordance with EN 868-8 Annex G, the tested (seal/gasket function and barrier properties) Steriset Containers showed no functional changes. The end of serviceable life tends to be determined by mechanical wear and tear and damage rather than frequent use of the container for its intended purpose and can be identified by the required visual and functional inspection before reuse (see Visual & Functional Inspection section).

Maintenance:

Only use Steriset containers following the instructions in this user manual and make sure this manual is easily accessible for relevant personnel.

Before using the container or opening the sterile container, or after any accident (i.e. dropped), it is essential that the sterile container undergoes a thorough visual and functional check for damage to ensure its intended use can still/has been met.

⚠

Careless handling or the use of inappropriate chemicals can cause damage to the sterilization containers, thereby putting at risk the ability to attain and preserve sterility.

Example: never use disinfection solutions containing of halogen or chloride – there is a danger of corrosion even for stainless steel containers. Steriset containers therefore require regular visual, and if necessary, functional checks.

Visual & Functional Inspection:

• Shape

The upper rim of the container bottom, where the seal sits, and the inner lid, where the seal is inserted, should be free of dents and visible deformations. The outer lid, inner lid or base should not show noticeable buckling, holes, dents or dimensional distortion.

• Seal

The seal of the inner lid must be completely inserted and undamaged.

• Handles and Latches

Handles, latches and similar fittings (e.g. rivets) must not be loose. The latches should lock the lid firmly to the bottom of the container.

• Inner Lid Valve & Valve Covers

Valve covers (the perforated disc section in image 3) should not show visible malformations. Use thumb to unlatch the inner lid from outer lid (see images 1-2). If necessary, use inner lid knob(s) to remove inner lid from outer (see image 3). Valves should not show visible malformations. Press knobs on valves and let go suddenly to check for spring back (see image 4-5). Use finger to pull upward on valve and let go suddenly to check for spring back (see image 6).

• Condensate Drain Valve (where applicable)

Remove drain by turning counter clockwise and lifting. Visually inspect drain for undamaged seal, cleanliness and foreign matter (see image 7). Re-Insert the drain by pressing vertically downwards while turning clockwise until it clicks into place and can no longer turn (see image 8). Check the drain's spring tension by thumb pressure while in the screwed in position. An audible and/or visual spring back should be observed.


Tamper Proof Sealing:

In accordance with EN ISO11607-1, Steriset containers offer two tamper-proof sealing methods:

1. Disposable Plastic Seals

Applicable only for the Taskit (mini) valve containers, but once attached, can be opened only by breaking. They are inserted through the horizontal hole in the closure before sterilization, and fastened between thumb and index finger. To unseal/open: Insert, for example, an index finger and twist (do not pull – see images 9– 11).

If seals purchased from another company other than Medline are to be used, the user must make sure that the containers are sealed in such a way as to prevent inadvertent opening of containers and to ensure that the tampering is evident.

2. Automatic Seal (ThermoLoc)

Containers with this type of closure indicate sealing automatically as a result of the heat during sterilization procedure. In the sealed condition, a red block appears below the closure. To unseal, push the block vertically upwards until the red block has disappeared into the housing. Only then can the clasp be released (see images 12 –14).

The sealing function of ThermoLoc can be "reset" (pushed upward without automatically reappearing downward) after cool down (approximately 10 minutes after removal from the sterilizer). Attempts to reset the seal before sufficient cool down will require an application of force and gloves to be worn due to the container still being hot. This may result in the ThermoLoc reappearing after the lid latch has been opened to remove the contents of the container.

⚠

The Thermoloc is not a substitute for external chemical indicator monitoring (see External/Internal Indicators under Sterilization section).

Internal Packaging

Use of absorbent surgical towels/linen or disposable equivalent (e.g. tray liner) is recommended to assist in the final drying stage.

The recommended minimum drying time for pre-vacuum steam sterilization is 20 minutes.

Recommended absorbent tray liner configuration with surgical instruments:

- Unfold tray liner and place within inner basket (see image 15).
- Place instruments over the liner within the basket. If necessary, fold the excess liner material over the instrumentation (see image 15).

Recommended absorbent towel/linen configuration with surgical instruments:

- Absorbent towel/linen can be utilized in the same manner as described with the tray liner.
- Absorbent towel/linen can be utilized as a "topper" (laid over, or on "top," of the surgical instruments). This method is the recommended configuration for absorbent towels/linens for Steriset containers with a condensate drain valve as the condensate that forms prior to the drying phase should be able to drain unimpeded (see Sterilization section). Steriset valve containers are containers with a closed (unperforated) base or a thermostatic condensate drain with a permanent stainless-steel pressuresensitive valve in the inner lid. These valves open and close in response to the changes in pressure during vacuum steam sterilization cycle.

⚠

There is a risk of non-sterility if the container is overloaded or protruding corners of internal packaging prevent the container from closing properly.

⚠

Never use additional external packaging of the container (e.g. wrap) during sterilization as the resulting increase in flow resistance may affect sterilization or damage the container (e.g. implosion).

Sterilization

Operational Limits:

In the case of instrument sterilization, the inner content weight (including perforated tray) should not exceed 10 Kg. If exceeded, residual moisture may remain even if internal packaging materials to assist in drying are utilized. In the case of cloth loads (or similar), the load weight should not exceed 6.5 - 7 Kg.

- To prevent damage to the parts of the container or inner contents, it is recommended that the container be transported with its lid latched (sealed).
- As a general rule, containers should only be handled after cool down (DIN 58953/9 recommends 30 minutes).

Due to the double-lid design, Steriset containers can be sterilized while stacked.

- Stack height should be less than or equal to 60 cm.
- To prevent accidents or mechanical damage, do not handle stacked containers.


The containers should be placed horizontally in the sterilizer to decrease probability of retained moisture due to condensation build on one side of the container.

⚠

If the sterilization procedure causes sterilization containers to become deformed in any way, then there is no guarantee of sterility. In such cases, the entire load should not be used and an investigation initiated to determine the root-cause (e.g. analysis of sterilization records, inspection of sterilizer and other sterilization packaging, functional testing/review of damaged Steriset).

External/Internal Indicators

External process indicators are used to indicate that the container has been exposed to the sterilization process and distinguishes between processed and unprocessed containers. Use of external indicators should be in accordance with the facility's policies and procedures (See also DIN 58953/9). Medline indicator label (see image 16), is properly sized to be inserted into the slot on the outside of the Steriset container.

The internal biological and/or chemical indicator/integrator should be placed where the user has determined the most challenging location within the container. Use of internal indicators/integrators should be in accordance with the facility's policies and procedures.

⚠

If external or internal indicators are not utilized, then other organizational measures should be taken to ensure non-sterile Sterisets are not inadvertently released.

Condensation Drain Valve Bottom

Many factors contribute to the effectiveness of drying time with some being unique to the individual facility. Steriset containers fitted with a condensate drain valve in the bottom can assist in drying effectiveness particular situations where retained moisture may be of concern to the user (e.g. heavy, dense instruments). The drain valve is thermostatic, only opening when the temperature rises above 130°C, thereby allowing condensate to drain through its central outlet during sterilization. During the drying phase, while still inside the sterilizer, the drain closes when the temperature first drops below 110°C (typical removal temperatures from the sterilization chamber are between 75–80°C.

For optimal effectiveness of the condensate drain valve, the condensate should be able to drip to the container base unimpeded (see Internal Packaging section).

⚠

121°C program: standard condensation valves do not open during the 121°C sterilization program for the above reasons.

⚠

Steriset containers with drain valves must be removed from the sterilizer within 1 hour of program end. If the chamber door remains closed for a prolonged period of time, the temperature within may increase enough to reopen the drain valve, presenting a sterility risk upon the containers removal from the sterilizer. In the case of small load sterilization (<2 Kg), the container must be removed from the sterilizer promptly at program end.

Post-Sterilization

Cooling

To safeguard against accidents (burns, dropping, etc.), terminally sterilized items should not be touched during the cooling process. This should be followed despite the Steriset system's thermally insulated handles (silicone coating).

Steps should be taken to ensure cooling does not take place too rapidly after removal from the sterilizer (e.g. containers should remain on sterilizer cart during cooling and away from air conditioning and cold air-vents). Rapid cooling can lead to re-condensation of water vapor inside the container resulting in accumulated retained moisture.

Storage/Transport

The Steriset system is intended to allow sterilization of the enclosed medical device and also to maintain sterility until used.

Sterilized items should be stored under environmentally controlled conditions with appropriate handling practices to reduce the potential for contamination as well as kept clean and dry.

According to scientific research, (see also DIN 58953/9;2000), the potential contamination risk can be reduced by the following measures.

- The use of internal packaging.
- Storage under (dust)-protected conditions.
- Limitation of storage period.

We recommend to follow DIN 58953/9;2000 - operation procedure for sterilization containers.

Aseptic Presentation

The contents of Steriset can only be considered sterile if the container has been correctly sterilized and not tampered with (as indicated by the applicable tamperproof seal) prior to its intentional opening at the point of use.

If Steriset containers are being opened after a long period of storage or after storage under non-controlled environmental conditions, then it is recommended that the outer bio-barrier be wiped with a disinfectant before handling (see Automatic & Manual Cleaning/Disinfection section for considerations).

Inspection should also be conducted at the point of use to ensure that no visible mechanical damaged occurred during handling (See Visual & Functional Inspection section).

⚠

The outer lid latches of the container must be opened so that lid is opened parallel with the container bottom. Opening only one side/latch with the other still closed risks damage to the latch that is still closed and the integrity of the container (see images 17-18).

Decontamination and Cleaning

Initial Receiving

Inspect, thoroughly clean, and rinse all reusable components before placing into service. Not properly preparing your Steriset system per the following cleaning instructions may adversely affect the containers anodized finish and performance.

Steriset Disassembly

After aseptic presentation, direct contamination of containers typically does not occur if removed from the operating room before the patient enters. For these instances, the facility should determine best cleaning practice based on its established policy and procedures.

If contamination cannot be excluded containers must be fully cleaned/decontaminated in consideration of the below instructions.

Automatic & Manual Cleaning/Disinfection

Steriset TVA and KVA container bottoms, inner and outer lids are made of anodized aluminum. The outer lids of TVE container, fittings and valves of all types of Steriset containers are made of electrolytically polished, chemically resistant stainless steel.

When selecting cleaning and disinfection agents and methods it should be ensured that they recommended/compatible with anodized aluminum as well as the following considerations should be made:

- Do not use lathering cleaning substances (e.g. powder) or abrasive metal brushes or wipes on container surfaces.

· Do not use substances containing phosphates, halogens or chlorides as there is a risk of corrosion to container surfaces and parts.

· Rinsing should remove all cleaning agent residues. Desalinated or Critical water should be used for the final rinse.

· Container parts must be thoroughly dried and stored in a dry, clean area following cleaning.

· Selected cleaning agents should be appropriate for the available water quality (e.g. pH, temperature, hardness, etc.) as a highly acidic or alkaline pH (detergent and water solution) could damage the container surfaces or parts and/ or reduce the effectiveness of the cleaning process.

For thermal cleaning and disinfection (recommended):

· **In the case of desalinated or critical water:** pH-neutral enzymatic cleaner up to 60°C and subsequent thermal disinfection during final rinsing by preferred use of desalinated or critical water.

· **In case of softened water:** mild alkaline, non-chlorous detergent at approximately 55°C followed by thermal disinfection (desalinated or critical water for thermal disinfection during final rinsing is preferred).

For chemical cleaning and disinfection, the use of pH-neutral or weak acidic cleaners and:

· **In the case of desalinated or critical water:** pH-neutral final rinse aid.

· **In the case softened water:** mild alkaline, soft water compatible final rinse aid.

Inner and outer Steriset lids should be separated (see image 19), and placed diagonally within the washer. The container bottom should be inverted and face downwards. If the Steriset is fitted with a condensation drain valve, the valve should removed from the bottom during washing (see image 20).

⚠

When disinfectant substances are applied manually, the instructions for use should be followed to ensure that inadequate effect and damage to the Steriset are avoided.

Specific attention should be paid to ensure the substances selected are tolerated by the Steriset materials (e.g. anodized aluminum), while following the specified dilution rate, exposure time, and compatibility in combination with other cleaning agents/methods.

⚠

If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

⚠ Water Quality

Please note:

a) Regular tap water can only be used cold, and thus is unsuitable for washing.

b) Softened water should only be used up to max. 60 °C:
- Rinsing: only << 60°C
- Washing: up to max. 65°C – when suitable cleaner selected (see recommendations)
- Thermodisinfection: NOT applicable! (hot softened water may degrade aluminum-white deposit!)

c) Desalinated water (EN 285 quality, appendix B)
- Rinsing: suitable up to 95°C
- Washing: in combination with a suitable cleaner (see recommendations)
- Thermodisinfection: suitable up to 95°C

The brochure "Proper maintenance of instruments" of the German working group "Instrument preparation now" gives further valuable guidance for re-processing of reusable medical devices: See i.e. PDF-brochure under www.a-k-i.org

Suitable cleaners / disinfection solutions are i.e. offered by:

· Chemische Fabrik Dr. Weigert: the MediClean - Medizym - Neodisher family: **www.drweigert.de**

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· Do not use lathering cleaning substances (e.g. powder) or abrasive metal brushes or wipes on container surfaces.

· Do not use substances containing phosphates, halogens or chlorides as there is a risk of corrosion to container surfaces and parts.
· Rinsing should remove all cleaning agent residues. Desalinated or Critical water should be used for the final rinse.
· Container parts must be thoroughly dried and stored in a dry, clean area following cleaning.
· Selected cleaning agents should be appropriate for the available water quality (e.g. pH, temperature, hardness, etc.) as a highly acidic or alkaline pH (detergent and water solution) could damage the container surfaces or parts and/ or reduce the effectiveness of the cleaning process.

For thermal cleaning and disinfection (recommended):
· **In the case of desalinated or critical water:** pH-neutral enzymatic cleaner up to 60°C and subsequent thermal disinfection during final rinsing by preferred use of desalinated or critical water.
· **In case of softened water:** mild alkaline, non-chlorous detergent at approximately 55°C followed by thermal disinfection (desalinated or critical water for thermal disinfection during final rinsing is preferred).

For chemical cleaning and disinfection, the use of pH-neutral or weak acidic cleaners and:
· **In the case of desalinated or critical water:** pH-neutral final rinse aid.
· **In the case softened water:** mild alkaline, soft water compatible final rinse aid.

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19

20

⚠
When disinfectant substances are applied manually, the instructions for use should be followed to ensure that inadequate effect and damage to the Steriset are avoided.

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⚠
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

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b) Softened water should only be used up to max. 60 °C:
· Rinsing: only << 60°C
· Washing: up to max. 65°C – when suitable cleaner selected (see recommendations)
· Thermodisinfection: NOT applicable! (hot softened water may degrade aluminum-white deposit!)

c) Desalinated water (EN 285 quality, appendix B)
· Rinsing: suitable up to 95°C
· Washing: in combination with a suitable cleaner (see recommendations)
· Thermodisinfection: suitable up to 95°C

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Suitable cleaners / disinfection solutions are i.e. offered by:
· Chemische Fabrik Dr. Weigert: the MediClean - Medizym - Neodisher family: **www.drweigert.de**
· Henkel-Ecolab : the Sekumatic family: **www.ecolab.com**

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EC REP Medline International France SAS,
Medline Industries Ltd, WA4 6HL, UK
Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands