Medium Concentration Oxygen Mask

Instructions for use

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RB26SWY 2026-02

REF	Tubing Included	Length	Tubing Color	Patient Size
1041S	✓	2.1 m	Clear	Adult
1049S	-	-	-	Adult
41042S	✓	2.1 m	Clear	Pediatric
41050S	-	-	-	Pediatric
41035S	✓	2.1 m	Clear	Pediatric
41040S	✓	2.1 m	Clear	Adult
1041SG	✓	2.1 m	Green	Adult
1042SG	✓	2.1 m	Green	Pediatric
1043SG	✓	3.0 m	Green	Adult

Tubing Included; Schlauch enthalten; Tubulure incluse; Inclusief slang; tubo incluso; Tubos incluidos; Tubo incluído; Vrátane hadičiek; Včetně kroužku s háčkem; Slang medföljer; Putket mukana; Slangen medfølger; Slange følger med; Περιλαμβάνεται σωλήνας; Rurka w zestawie; Priložena cevka; Tuburi incluse; Включени тръби; Voolikud kaasa pandud; Cső mellékelve; Pridedami vamzdeliai; Cijevi su sadržane; Hortum dahildir; الأنابيب مشمولة في المنتج.

Length; Länge; Longueur; Lengte; Lunghezza; Longitud; Comprimento; Délka; Dĺžka; Längd; Pituus; Længde; Lengde; Μήκος; Długość; Dolžina; Lungime; Дължина; Pikkus; Hosszúság; Ilgis; Duljina; Uzunluk; الطول

Tubing Colour – Clear, Green; Schlauch Farbe - Transparent,Grün; Couleur de tubulure - Transparent, vert, Slangkleur - doorzichtig, groen, Colore del tubo - verde chiaro, Color del tubo - Verde transparente, Color do tubagem - Transparente, Verde, Farba hadičiek – číra, zelená, Barva hadičky – čirá, zelená, Slangens färg – Klar, Grön, Letkun väri – Kirkas, Vihreä, Slangefarve – klar, grøn, Farge på slange – Klar, grønn, Χρώμα σωλήνων – Ανοιχτό, Πράσινο, Kolor rurki – przezroczysty, zielony, Barva cevke - prozorna, zelena, Culoare tubulaturã – Transparent, verde, Цвят на тръбата – Прозрачен, Зелен, Vooliku värv – läbipaistev, roheline; Tömlő színe – Átlátszó, Vamzdelio spalva – skaidri, žalia; Boja cijevi – prozirna, zelena; Zöld, Tüp Rengi – Şeffaf, Yeşil; لون الأنبوب - شفاف، أخضر.

Patient Size - Adult, Pediatric; Patienten Größe - Erwachsene, Kinder; Taille patient - Adulte, Pédiatrique; Patiëntenmaat - Volwassenen, Pediatrisch; Misura adulto, pediatrico; Tamaño - Adulto, Pediátrico; Tamanho Peciente - Adulto, Pediátrico; Veľkosť pacienta - Dospelý, Pediatrická; Velikost pacienta - Dospělý, Pediatr; Patientstorlek - Vuxen, Pediatrisk; Potilaan koko - aikuiset, lapset; Patientstørrelse – voksen, barn; Pasientstørrelse - Voksen, Pediatrisk; Μέγεθος ασθενούς - Ενηλίκων, Παιδιατρικό; Rozmiar pacjenta – dorosły, dziecko; Velikost bolnika – odrasli, pediatrični bolniki; Dimensiune pacient - adult, pediatric; Пациентски размер - възрастни, педиатрични пациенти; Patsiendi suurus – täiskasvanu, laps; Betegméret - felnőtt, gyermek; Pacientų amžius – suaugusieji, vaikai; Veličina pacijenta – odrasli, djeca; Hasta boyutu - Yetişkin, Pediatrik; حجم المريض- بالغ، أطفال.

en Medium Concentration Oxygen Mask

Intended use:

For administration of oxygen by inhalation for patients experiencing hypoxia.

These devices are intended for use by medical personnel, lay persons, or patients that have been prescribed oxygen.

Intended to be used for a maximum of 30 days.

Instructions for use:

1. For masks provided without tubing: attach separately acquired tubing to Mask grommet.
2. Attach supply end of tubing by pushing end connector directly onto compatible gas outlet.
3. Set oxygen to prescribed flow as directed by trained clinician.
4. Check for secure connections by confirming flow of gas through device.
5. Place the mask on the patient's face with the elastic strap below the ears and around the neck.
6. Gently pull the strap ends until the mask is secure.
7. Mold the metal strip on the mask to fit the face.

Disposal Instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warning:

Prior to use, ensure all connections are secure. Do not over tighten or force connections. Do not alter device. Improper connections may result in disconnection or leakage during use, which may result in performance degradation and increase the risk of respiratory distress.

Do not use in the presence of an ignition source.

Prior to use, ensure compatibility with gas outlet. Standard connectors are compatible with 5-7mm male fittings.

Do not use if device is soiled or damaged.

Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Intended to be used for a maximum of 30 days.

Max Pressure: 50 PSI. Max Flow: 15LPM

Non-Sterile. Do not sterilize prior to use.

If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Caution should be exercised when delivering oxygen via a face mask to patients with chronic obstructive pulmonary disease or facial trauma.

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