VAPREVENT® ORAL CARE

REF MDS096513MEU

en Suction Swabs Kit

Intended Use: Oral Care products are intended to clean oral secretions from the oral cavity, suction swabs and toothbrush also are intended to aspirate oral secretions, during oral care procedures. They can be used alone or in combination with Oral Care solutions to clean the oral cavity. These devices are not intended to be used by lay users.
1. With gloved hands, peel back film. Ensure foam on suction swabs and DenTips® is intact. Connect Suction swab to standard suction line. Ensure connection is secure. 2. Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands. 3. Open Medline mouth rinse and pour into basin. Saturate suction swabs with mouth rinse. 4. Gently clean teeth and gums for approximately 1-2 minutes. Provide suction by using the thumb port. 5.Ensure foam is intact before discarding suction swab. Inspect oral cavity and remove any particles, if present. 6.Open mouth moisturizer and squeeze onto DenTips®. Gently coat inside of the mouth as indicated on the product. 7.Ensure foam is intact before discarding DenTips®. Inspect oral cavity and remove any particles, if present.

Warning: Discontinue use and consult physician if irritation occurs. Keep out of the reach of children. Protect from freezing, Avoid excessive heat.

Directions: MOUTH RINSE: Adults and children over 12 years, use approximately 1 tablespoon (7 ml) of rinse unless otherwise directed by a healthcare professional. Rinse in mouth for 30 seconds and spit out. Do not swallow. **MOUTH MOISTURIZER**: Apply gel to inside of mouth and/or lips as needed. Not for use by children 12 years and under, unless directed by a healthcare professional. Do not swallow.

Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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Packaged in Mexico for: **Medline Industries, LP,** Three Lakes Drive, Northfield, IL 60093 USA.

2797

EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH Rep and importer:** Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

0086

UK representative and importer in UK: Medline Industries Ltd, WA4 6HL,UK

Medline International B.V., Catharijnesingel 47, 3511 GC Utrecht, The Netherlands

Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.
©2026 Medline Industries, LP