VAPREVENT® ORAL CARE

REF MDS09651300EU

en Suction Swabs Kit

Intended Use:
Oral Care products are intended to clean oral secretions from the oral cavity, suction swabs and toothbrush also are intended to aspirate oral secretions, during oral care procedures. They can be used alone or in combination with Oral Care solutions to clean the oral cavity. These devices are not intended to be used by Lay users.
1. With gloved hands, peel back film. Ensure foam on suction swabs and DenTips® is intact. Connect Suction swab to standard suction line. Ensure connection is secure.
2. Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
3. Saturate suction swab with water or desired solution.
4. Gently clean teeth and gums for approximately 1-2 minutes. Provide suction by using the thumb port.
5. Ensure foam is intact before discarding suction swab. Inspect oral cavity and remove any particles, if present.
6. Open mouth moisturizer and squeeze onto DenTips®. Gently coat inside of the mouth as indicated on the product.
7. Ensure foam is intact before discarding DenTips®. Inspect oral cavity and remove any particles, if present.

Warning: Do not use for tracheostomy care. Discontinue use and consult physician if irritation or excessive bleeding occurs. Keep out of the reach of children. Protect from freezing. avoid excessive heat.
Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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©2026 Medline Industries, LP
Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.
EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France
UKCA 0086 **UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL,UK
Medline International B.V., Catharijnesingel 47, 3511 GC Utrecht, The Netherlands