Adaptors and Connectors

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Besmed Health Business Corp. No. 5, Lane 116, Wu-Kong 2nd Road, Wu-Ku District, New Taipei City 24888, Taiwan

EC REP Casus Europe B.V.Lange Viestraat 2b, Utrecht, 3511 BK, The Netherlands

0344 CH REP Casus Switzerland GmbH Hinterbergstrasse 49 6312 Steinhausen Switzerland

CH importer: Medline International Switzerland Sàrl, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK representative:
Qserve Group UK Ltd. 282 Farnborough Road, Farnborough, GU14 7NA, United Kingdom

Importer in UK:
Medline Industries Ltd, WA4 6HL, UK

Medline International B.V.,
Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

506-CT1318-1007 Rev.0 202411

REF 1642B REF 1651B REF 1659B REF 1077B

en Adaptors and Connectors

Intended use:

Aerosol Adaptors and Connectors are intended to be connected in-line with a nebulizer system.
Adaptors and Connectors are intended for use in both acute and non-acute care settings by qualified medical personnel, home healthcare providers, or patients under the general supervision of qualified medical personnel.
Adaptors and Connectors are intended for use on general populations, which may include pediatric and adult patients.
These devices are intended to be used for a maximum of 7 days.

Instructions for use:

REF **1642B, 1651B, 1077B**
1. Connect the adaptor or connector in-line. Ensure all connections are secure.

REF **1659B**
1. Connect the adaptor in-line.
2. Turn the white arrow in direction medication is to be administered.
3. Insert the canister and dispense medication as prescribed.

Disposal instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warnings:

For devices with arrows, ensure all arrows are pointing in the direction of flow.
For devices with ports, close off port with cap when not in use to prevent leakage.
These devices are not intended to be used in a pressurized or closed system.
Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
Non-Sterile. Do not sterilize prior to use.
Do not use if device is soiled or damaged.
Do not use in the presence of an ignition source.
If a serious incident has occurred in relation to the device, it should be reported to Medline/Besmed and the competent authority of the Member State in which the user and/or patient is established.

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