CHEMOTHERAPY INFORMATION

Surgical Gloves Not Made with Natural Rubber Latex.

en

Intended use
Medline Surgical Gloves are single use, disposable devices intended for surgical procedures that is worn on the surgical staff's hands to prevent cross contamination between patient and surgical team and to protect a surgical wound from contamination. Offered for sale to end users, such as hospitals, clinics, and surgeons.
Indications for User
The surgeon's glove is a disposable device made of synthetic rubber intended to be worn by operating room personnel to protect a surgical wound from contamination. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978 Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs.
Device supplied sterile using gamma radiation.
The following chemicals have been tested with these gloves. Refer to Table 1.
Chemotherapy Drug Permeation (Shortest Detection Time) in Minutes.

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en

Caution: Gloves used for protection against Chemotherapy Drugs exposure must be selected specifically for the type of chemotherapy agent used. Users are recommended to review Safety Data Sheets (SDS National Repository) for Chemotherapy Drugs being used to determine the required level of protection.
Warnings: Do not use if package is damaged. Use of damaged gloves may result in potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient.If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal: Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for chemotherapy hazardous waste.
MSG90xx, MSG15xx, MSG10xx, MSG92xx, MSG96xx: Warnings: Do not use with Carmustine and Thiotepa.
MSG90xx, MSG15xx, MSG10xx, MSG92xx, MSG96xx: Warnings: Do not use with Cyclophosphamide.
MSG63xx: Caution: Testing showed a minimum breakthrough time of 36.3 minutes with Carmustine and 49.2 minutes for Thiotepa.
MSG63xx: Caution: Testing showed a minimum breakthrough time of 38.5 minutes with Carmustine and 38.6 minutes for Thiotepa.
MSG99xx: Caution: Testing showed a minimum breakthrough time of 34.6 minutes for Thiotepa.
MSG99xx: Caution: Testing showed a minimum breakthrough time of 32.6 minutes with Thiotepa.

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