HUDSON RCI®

AQUAPAK®

REF 400040, 000-40F, 000-40B STERILE R

REF 400301, 400501, 400601 STERILE A 0°C

REF 400340, 003-40F, 400540, 400640, 006-40F, 003-40B STERILE R STERILE A 0°C

Rx ONLY

Teleflex MEDICAL
3015 Carrington Mill Blvd
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

CE 1639 EC REP
Teleflex Medical
IDA Business and Technology Park
Dublin Rd.,
Athlone, Co. Westmeath
Ireland

L02551 R01

© Copyright 2007, 2009, 2011, 2014, 2015, 2016, 2018 Teleflex

en

Prefilled Sterile Humidifier

Indication for use: Adds humidity to a patient's inspired breathing gases.

Instructions for use:

NOTE: Use only Hudson RCI AQUAPAK® sterile humidifier reservoirs.

1. Thread the Humidifier Adaptor onto the reservoir and tighten.
2. To open the reservoir, snap the outlet trigger upwards.
 CAUTION: DO NOT TWIST.
3. Connect the Humidifier Adaptor to a flow-meter.
 CAUTION: Do not over tighten wing nut.
4. Verify audible alarm function prior to use.
 Occlude port and listen for alarm to sound.
 CAUTION: An alarm from the pressure relief valve indicates the reservoir has accumulated too much pressure. Inspect the system for excessive blockage or a kinked tube.
5. Attach supply tubing to the reservoir outlet.
6. Turn on the flow-meter and verify the flow of gas through the device.

Warning: Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Contents of reservoir for inhalation therapy only, not for parenteral administration. Contains no antimicrobial agents.

Do not use if package is damaged.

AQUAPAK and Hudson RCI are registered trademarks of Teleflex Incorporated or its affiliates.

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