HUDSON RCI® AquaPak®

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en Prefilled Sterile Nebulizer

For inhalation, USP · Non-pyrogenic

Intended use:
Intended to add sterile water or saline solution in aerosol form to a patient's breathing gases. Medline Hudson RCI AQUAPAK Pre-Filled Nebulizers are intended for use in general population, which may include pediatric (ages 2 years and above) and adult. To be used with min 50PSI of oxygen as their main gas source. It is a single use device used by trained clinicians in clinical environments. This product is intended to be used until water/saline depletion or for a maximum of 14 days.

Instructions for Use:
NOTE: Use only Hudson RCI AQUAPAK Nebulizer Adaptors with AQUAPAK Nebulizer Reservoirs. For heated nebulization, use with only Hudson RCI AQUATHERM.

1. Thread the Nebulizer Adaptor onto the reservoir and tighten.
 CAUTION: Do not over-tighten the adaptor.
2. Grasp the return tube behind the thumb flanges on the puncture pin and remove the pin sheath.
3. Insert the pin into the puncture port on the reservoir and push until the thumb flanges rest against the reservoir.
4. Thread the nebulizer onto the flow-meter.
 CAUTION: Do not over-tighten the wing nut.
5. Attach the patient tube to the nebulizer outlet.
6. Turn the entrainment collar on the nebulizer to the desired oxygen concentration.
 NOTE: The following table depicts the estimated FIO2 output at each adaptor setting:

LPM Setting	Adaptor Setting %	Estimated FIO2 output %
5LPM	28%	24.6% to 31.4%
8LPM	33%	29.0% to 37.0%
10LPM	98%	≥ 86.2%

7. Set the flow-meter to the desired setting and verify the flow of gas through the device.
8. The nebulizer is ready for use when liquid flows through the return tube.

Warnings:
- Single Use Only. Do not reuse, reprocess, or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
- Contents of reservoir are for inhalation therapy only and are not intended for parenteral administration.
- Contains no microbial agents.
- Sterile. Internal contents sterile if bottle intact. Non-Pyrogenic. No bacteriostatic agent or other preservative added.
- Do not use in the presence of an ignition source.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal instructions:
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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CE 2797 EC REP Medline International France SAS,
Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia
Medline International B.V.,
Medline Industries Ltd,