MEDLINE

Instructions for Use

de Gebrauchsanweisung
fr Mode d'emploi
nl Gebruiksaanwijzing
it Istruzioni per l'uso
es Instrucciones de uso
pt Instruções para uso
sk Návod na použitie
cs Návod k použití

sv Instruktioner för användning
fi Käyttöohjeet
da Instruktion
no Bruksanvisning
el Οδηγίες χρησης
pl Instrukcja stosowania
sl Navodila za uporabo
ro Instructiuni de folosire

bg Инструкции за употреба
et Kasutusjuhised
hu Használati utasítás
lv Lietošanas instrukcijas
lt Naudojimo instrukcija
tr Upute za korištenje
hr Kullanim talimatlari
ar تعليمات الاستخدام

en

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and tests only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For correct use, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defects or imperfections. Discard if glove is damaged.

Approved for use with chemotherapeutic drugs per ASTM D6978. Please contact Medline for detailed information.

Refer to the guidelines provided by your healthcare facility for donning, doffing and hand hygiene when using the gloves.

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Restore Sense:
OATN32XS, OATN32S, OATN32M, OATN32L, OATN32XL

Mediguard Special Nitrile:
MGSNXS, MGSNS, MGSNM, MGSNL, MGSNXL

Mediguard Special Nitrile 100:
MGSN100XS, MGSN100S, MGSN100M, MGSN100L, MGSN100XL

Sensicare Extended Cuff:
MDSEXTXS, MDSEXTS, MDSEXTM, MDSEXTL, MDSEXTXL

Sensicare Ice:
486800, 486801, 486802, 486803, 486804

Curad 3G:
CUR8233, CUR8234, CUR8235, CUR8236, CUR8237

Accutouch:
MDS192073, MDS192074, MDS192075, MDS192076, MDS192077

Eco vinyl:
MGE323, MGE324, MGE325, MGE326, MGE327

Curad latex:
CURE8103, CURE8104, CURE8105, CURE8106, CURE8107

Sensicare Silk:
MDS9083, MDS9084, MDS9085, MDS9086, MDS9087

Sensicare Free:
MDSFREE-XS, MDSFREE-S, MDSFREE-M, MDSFREE-L, MDSFREE-XL

Sensicare Pro:
PRO400XS, PRO400S, PRO400M, PRO400L, PRO400XL

Generation Pink:
PINKE2500, PINKE2501, PINKE2502, PINKE2503, PINKE2504

Sensicare 300 long:
MDSMID-XS, MDSMID-S, MDSMID-M, MDSMID-L, MDSMID-XL

Mediguard White:
WMGSNXS, WMGSNS, WMGSNM, WMGSNL, WMGSNXL

Personal Protective Equipment of Complex Design Category III, in compliance with Regulation (EU) 2016/425, as amended to apply in GB, type tested to EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-2:2019, EN 16523-1:2015+A1:2018, EN ISO 374-4:2019, ISO 16604:2004, EN ISO 374-5:2016 & ASTM D6978-05:2023; Class 1 Medical Device in compliance with Regulation (EU) 2017/745 and UK Medical Device Regulation, type tested to EN455-1:2020+A2:2024, EN455-2:2024, EN455-3:2023 & EN 455-4:2009

Declaration of conformity/Konformitätserklärung/Déclaration de conformité/Conformiteitsverklaring/Dichiarazione di conformità/Declaración de conformidad/Declaração de conformidade/Vyhlásenie o zhode/Prohlášení o shodě/Försäkran om överensstämmelse/Vaatimustenmukaisuusvakuutus/Overensstemmelseserklæring/Konformitetserklæring/δήλωση συμμόρφωσης/Deklaracja zgodności/Izjava o skladnosti/Deklaratija de conformitate/декларация за съответствие/Vastavusdeklaratsioon/Megfelelőségi nyilatkozatot/Atbilstības deklarāciju/Atitikties deklaracija/Izjava o sukladnosti/Uygunluk beyanına/إعلان المطابقة/**www.medline.eu/document-library**

Medical Device class I under Regulation (EU) 2017/745,
PPE category III under Regulation (EU) 2016/425, certified by BSI Group The Netherlands B.V. 1066 EP, Amsterdam, Netherlands (Notified Body number 2797)

UK approved body:
BSI Assurance UK Ltd, MK5 8PP, UK
UK Representative and importer in UK:
Medline Industries Ltd, WA4 6HL, UK

Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH Rep and importer:**
Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

©2026 Medline Industries, LP
©2026 Medline, Curad, Sensicare and Restore are registered trademarks of Medline Industries, LP

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