Instructions for Use

- de Gebrauchsanweisung
- fr Mode d'emploi
- nl Gebruiksaanwijzing
- it Istruzioni per l'uso
- es Instrucciones de uso
- pt Instruções para Uso
- sk Návod na použitie
- cs Návod k použití
- sv Instruktioner för användning
- fi Käyttöohjeet
- da Instruktion
- no Bruksanvisning
- el Οδηγίες χρήσης
- pl Instrukcja obsługi
- sl Navodila za uporabo
- ro Instructiuni de folosire
- bg Инструкции за употреба
- et Kasutusjuhised
- hu Használati utasítás
- lv Norādījumi par lietošanu
- lt Naudojimosi instrukcijos
- hr Upute za uporabu
- tr Kullanim talimatlari
- ar تعليمات الاستخدام

REF SBA5CE
REF SBB5CE
REF SBP5CE
REF SBR5CE
REF SBW5CE
REF SBY5CE

MD STERILE EO

2 STERILIZE

LATEX

Latex free, Latexfrei, Sans latex, Latexvrij, Senza lattice, Sin látex, Sem látex, Bez latexu, Bez Latexu, Latexfri, Lateksiton, Latexfri, Lateksfri, Χωρίς λάτεξ, Bezlateksowe, Brez lateksa, Fără late, Несъдържащ латекс, Lateksivaba, Latex mentes, Nesatur lateksu, Sudētyje nėra latekso, Bez lateksa, Lateks içermez, خالي من اللاتكس

Medline International France SAS,
5 rue Charles Lindbergh, 44110 Châteaubriant, France.

CE 2797

CH REP

CH Rep and importer:
Medline International Switzerland SARL,
c/o MN & Associés SA, 1 Place de Longemalle,
1204 Genève, Switzerland

UK CA 0086

UK representative and importer in UK:
Medline Industries Ltd, WA4 6HL,UK

RA26SJP 2026-01

en Suture Boots

Intended use:
Designed to be placed over the tip of the instrument, suture boots protect, tag and help locate sutures during surgical procedures.

Intended user:
Intended to be used by healthcare professionals.

Instructions for use:

1. Select the appropriate size: pink boot for large instruments, or other colored boots for standard instruments.
2. Remove adhesive peel-off liner from cartridge, and place onto any flat surface within the sterile field.

3. Insert tip of instrument (artery forceps, vascular clamps, or similar surgical instruments) into suture boot.
4. Remove from cartridge and rinse with sterile water before use.
5. Once surgical procedure is done, place the used suture boot back into the cartridge for ease of counting.

Precaution:
Potential for complication includes but not limited to irritation. Single use only. Do not re-use. Re-use may cause a risk of infection and/or compromise functional reliability which could lead to patient harm. Do not use if the package is open or damaged.

Disposal:
Precaution should be taken when discarding this device and disposal of the device will be made in accordance with applicable hospital or national guidelines/ regulations for biologically hazardous waste.

Notice to Users and/or Patients:
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

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