de Gebrauchsanweisung
fr Notice d'utilisation
nl Gebruiksaanwijzing
it Istruzioni per l'uso
es Instrucciones de uso
pt Instruções para Uso
sk Návod na použitie
cs Návod k použití
sv Instruktioner för användning
fi Käyttöohjeet
da Instruktion
no Bruksanvisning
el Οδηγίες χρήσης
pl Instrukcja obsługi
sl Navodilo za uporabo
ro Instructiuni de folosire
bg Инструкции за употреба
et Kasutusjuhised
hu Használati utasítás
lt Naudojimosi instrukcijos
hr Upute za uporabu
ru Инструкция по применению
tr Kullanim talimatlari

Instructions for use

Medline sock size	EU size	Foot length (cm)
S	< 34	< 21.5
M	34 - 36	21.5 - 22.5
L	36 - 38	23.0 - 24.0
XL	39 - 43	24.5 - 27.5
XXL	44 - 46	28.0 - 29.5
XXXL / Bariatric	> 46	> 29.5
One size fits most	36 - 43	23.0 - 27.5

en The cuff must sit above the ankle bone and the tread must extend past the heel. The fit must be close enough to hold the slipper sock in position in order to minimise unnecessary movement that could cause friction. However, the sock must not compress the foot as this may impede circulation. Reassess sizing at regular intervals as the patient may require different sizes during their stay. The slipper socks should be changed on a daily basis to ensure maximum protection or when the anti-slip tread has been visibly damaged, or once soiled or wet. Disposal of the device shall be made in accordance with applicable hospital or national regulations for textiles. Not washable. Single use.

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PPE category II under Regulation (EU) 2016/425
Type tested to specification 2 PfG 2682: 2019-02 created by TÜV Rheinland LGA Products GmbH
(Notified body N°0197) located at Tillystraße 2, D - 90431 Nürnberg, Germany.

Declaration of conformity/Konformitätserklärung/Déclaration de conformité/Conformiteitsverklaring/Dichiarazione di conformità/Declaración de conformidad/Declaração de conformidade/ Vyhlásenie o zhode/Prohlášení o shodě/Försäkran om överensstämmelse/Vaatimustenmukaisuusvakuutus/Overensstemmelseserklæring/Konformitetserklæring/Δήλωση συμμόρφωσης/ Deklaracja zgodności/Izjava o skladnosti/Declarația de conformitate/Декларация за съответствие/Vastavusdeklaratsioon/Megfelelőségi nyilatkozatot/Atitikties deklaracija/ Izjava o sukladnosti/Декларация о соответствии/Uygunluk beyanı: **www.medline.eu/document-library**

Medline International France SAS,
5 rue Charles Lindbergh, 44110 Châteaubriant, France

PPE category II under Regulations (EU) 2016/425 on personal protective equipment, as amended to apply in GB
UK approved body (n°2571): TUV Rheinland UK Ltd., Friars Gate (Third Floor), 1011 Stratford Road, Shirley, Solihull, B90 4BN, UK
UK representative and importer in UK:
Medline Industries Ltd, WA4 6HL, UK