VAPREVENT® ORAL CARE

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en DISPOSABLE ORAL SWABS

1. With gloved hands, remove Dentips® from packaging. Ensure foam on Dentips® is intact.
2. Use of bite block/Berman airway strongly recommended.
3. Saturate foam with desired mouthwash. Gently clean oral cavity for approximately one minute.
4. After use, ensure foam is intact before discarding Dentips®. Inspect oral cavity and remove any particles, if present.

Warning: Do not use for tracheostomy care. Discontinue use and consult physician if irritation occurs. Keep out of the reach of children.

Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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©2020 Medline Industries, Inc. RJ20JYU - 10/2020
Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, Inc.
Packaged in USA for: **Medline Industries, Inc.,** Three Lakes Drive, Northfield, IL 60093 USA.
EC REP **Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France**