Impregnated dentifrice, ingredients / Imprägniertes zahnputzmittel, inhaltsstoffe / Dentifrice complet, ingrédients / Geïmpregneerd product voor mondreiniging, ingrediënten / Dentifricio imbevuto, ingredienti / Dentífrico impregnado, ingredientes / Dentífrico impregnado, ingredientes / Impregnovaný zubný prípravok, zložky / Impregnovaná zubní pasta, složení / Impregnerad tandkram, ingredienser / Impregnoitu hammastahna, ainesosat / Imprægneret tandpasta, ingredienser / Mettet tannkrem, innholdsstoffer / Εμποτισμενο με προϊον καθαρισμου δοντιων, συστατικα / Impregnowane środkiem do mycia zębów, składniki / Impregnirana sredstva za čiščenje zob, sestavine / Pastă de dinți impregnată, ingrediente / Импрегниран дентифрис, съставки / Immutatud hambapulber, koostisosad / Fogkrémmel impregnálva, összetevők / Impregnuota dantų pasta, sudedamosios dalys / Impregnirana pasta za zube, sastojci / Emprenye edilmiş diş temizleme maddesi, içindekiler

Water (Aqua), Calcium Carbonate, Glycerin, Peppermint flavor, Sodium Saccharin, Carboxymethylcellulose, Sodium Tri-poly phosphate, Xanthan Gum, Sodium benzoate, Methyl Paraben, Green 5

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Packaged in Mexico for: **Medline Industries, LP,** Three Lakes Drive, Northfield, IL 60093 USA

EC REP Medline International France SAS,
5 rue Charles Lindbergh, 44110
Châteaubriant, France

UK representative and importer in UK:
Medline Industries Ltd, WA4 6HL,UK

CH REP **CH Rep and importer:** Medline International
Switzerland SARL, c/o MN & Associés SA, 1 Place de
Longemalle, 1204 Genève, Switzerland

Medline International B.V.,
Catharijnesingel 47, 3511 GC
Utrecht, The Netherlands

Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.
©2025 Medline Industries, LP

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en Disposable Oral Swabs

Oral Care products are intended to clean oral secretions from the oral cavity during oral care procedures. They can be used alone or in combination with Oral Care solutions to clean the oral cavity.

1. With gloved hands, remove Dentips® from packaging. Ensure foam on Dentips® is intact.
2. Use of bite block/Berman airway strongly recommended.
3. Saturate foam with desired mouthwash. Gently clean oral cavity for approximately one minute.
4. After use, ensure foam is intact before discarding Dentips®. Inspect oral cavity and remove any particles, if present.

Warning: Do not use for tracheostomy care. Discontinue use and consult physician if irritation occurs. Keep out of the reach of children. Protect from freezing. Avoid excessive heat.

Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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