VAPREVENT® ORAL CARE

REF MDS096525EU

en Suction Swabs Individually Wrapped

Intended Use:

Oral Care products are intended to clean oral secretions from the oral cavity, suction swabs and toothbrush also are intended to aspirate oral secretions, during oral care procedures. These devices are not intended to be used by lay users or with freestanding suction vacuum equipment. They can be used alone or in combination with Oral Care solutions to clean the oral cavity.

1. With gloved hands, remove suction swab from packaging. Ensure foam on suction swab is intact. Connect Suction swab to standard suction line. Ensure connection is secure.
2. Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
3. Saturate suction swab with water or desired solution.
4. Gently clean teeth and gums for approximately 1-2 minutes. Provide suction by using the thumb port.
5. Ensure foam is intact before discarding suction swab. Inspect oral cavity and remove any particles, if present.

Warning: Discontinue use and consult physician if irritation occurs. Keep out of the reach of children.

Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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Packaged in Mexico for: **Medline Industries, LP,** Three Lakes Drive, Northfield, IL 60093 USA

2797

EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH Rep and importer:** Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK CA 0086 **UK representative and importer in UK:** Medline Industries ltd, WA4 6HL, UK

Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.
©2024 Medline Industries, LP