297 mm

210 mm

ELECTROSURGICAL PENCIL

EN INSTRUCTIONS FOR USE
DE GEBRAUCHSANWEISUNG
IT NOTICE D'UTILISATION
ES MANUAL DE INSTRUCCIONES
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR KULLANIM TALİMATLARI
PT INSTRUÇÕES PARA USO
NL INSTRUCTIEBOEKJE VOOR GEBRUIK
SV ANVÄNDARINSTRUKTION
DA BRUGSANVISNING
NO BRUKSANVISNING
FI KÄYTTÖOHJEET
IS NOTKUNARLEIÐBEININGAR
RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
PL INSTRUKCJA OBSŁUGI
CS NÁVOD K POUŽITÍ
SK NÁVOD NA POUŽITIE
HU HASZNÁLATI UTASÍTÁS
RO INSTRUCȚIUNI DE UTILIZARE
SL NAVODILA ZA UPORABO
HR UPUTE ZA UPORABU
SR УПУТСТВО ЗА УПОТРЕБУ
BG ИНСТРУКЦИИ ЗА УПОТРЕБА
UK ІНСТРУКЦІЯ З ЗАСТОСУВАННЯ
ET KASUTUSJUHEND
LV LIETOŠANAS INSTRUKCIJA
LT NAUDOJIMO INSTRUKCIJOS
MT ISTRUZZJONIJIET GĦALL-UŻU
GA TREORACHA DON ÚSÁID
SQ UDHËZIME PËR PËRDORIM
MK УПАТСТВО ЗА УПОТРЕБА
AR تعليمات الاستخدام
FA دستورالعمل استفاده
HE הוראות שימוש
ZH 使用说明书
JA 使用説明書
KO 사용 설명서
TH คู่มือการใช้งาน
VI HƯỚNG DẪN SỬ DỤNG
ID PETUNJUK PENGGUNAAN
MS ARAHAN PENGGUNAAN
HI उपयोग के निर्देश

REF ESP0300LQ
REF ESP0300LQ
REF ESP0400LQ
REF ESP0500LQ
REF ESP0400DLQ
REF ESP0400LQ
REF ESP0400SLQ
REF ESP0500LQ
REF ESP0400DLQ
REF ESP0400LQ
REF ESP0400SLQ
REF ESP0500LQ

Do not use if the package is opened or damaged
Store between +5°C and +40°C in a dry place
Keep away from sunlight
Keep away from heat
Keep away from rain
Do not reuse
Do not resterilize
Do not use if the package is damaged
Consult instructions for use
Non-sterile
Sterilized using ethylene oxide
Sterilized using irradiation
Sterilized using steam or dry heat
Manufacturer
Date of manufacture
Use by date
Batch code
Catalogue number
Serial number
Quantity
Single sterile barrier system
Single sterile barrier system with protective packaging inside
Double sterile barrier system
Single sterile barrier system with protective packaging outside
Double sterile barrier system with protective packaging inside
Double sterile barrier system with protective packaging outside
Caution
Medical device
Authorized representative in the European Community
Unique device identifier
Do not use if package is damaged and consult instructions for use
Do not use if sterile barrier is damaged
Do not use if sterile barrier is damaged and consult instructions for use
Keep dry
Keep away from sunlight
Keep away from heat
Keep away from rain
Do not reuse
Do not resterilize
Do not use if the package is damaged
Consult instructions for use
Non-sterile
Sterilized using ethylene oxide
Sterilized using irradiation
Sterilized using steam or dry heat
Manufacturer
Date of manufacture
Use by date
Batch code
Catalogue number
Serial number
Quantity
Single sterile barrier system
Single sterile barrier system with protective packaging inside
Double sterile barrier system
Single sterile barrier system with protective packaging outside
Double sterile barrier system with protective packaging inside
Double sterile barrier system with protective packaging outside
Caution
Medical device
Authorized representative in the European Community
Unique device identifier

STERILE R

Medical Device

CE 0482

Dongguan Gaomifeid Equipment Mfg Co., Ltd.
Building 1, Dongguan Branch of Gaomi Hi-tech Cuck Management Zone, Keqiang Town, Dongguan, Guangdong, China
Tel:(0769)82000123 Fax:(0769)82000183

EC REP
Medical Device Borkstrasse 10 48163 Muenster, Germany
Tel.: (49)251 32266-0 Fax: (49)251 32266-22

CH REP
Medfair SWISS Gmbh Bäckerstrasse 15, 5400 Baden, Switzerland

Medline International B.V., Carlserlaangen 49 2011 DC Utrecht, The Netherlands

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

2026-02

EN ELECTROSURGICAL PENCIL

The Electrosurgical pencils are intended for monopolar cutting and coagulation with various types of RF electrosurgical generators. They are designed to accept standard electrodes with a 2.4mm (3/32") diameter shaft.

HOW TO USE
- Use aseptic techniques when you remove the product from its packaging.
- Remove tip protection from electrode. Check that the electrode is firmly sealed in the hand piece. If you change the electrode, make sure that it is inserted fully into the hand piece.
- Plug in to the RF generator.
- The design and construction of this product does not allow its reuse.
- Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.

SAFETY INSTRUCTION WARNING
- Sterility is only guaranteed unless package is damaged or seal is broken. DO NOT RESTERILIZE.
- Do not use after expiration date indicated on the packaging.
- Broken cable may cause arcing or patient burns.
- Check instrument before use. Do not use damaged pencils or electrodes.
- Do not use in the presence of flammable or explosive materials.
- Keep the electrode clean. Do not use electrodes not in use.
- The design and construction of this product does not allow its reuse.
- Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.
- Medline explicitly warns against making any changes to the instrument. Any such change will void Medline from all liability.
- Always place the active Electrode in a clean, dry, insulated safety holster when not in use.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place away from the patient and drapes.
- Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.

ELECTRICAL CAPACITY: The rated voltage is 5kV.

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