297 mm

210 mm

297 mm

ELECTROSURGICAL PENCIL

EN INSTRUCTIONS FOR USE
DA GEBRAUCHSANWEISUNG
FI NOTICE D'UTILISATION
IT GEBRUIKSAANWIJZING
IS ISTRUZIONI PER L'USO
ES INSTRUCCIONES DE USO
PT INSTRUÇÕES PARA USO
SV INSTRUKTIONER FÖR ANVÄNDNING
EL ΚΑΫΤΟΟΗΙΟΕΙ
CS INSTRUKTION
NO BRUKSANVISNING
SK ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
PL INSTRUKCJA OBSŁUGI

REF ESP83001Q	REF ESP84001LQ
REF ESP84000Q	REF ESP84001ERBQ
REF ESP84001Q	REF ESRK4001Q
REF ESP84000LQ	REF ESRK4001LQ
REF ESP84002Q	

EN Do not use if the package is opened or damaged
DA Nicht benutzen, wenn die Packung offen oder beschädigt ist
FI Ne pas utiliser si l'emballage est ouvert ou endommagé
IS Niet gebruiken indien verpakking geopend of beschadigd is
IT Non utilizzare se la confezione non e' integra
ES No utilizar si el paquete está abierto o dañado
PT Não usar se a embalagem estiver aberta ou danificada
SV Används ej om förpackningen är öppnad eller skadad
EL Μη χρησιμοποιείτε εάν η συσκευασία είναι ανοιχτή ή κατεστραμμένη
CS Nepoužívejte, jestliže je obal otevřený nebo poškozený

STERILE R

Dongguan Queen Med Equipment Mfg Co., Ltd.,
Building 1, Dongyuan Industrial Centre, Haibian Village,
Caojie Management Zone, Xiagang Town, Dongguan,
Guangdong, China
TEL:86-769-82023333 Fax:86-769-82023111

Medical Device Limited
Unit3, Abbey House, Shannon,
Town Centre,County Clare,Ireland
TEL:(353)61 771888Fax:(353)61 771944

0482 EC REP

201×26

PMS COLORS USED CMYK
Black

ELECTROSURGICAL PENCIL

The Electrosurgical pencils are intended for monopolar cutting and coagulation with various types of RF electrosurgical generators. They are designed to accept standard electrodes with a 2.4 mm (3/32") diameter shaft.

HOW TO USE
- Use aseptic techniques when you remove the product from its packaging.
- Remove tip protection from electrode. Check that the electrode is firmly sealed in the hand piece. If you change the electrode, make sure that it is firmly sealed in the hand piece.
- Plug in to the RF generator.
- The design and construction of this product does not allow its reuse.
- Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.

SAFETY INSTRUCTION-WARNING
- Sterility is only guaranteed unless package is damaged or seal is broken. DO NOT RESTERILIZE!
- Do not use after expiration date indicated on the packaging.
- Broken cable may cause arcing or patient burns.
- Check instrument before use. Do not use damaged pencils or electrodes.
- Do not use in the presence of flammable or explosive materials.
- Place the instrument away from the patient when not in use.
- The design and construction of this product does not allow its reuse.
- Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.
- Medline explicitly warns against making any changes to the instrument. Any such change will exempt Medline from all liability.
- Always place the active Electrode in a clean, dry, insulated safety holster when not in use.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place away from the patient and drapes.
- Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.

ELECTRICAL CAPACITY: The rated voltage is 5kVp.

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