ELECTROSURGICAL ELECTRODES

Read Instructions before use.

Do Not Reuse. This is a single use device.

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REF:

Item code	Description
WN01-FG003001	PTFE, Blade electrode, extended insulation, 7 cm
WN01-FG000602	PTFE, Blade electrode, extended insulation, 10 cm
MN01-FG001401	PTFE, Blade electrode, 7 cm
WN02-FG002202	PTFE, Needle electrode, 15 cm
WN01-FG003201	PTFE, Blade electrode, 15 cm
MN01-FG001501	PTFE, Needle electrode, 7 cm

en INSTRUCTIONS FOR USE

ELECTROSURGICAL ELECTRODE

The Electrosurgical electrode are intended for monopolar cutting and coagulation with various types of RF electrosurgical generators.

HOW TO USE:

.Use aseptic techniques when you remove the product from its packaging.
.Remove tip protection from electrode. Check that the electrode is firmly sealed in the hand piece. If you change the electrode, make sure that it is inserted fully into the hand piece.
.Plug in to the RF generator.
.The design and construction of this product does not allow its reuse.
.Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.

SAFETY INSTRUCTION WARNING:

. Sterility is only guaranteed unless package is damaged or seal is broken. DO NOY RESTERLIZE
. Do not use after expiration date indicated on the packaging.
. Check instrument before use. Do not use damaged pencils or electrodes.
. Do not use in the presence of flammable or explosive materials.
. Place the instrument away from the patient when not in use.
. The design and construction of this product does not allow its reuse.
. Always refer to the instruction manual of the RF generator and to its general information about electrosurgical procedures.
. Medline explicitly warns against making any changes to the instrument. Any such change will exempt Medline from all liability.
. we recommend to always place the active Electrode in a clean, dry, insulated safety holster when not in use.
. Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes.

Place away from the patient and drapes.

. Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
. Improper electrode installation may result in injury to the patient or operating room personnel.
. Refer to generator's manual of instruction for power setting recommendations for various electrosurgical electrode configurations and specifications on voltage output characteristics to ensure that all safety precautions are followed.
. Confirm proper electrosurgical generator setting before proceeding with surgery. Power setting guidelines may vary due to differences in surgical technique, patient, electrodes and surgical set-up.
. Before inserting or changing an electrode be sure that the active accessory is not connected to an electrosurgical generator or make sure that the generator is off (standby).

ELECTRICAL CAPACITY: The rated voltage is 3kVp.

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CMC Medical Devices GmbH
Luxus Lebenswelt GmbH Kochstr.1,47877,Willich,Germany
Ltd.3689# , East Yangguang Avenue,Dipu
Medline International B.V., Nieuwe Stationsstraat 10,