

CLOSED WOUND DRAINAGE SYSTEMS

en **INSTRUCTIONS FOR USE** (p. 2–6)
de **GEBRAUCHSANWEISUNG** (p. 7–9)
fr **NOTICE D'UTILISATION** (p. 10–12)
nl **GEBRUIKSAANWIJZING** (p. 13–15)
it **ISTRUZIONI PER L'USO** (p. 16–18)
es **INSTRUCCIONES DE USO** (p. 19–21)
pt **INSTRUÇÕES PARA USO** (p. 22–24)
sv **INSTRUKTIONER FÖR ANVÄNDNING** (p. 25–27)
fi **KÄYTTÖOHJEET** (p. 28–30)
da **INSTRUKTION** (p. 31–33)
no **BRUKSANVISNING** (p. 34–36)
ja 取扱説明書 (37〜39ページ)

ExuFlow ExuFlow

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Single Use Only. 再使用禁止

Do not (re) sterilize.

STERILE EO Sterilized by Ethylene Oxide.

MD Medical Device

Single sterile barrier system with protective packaging inside.

Caution

Caution: Mandatory action to read the instructions for use.

Keep Dry

Keep away from Sunlight

℞ Only Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Storage: Protect from freezing. Avoid excessive heat.

Manufactured in China for: Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093 USA.

CE 2797 EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

UK CA 0086 **UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL, UK








CH REP **CH Rep and importer:** Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia.

Made in China - Hergestellt in China - Fabriqué en Chine - Geproduceerd in China - Fabbricato in Cina - Fabricado en China - Fabricado na China - Tillverkad i Kina - Valmistettu Kiinassa - Fremstillet i Kina - Laget i Kina - 中国製

ExuFlow is a trademark and Medline is a registered trademark of Medline Industries, LP.
V2 RE24NID 2024-05

				
	REF DYNJWE1320 DYNJWE1321A DYNJWE1323A DYNJWE1325A	REF DYNJWE0321 DYNJWE0323 DYNJWE0325		REF DYNJWE1010 (400 cc/mL)
	REF DYNJWE1308A DYNJWE1309 DYNJWE1310 DYNJWE1311	REF DYNJWE1410 DYNJWE1411 DYNJWE1409	REF DYNJWE1348 (100 cc/mL) DYNJWE1349 (100 cc/mL) DYNJWE1360 (100 cc/mL) DYNJWE1361 (100 cc/mL)	REF DYNJWE1481 (100 cc/mL)
	REF DYNJWE2186 DYNJWE2188 DYNJWE2190 DYNJWE2234 DYNJWE2186NH	REF DYNJWE2187 DYNJWE2189 DYNJWE2191 DYNJWE2187NH		REF DYNJWE2264 (100 cc/mL) DYNJWE2265 (100 cc/mL)
	REF DYNJWE2211 DYNJWE2213 DYNJWE2214	REF DYNJWE2212 DYNJWE2215		
 ExuFlow	REF ORHUR100 ORHUR150 ORHUR190 ORHUR860 ORHUR880 ORHUR900	REF ORHUR101 ORHUR151 ORHUR195 ORHUR870 ORHUR890 ORHUR910		
 ExuFlow	REF ORHUF100 ORHUF070	REF ORHUF071 ORHUF074 ORHUF101 ORHUF104		

en CLOSED WOUND DRAINAGE SYSTEMS

Description
Medline's silicone closed wound drainage systems are single use sterile products. The closed wound drainage system consists of two primary components: a surgical wound drain and a fluid collection reservoir. The wound drains are comprised of clear silicone tubing and are manufactured to be radiopaque. Medline's wound drains come in three styles: perforated, channeled, and ExuFlow. They also include optional trocar and adaptor offerings. The wound drain may be surgically placed into the surgical wound site using the optional pre-attached trocar. After wound drain placement, the trocar should be removed and disposed of. Exudate is collected through a wound drain and into the fluid collection reservoir.

Indications/Intended Purpose
Medline Closed Wound Drainage System devices are used to remove body fluid/exudate from the patient's body. The Wound Drain is used to drain fluids/exudate from body cavities. The Reservoir is intended to provide a vacuum for suction drainage of surgical wounds. The (optional) Trocar is used to access the suction site. Medline's closed wound drainage system is commonly used post-operatively for collecting exudate from surgical sites.

Intended for use on patients with surgical wounds requiring wound drainage, including adult and geriatric patients

Intended to be used by healthcare professionals in a surgical/operating room, procedural room, acute care, or post-acute care environment. Patients sent home with the device in place may empty and reactivate the reservoir, as directed by a healthcare professional.

Warnings/Precautions
1. Medline's closed wound drainage system is not a measuring device and is not intended to be used for accurate fluid measurements.
2. Repeated application of Medline's closed wound drainage system is not recommended.
3. Medline's closed wound drainage system cannot be used longer than 30 days.
4. Medline's closed wound drainage system and its components are for single use only and should be discarded after use.
5. Wound drains should not be handled with any instruments. Sharp objects can lead to tearing, weakening, or breakage of the drain.
6. In order to allow for easy removal, the wound drain's tubing should not be curled, stretched or sutured into place.
7. Wound drains should be placed and removed with care, avoiding fast movements, while applying slow and steady pressure.
8. Utilize sterile technique when making wound drain insertions and connections.
9. The product can be used with medical grade adhesive tape.
10. During wound drain placement, take care to not cut, tear or damage the wound drain, as this can lead to tearing or complete breakage.
11. Careful attention should be directed to make sure the drains do not clog or become occluded when they are still in use. If occlusion occurs, squeeze and pump the reservoir to generate suction and relieve the occlusion.
12. The reservoir should not be allowed to fill completely, this can lead to increased likelihood of back-flow contamination.
13. All drainage holes, channels, or ducts must be located within the wound site, creating an air-tight seal to maintain vacuum. An air-tight seal must be maintained to prevent the reservoir from filling with air, which can increase the likelihood of backflow.
14. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.
15. Patients sent home with the device should consult a healthcare professional if they experience inadequate/failure to drain exudate/bodily fluids, the Wound Drain becomes dislodged, occluded or damaged, or exudate/bodily fluid backflows to the wound site.

Contraindications
1. Do not use with any other closed wound drainage system or components on the market.
2. Do not reuse Medline's closed wound drainage system.
3. Do not use Medline closed wound drainage system for longer than 30 days.
4. Do not suture through, cut, nick, tear or damage the drain, as this may lead to breakage.
5. Do not use with any other suction vacuum source other than Medline's 100 cc/mL, 200 cc/mL, or 400 cc/mL reservoir.
6. Do not use if package has been opened or if the seal has been compromised.
7. Do not resterilize the product.

Risks/Side Effects
Possible risks include infection and inflammation. Prior to the surgery, the patient should be informed of such risks and complications that may occur, depending on the patients degree of intolerance to any object placed into the body.


A drain left in a patient for any period of time allows for tissue ingrowth around and into the drain, which can cause breakage and discomfort upon removal.

Reservoir Initial Testing
For Medline's reservoir, please follow the instructions / diagram below.
1. Remove Reservoir (D) from packaging and test for efficacy as follows:
 a. Fully compress Reservoir (D).
 b. Close Disposal Port (C) by tightly securing Screw Cap (A) to Disposal Port (C).
 c. Occlude Patient Port (B) and release Reservoir (D).
 d. If vacuum is effective, Reservoir (D) should remain compressed. If Reservoir (D) expands, it should be discarded.


Drain Connection/Drainage Placement
Place drain in wound site, ensuring that all holes, channels, or ducts are located within the wound site.
On drains marked with indicators, pull completely through the insertion site until the initial set of indicators appear at the skin surface.
1. Attach the drain or adaptor to the reservoir via Patient Port (B).
2. For dual drain connections (400 cc/mL bulb only):
 a. Cut tip from Patient Port 2 (E).
 b. Attach second drain or adapter to Patient Port 2 (E).
3. Firmly compress Reservoir (D) to evacuate air.
4. Close Disposal Port (C) with Screw Cap (A).
5. Release Reservoir (D) – the closed wound drainage system is now active.

Reservoir should be located at patient level or below drain exit for adequate suction to be maintained.
Drainage should be monitored such that Reservoir is not filled to excess.


Silicone Reservoir
100 cc/mL
DYNJWE1305

Silicone Reservoir
200 cc/mL
DYNJWE2000

Silicone Reservoir
400 cc/mL
DYNJWE1000

Storage/Disposal
Protect from freezing, avoid excessive heat. Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Risks Associated with Reuse
Medline's closed wound drainage system and its components are not validated for reuse. The design of the device does not facilitate cleaning and can result in an increased risk of infection if not properly decontaminated.
The device may not perform as intended by the manufacturer if reused, resterilized, or decontaminated. Reuse of the product may result in material degradation resulting in a device failure. Reuse of single use devices creates potential risk for the patient or user. It may cause contamination and/or impairment of function that may lead to injury, illness and/or death of a patient.

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