ExuFlow® Drain

Instructions for Use (p. 3–4)

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



en ExuFlow® Drain

Description
ExuFlow Drain is a single use sterile products. ExuFlow perforated channel drains combine a structure with multiple narrow ducts and internal perforations and are manufactured to be radio-opaque. The Drain is surgically placed into the surgical wound site using a trocar. Exudate is collected through a wound drain into a fluid collection Reservoir.

Intended Use
To withdraw fluids/gases from the wound or surgical area.

Indications for Use
The ExuFlow Drain are intended to be used in clinical cases which require the removal of fluids and exudates accumulated in body cavities or wounds following surgery.

Intended User
ExuFlow Intended to be used by healthcare Professionals only (physicians or trained Nurses) in hospital settings such as operating rooms, intensive care units (ICUs), and emergency departments.

Patient Population
Depending on the size of the tube, this device can be used in adult, pediatric and neonatal patients undergoing surgical procedures requiring drainage.

Warnings/Precautions
1. Drain, and Trocars (optional) are for single use only and should be discarded after use.
2. Drain should not be handled with any instruments. Sharp objects can lead to tearing or weakening and breakage of the drain.
3. Do not use if package has been opened or if the seal has been compromised.
4. In order to allow for easy removal, the drain should not be curled, stretched or sutured into place. Drains should be placed and removed with care avoiding fast movements applying slow and steady pressure.
5. Use an aseptic technique during placement and handling
6. The product can be used with medical grade adhesive tape.
7. During drain placement, take care to not cut, tear or damage the drain as this can lead to tearing or complete breakage.
8. Careful attention should be directed to make sure the drains do not clog or become occluded when they are still in use.
9. Regularly inspect the drain for clogging or occlusion; in case of obstruction, the drain can be aspirated by connecting suction to the reservoir outlet or temporarily disconnecting the drain from the reservoir and applying suction directly to the drain.
10. All drainage holes, channels, or ducts must be located within the wound site creating an air tight seal to maintain vacuum. An air tight seal must be maintained to prevent the reservoir from filling with air which can increase the likelihood of backflow.
11. Do not resterilize the product.

Contraindications
The Surgical Drains have no absolute contraindication.

Possible Complications
Possible risks include infection and inflammation. Prior to the surgery, the patient should be informed of such risks and complications that may occur depending on the patients degree of intolerance to any object placed into the body. A drain left in a patient for any period of time allows for tissue ingrowth around and into the drain which can cause breakage and discomfort upon removal.

Drain Connection/Drainage Placement
Place drain in wound site, ensuring that all holes, channels, or ducts are located within the wound site.
Connect drain to the drainage system, according to the methods indicated in the respective instruction sheet.

Storage/Disposal
Protect from freezing, avoid excessive heat.
Disposal should be handled according to local environmental procedures.

Risks Associated with Reuse
The wound drainage system devices are not validated for reuse. The design of the device does not facilitate cleaning and can result in an increased risk of infection if not properly decontaminated.
The device may not perform as intended by the manufacturer if reused, resterilized, or decontaminated. Reuse of the product may result in material degradation resulting in a device failure.

CONTACT
DEGANIA SILICONE Ltd.
Any serious incident that has occurred in conjunction to this device, must be reported to Degania Silicone Ltd (at the address listed below) and to the relevant regulatory authority and/or competent authority of the Member State in which the user/patient is established.

Glossary	
	Manufactured by Degania Silicone LTD., Degania Bet 1513000, Israel
EU REP	Authorized - Representative in the European Community/ European Union. Biometrix S.R.O. Vicenzy 16, Samorin 93101, Slovakia
CH REP	Swiss authorised representative MedNet Swiss GmbH, D4 Platz 4, 6039 Root D4, Switzerland
	Importer Biometrix S.R.O. Vicenzy 16, Samorin 93101, Slovakia
Rx Only	Caution: This device is for sale by or on order of a physician
LATEX	Not made with natural rubber latex
PHT DEHP	Not made with DEHP
BPA	Not made with BPA
MR	This device is MR Safe
MD STERILE EO REF LOT	

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Produsert av Degania Silicone LTD., Degania Bet 1513000, Israel
CH REP MedNet Swiss GmbH
Importer in UK: Medline Industries Ltd, WA4 6HL, UK