Equipment Drape - Fluoroscope
Hybrido™ Adjustable Clip Drape for C-Arm
Instructions for Use

MEDLINE MICROTEK

Adjustable Clip Drape for C-Arm
REF 17685

Adjustable Clip Drape Complete Set for C-Arm
REF 17686
REF 17687

Intended Use / Indications for Use
Equipment Drapes are intended to cover/wrap various surgical and/or nonsurgical instruments/equipment [e.g. robotic arms, microscopes, tables, x-ray systems, light handles, etc]. Equipment Drapes are used to establish a barrier that minimizes the passage of microorganisms, body fluids, and other contaminants. They are typically made of natural and/or synthetic flexible materials and shaped to fit the type of instrument/equipment. This is a single-use device [not to be reprocessed].

Contraindications
There are no known contraindications for these devices.

Warnings
- **Visually inspect for breaches of packaging integrity (such as tears, cuts, open seals, etc.)**
- Contents are supplied sterile using an ethylene oxide (EO) process. Do not use if sterile barrier package is damaged. If damage is found, call your representative.
- Follow aseptic technique during the draping process. If at any time during the draping process the drape becomes damaged, remove the drape and dispose of it. Use a new drape and restart the draping process.
- The drapes are single-use devices. Do not reuse or re-sterilize single-use devices. Reuse of single use devices creates a risk of patient or operator infection.
- Caution: Ensure film edge is not between the clips and the equipment.

Precautions
There are no known precautions for this device.


Product / Device Description
The drape is composed of synthetic materials and is intended to protect the C-Arm from contamination during various surgical procedures.


How Supplied
- Do not use the drape if the sterile barrier package is damaged.
- This product may be supplied with rubber bands and/or ties.


Report serious incidents to legal manufacturer/importer and competent authority.


Reuse Warning Statement
For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.



Handling and Storage
- After use, dispose of product and packaging in accordance with hospital,administrative and/or local government policy.

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©2024 Medline Industries, LP. All rights reserved.
Microtek Medical Malta Ltd.