Rotation Anesthesia Mask, Valve

Instructions for use

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en Rotation Anesthesia Mask, Valve

Model & Specification

No.	Spec.	Size	Model
1	AMR53020H	0#	Model A
2	AMR53021H	1#	
3	AMR53022H	2#	
4	AMR53023HH	3#	
5	AMR53024HH	4#	
6	AMR53025HH	5#	
7	AMR53026HH	6#	

Parts Description

As is shown in Fig.1 (Model A), the Disposable Anesthetic Mask is consisted of shell, air cushion, inflation port and retaining ring.


Fig.1

Product Description

During clinical use, the clinician should first unpack the mask, and connect the connector of the mask to the other one of breathing circuit that is connected with respirator, ventilator or anesthetic machine, then open the machine and let the gas pass through the circuit and the mask, finally be absorbed by patient.

The device itself does not include the breathing circuit and the anesthetic gas, it is put into market in a non-sterile state for single use

Performance

1. Inflation and deflation
 For mask with inflation port, the inflation port shall allow users to inflate and deflate the air cushion using syringe without needle.
2. Security of the retaining ring attachment
 Overhang the anesthetic mask with conical connector upwards, the retaining ring and mask shall not detach.

Intended Use

The Disposable Anesthetic Mask is used as a passage for anesthetic gas to enter into the patient body, which is usually used with simple respirator, ventilator or anesthetic machine in medical institution.

Indications

Used in clinical anesthesia to introduce anesthetic gas into patients' upper airway.

Intended users

Professional medical staff who has been trained by anesthesia and first-aid knowledge.

Intended patient group

No.	Size	Intended patient group
1	0#	Neonate
2	1#	Infant
3	2#	Pediatric
4	3#	Adult (S)
5	4#	Adult (M)
6	5#	Adult (L)
7	6#	Adult (XL)

Contraindication

Patients who are allergic to PVC shall use with caution.

⚠ **Attention:**

1. Excessive pressure applied to face mask can lead to facial and/or optic nerve damage.
2. The device is designed to mate with standard 15 mm/ 22 mm connector.
3. Disposal of Disposable Anesthetic Masks shall be in accordance with local regulations.
4. Refrain from pressing rigid parts of the mask against the face, nose or eyes.

⚠ **Warning**

1. Do not use if device or package is damaged.
2. Do not use if device has past its expiry date .
3. Do not use in the presence of flammable anesthetics.
4. The device is for single use and reuse is forbidden, or it may cause cross infection.
5. Only to be used by medical professionals with the required knowledge and training.
6. The device shall not be contacted with oil, or it may cause contamination and discomfort .
7. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of Member State in which the user and/ or patient is established.
8. After use, the device should be isolated and managed as waste, not stored at will. It should be destroyed and disposed in accordance with the local laws and regulations. Do not discard at will to avoid infection or microbial hazards
9. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of Member State in which the user and/ or patient is established.

Operation Method

Step 1. Unpack the mask and check air cushion for proper inflation pressure.

Step 2. Connect mask to a breathing circuit that is connected with respirator, ventilator or anesthetic machine.

Step 3. Be sure all connections are secured, then open the machine and let gas pass through the circuit and mask, finally be absorbed.

Step 4. If there's a concern of mask failure, please replace it immediately.

Placement

1. Place the air cushion on the patient's chin first.
2. Roll air cushion up the face with gentle, but firm pressure.

Storage and Transportation Condition

1. The device should be stored in a non-corrosive gas, relative humidity NMT 80% and clean environment.
2. Avoid hard pressure, direct sunlight and rain-snow exposure.

Packaging marks

Symbol	Meaning	Symbol	Meaning
	Country and date of manufacture	LOT	Batch Code
	Do not reuse		Keep away from sunlight
	Caution		Consult instructions for use
	Manufacturer		Use by
EU REP	Authorised representative in the european community		Do not use if the package is damaged
CE 0197	The product has passed CE certification		One-way valve
	Keep dry	MD	Medical device
	Fluid path	REF	Model number
A 文	Translation	%	Humidity limitation
	Temperature limit	UDI	Unique device identifier

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Ningbo Xinwell Medical Technology Co., Ltd.,
No.188 Binjiang Road, Cixi High-tech Industrial Development Zone, Cixi City, 315301 Zhejiang P.R. China

CE 0197

EU REP Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

CH REP SUNGo Technical Service GmbH Bahnhofstrasse 21,6300 Zug, Switzerland

CH importer: Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK representative: Lotus Global Co Ltd, 23 Maine Street, Reading, RG2 6AG, England, United Kingdom
Importer in UK: Medline Industries Ltd, WA4 6HL, UK

Medline International B.V.,
Catharijnesingel 47, 3511 GC Utrecht, The Netherlands

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia