Instruction for use

SOLUNUM SİSTEMİ KULLANIM TALİMATI (STERIL OLMAYAN) (MEDLINE)
BREATHING SYSTEM INSTRUCTION FOR USE (NON STERILE) (MEDLINE)

IFU ID: AW.04.ALT.MED.0013-REV.02 (05/2024)

Contents

CS	čeština (Návod k použití)	2/36
DA	Dansk (Brugervejledning)	4/36
DE	Deutsch (Bedienungsanleitung)	6/36
EL	Ελληνικά (Οδηγιεσ Χρησησ)	8/36
EN	English (Instructions For Use)	10/36
ES	Español (Instrucciones De Uso)	12/36
FR	Français (Mode D'emploi)	14/36
HU	Magyar (Használati utasítás)	16/36
IT	İtaliano (Istruzioni Per L'uso)	18/36
NL	Nederlands (Gebruiksaanwijzing)	20/36
NO	Norsk (Bruksanvisning)	22/36
PL	Polski (Instrukcje Użytkowania)	24/36
PT	Português (Instruções de uso)	26/36
RO	Romana (manual de utilizare)	28/36
SK	Slovenčina (Návod na použitie)	30/36
SV	Svenska (Användarmanual)	32/36

Manufactured by:
Meditera Tıbbi Malz. San. ve Tic. A.Ş.
İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
35900 Tire İzmir
TÜRKİYE
Tel: (+90 232) 513 51 10 (Pbx) www.meditera.com.tr

Medline International BV.
Catharijnesingel 47, 3511 GC
Utrecht, Netherlands

CH importer: Medline International
Switzerland SARL, c/o MN & Associés
SA, 1 Place de Longemalle,
1204 Genève, Switzerland.

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Definition:
The breathing circuit is a device consisting of accessories such as connectors that enable the transportation of respiratory gases, air and/or anesthetic gases between the anesthesia device or ventilator and the patient. It is connected between the mechanically ventilated patient, who is connected to the breathing/anesthesia machine, and the ventilation/anesthesia device. These products are manufactured according to the requirements of ISO 5367 and ISO 5356 standards.

Accessory Description:
allow hoses to stay together for ease of use in the breathing circuit.

Nebulizer Kit: The nebulizer kit, connected to the breathing circuit with the help of a T connector, allows medication to be given to the patient through the breathing circuit.

Connectors: These are accessories that enable breathing circuits to be connected between the patient and the device in different shapes and positions depending on the intended use.

Cap: It is an accessory used to cover the exposed ends of the breathing circuits. Luerlock Syringe Without Needle: It is an accessory that is attached to the water trap section of the breathing circuit and provides a locked connection. It drains the fluid inside the respiratory circuit by sucking or pressing.

Limb: It provides a connection between the anesthesia balloon and the anesthesia device, which is used for manual ventilation of the patient in inhalation anesthesia and general anesthesia applications and for visual and tactile feedback when mechanical breathing is subsequently switched to.

Circuits with smoothbore tubes on the inside:
This type of hose contributes to mobility (flexibility) with its flexible structure, creates low resistance to air flow, is resistant to breakage and bending thanks to its external spirals, and with its flat internal structure, it does not accumulate water, body fluid, dirt and reduces the risk of bacterial growth. Water that may accumulate due to condensation is also ensured to flow easily into the water trap. Since they are transparent, their inner parts are easy to observe. They have low compliance. They are preferred for long-term use.

Circuits with corrugated tubes:
Since the inner surfaces are serrated, there may be more secretion, fluid and accumulation. They have a more matte structure. However, they are preferred for short-term use because they are light and useful.

Circuits with extendable tubes:
They can grow up to 3 times their own length. Thanks to this structure, they provide positioning optimization. In this way, the pulling effect on the airway can be reduced. They are especially preferred in extension (limb) lines. Within the framework of the above definitions, the breathing circuit is defined into subgroups as follows.

Dual Line Circuit Breathing Circuit:
It is the breathing circuit that provides the transportation of respiratory gases between the anesthesia device or ventilator and the patient.

Dual Lumen Breathing Circuit:
Thanks to a hose divided in half by a membrane, respiratory gases, air and/or anesthetic gases are supplied from a single line between the anesthesia device or ventilator and the patient. The dry and cold air passing through the inspiration line (going to the patient) is heated by the hot and humid air passing through the expiratory line (exhaled by the patient). The exhaled air heats the inspired air and thus natural heating is achieved. Since these circuits are single line and lighter, they reduce crowding on the patient, eliminate clutter in the environment, reduce weight-related shrinkage and separation, and increase mobility. Although they are frequently preferred as an anesthesia circuit, they can also be easily used for ventilation in intensive care units. There are types with water traps.

Expected Clinical Benefit:
Supporting the respiratory functions of patients with mechanical ventilation. Mechanical ventilation of patients and application of anesthetic gases in general anesthesia applications.

Intended Use:
The breathing circuit is a device consisting of components such as tubes and connectors which enable the transportation of respiratory gases, air and/or anesthetic gases between the anesthesia device or ventilator and the patient. It is connected between the mechanically ventilated patient, who is connected to the breathing/anesthesia machine, and the ventilation /anesthesia device. These disposable breathing circuits enable the respiratory gases of the mechanical ventilator or the respiratory and anesthesia gases of the anesthesia device to be delivered to the patient. Select the appropriate breathing circuit according to the patient. Higher resistance and compliance values may lead to inappropriate ventilation and injury to the patient.

Intended User of the Product and the Patient Group It Will Be Used : Neonatal, pediatric and adult patients who are mechanically ventilated or mechanically ventilated during general anesthesia and inhalation anesthesia applications.

Users:
They are people who use the medical device in accordance with its intended purpose. Users for this medical device are all healthcare professionals including doctors, nurses, anesthesia technicians, paramedics. Designated users have the necessary expertise, training and knowledge to use and install this medical device. This medical device should only be used or installed by identified users.

WARNING: This product should only be administered by qualified healthcare personnel.

Indications: It is indicated for use in mechanically ventilated neonatal, pediatric and adult patients, including general anesthesia and inhalation anesthesia. They are used between mechanical ventilators and the patient in mechanical ventilation. They are used between anesthesia machines and the patient for inhalation anesthesia and general anesthesia. Thanks to various ports on the connectors, values such as temperature, humidity, pressure, air flow can be measured and gas sampling can be performed. They make these transactions easier.

Contraindications: It has no known contraindications.

Usage: This device is used to deliver air, a combination of air and/or anesthetic gas and oxygen.

WARNING: If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

WARNING: Do not modify the breathing circuit parts or design.

WARNING: Install the breathing circuit in accordance with the instructions and intended use. The addition of other components or incompatible components may adversely affect the performance of the ventilator by increasing inspiratory and expiratory resistance.

WARNING: Do not extend the circuit, allow it to bend or remain under mechanical stress.

WARNING: Connections of all circuit parts must be checked before use.

WARNING: Before use on the patient, check the tightness and durability of the entire system and perform leakage and pressure drop tests. Set and review the relevant alarms of the anesthesia device.

WARNING: Blockages, damage and foreign objects may cause malfunction. Before installation, check all components as well as the breathing circuit for blockages.

WARNING: Once the patient is connected, observe the patient appropriately during the ventilation process. The patient, ventilator, anesthesia device and breathing circuit should be regularly monitored according to the standards specified by the hospital. Apply necessary monitoring.

WARNING: During use, frequently check that all connections are properly and securely connected and that there are no leaks or blockages.

WARNING: During use, make sure that the connection is correct for all alarms occurring on the device to which the circuit is connected, and check that the patient's condition and patient connections are appropriate, and that the circuit airway is open (kink, kink, obstruction, water accumulation, etc.).

Explanation of Symbols	
REF	StockCode
LOT	LotNumber
QTY	Quantity
	ExpireDate
	Manufacturer Information
	Single Use
	Read the IFU
	Do not use if package is damaged
	Keep away from sunlight
	Caution
MD	Medical Device
+40 +5 °C	Temperature Limits
70% 15%	Humidity Limits
NON STERILE	Non sterile.
	Keep dry.
LATEX	Do not contain Latex
PHT	Do not contain phthalate.
CE 0476	Notified Body Number
CH REP	Swiss Representative
UK REP	United Kingdom Representative
UDI	Unique Device Identifier

WARNING: Check frequently that there is no condensation in the circuit, as this may increase resistance to flow and cause the ventilator to alarm or cause errors in respiratory measurements.

WARNING: If the breathing circuit contains a filter;

a. Active humidification or application of different medications via a nebulizer may cause increased resistance in the filter;
b. Avoid nebulizing towards the bacteria/virus filter and using the filter as a water trap;
c. Do not connect the bacterial/virus filter directly to the active humidifier device output;
d. If a gas sampling line is added to the filter, make sure its cover is closed.

WARNING: If the anesthesia bag is placed separately in the circuit bag in the anesthesia circuit, connect the balloon tightly to the extension line and check the integrity of the connection before use.

WARNING: In breathing circuits containing a water trap with a fully enclosed, needle-free connector, the water trap body and cup are glued together. Please do not try to open it. You can empty this type of water trap with the help of a luer syringe that you will connect to the needle-free connector.

WARNING: In circuits containing water traps, the position of the water trap must be perpendicular to the ground plane.

WARNING: If the water trap is installed too high in the breathing circuit, condensed water may flow into the breathing hoses. Install the water trap at the lowest point of the breathing circuit and below the patient level.

WARNING: If the water traps are not emptied regularly, condensed water may flow into the breathing tubes. The breathing circuit may be partially or completely blocked. Check the water trap regularly for condensation and empty the breathing circuit if necessary.

WARNING: In circuits containing standard water traps, incorrect installation may cause incorrect operation of the water trap. Install the water trap correctly.

WARNING: The contents of the water trap may be contaminated. Dispose of the contents of the water trap in accordance with hospital regulations. Be careful when using the water trap.

WARNING: When removing the water trap bowl from the water trap, there is a higher risk of leakage. During leakage, PEEP may decrease and VT may not be achieved. After emptying, install and close the water trap container as quickly as possible.

WARNING: Using an improper suspension arm may damage the breathing circuit.

WARNING: Do not use products with damaged packaging.

WARNING: Avoid contact with chemicals, cleaning agents or hand sanitizer.

WARNING: Check the interior of the medical device regularly for accumulation of fluid or secretions. Liquids or secretions may cause increased resistance or enter the patient's airways.

WARNING: Apply this disposable product immediately after opening. The user must replace the product regularly in accordance with hygiene rules. Change the product in accordance with the hospital infection committee and clinical rules.

USAGE OF TEST ADAPTER:

1. Remove the elbow connector.
2. Connect the coaxial breathing circuit to the expiratory port of the device with the coaxial test adapter (A) (Altech Test Adapter).
3. To detect a leak in the internal tubing, perform a leak test in accordance with the device's operating instructions.
4. Do not use the breathing circuit in case of leakage.
5. Remove the coax test adapter and reconnect the elbow connector.
6. Attach the complete breathing circuit (hose, filter/HME, etc.) to the device.
7. Perform an appropriate leak test to detect leakage of the complete breathing circuit.

Product Usage Duration on the Patient :
The product should be replaced in accordance with hospital infection control and intensive care/anesthesia protocols. The product can be used in intensive care units for up to 7 days.

Disposal of the Product:
After use, the medical device must be disposed of in accordance with applicable hospital hygiene waste disposal regulations.

WARNING: Do not reuse this disposable product. Reuse may cause serious harm, including contamination with infectious substances or interruption of treatment.

WARNING: Do not immerse the product in water, wash, clean, disinfect or sterilize it. Reprocessing of disposable medical devices may result in decreased performance or loss of function.

Fire Risk:
It is not suitable for use with electrosurgery or laser when used with conductive and easily flammable substances such as oxygen and nitrous oxide, and flammable (non-AP/APG) anesthesia gases. If used in these situations, the device poses a risk of fire. To protect the patient and user, oxygen and nitrous oxide etc. Eliminate leaks in the carrying hoses. Before starting electrosurgery and laser surgery, purge the area around gas-carrying parts (endotracheal tube, mask, Y-piece, tubing, filter and breathing bag) with sufficient air (<25% O2); Also clean surgical drapes.
To avoid damage to the breathing circuit, ensure a distance of at least 200 mm (7.9 inches) between hoses carrying oxygen or nitrous oxide and potential sources of flame (e.g., electrosurgery and laser surgery devices).


General information about the product:
The product does not contain latex and phthalates.

Storage Conditions:
Store in a cool and dry place away from direct sunlight. Before use, check the parcels and packages for any damage. Store in accordance with the storage conditions stated on the label.

Shelf life:
If stored under appropriate storage conditions, the product can be used until the expiration date stated on the label.

Shipment of the Product:
The products are placed in a box and stretched; It is shipped on pallets or in bulk, according to the agreement made with the customer.

Label Information:
The product label has been prepared with reference to 93/42/EEC Medical Device Directive, ISO 15223 and ISO 15986 Standard. An explanatory table regarding the symbols on the label is given below.

Explanation of Symbols	
	StockCode
	LotNumber
	Quantity
	ExpireDate
	Manufacturer Information
	Single Use
	Read the IFU
	Do not use if package is damaged
	Keep away from sunlight
	Caution
	Medical Device
	Temperature Limits
	Humidity Limits
	Non sterile.
	Keep dry.
	Do not contain Latex
	Do not contain phthalate.
CE 0476	Notified Body Number
CH REP	Swiss Representative
UK REP	United Kingdom Representative
UDI	Unique Device Identifier

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New codes proposal	Product Name
212170M	CIRCUIT,SMOOTH,170CM,2LT,2WT,AT150CM,ADULT,NS
212240M	CIRCUIT,SMOOTH,240CM,2LT,2WT,AT150CM,ADULT,NS
111200M	CIRCUIT,SMOOTH,200CM,2LT,W/O WT,AT150CM,ADULT,NS
111300M	CIRCUIT,SMOOTH,300CM,1LT,2WT,AT300CM,PED,NS
213240M	CIRCUIT,SMOOTH,240CM,1LT,2WT,AT150CM,PED,NS
505200M	CIRCUIT,SMOOTH,200CM,1LT,W/O WT,AT150CM,PED,NS
191431M	T CONNECTOR, 3 X 22M/15F, NS
191423M	ELBOW CON.,PRESS. PORT,15M - 22M/15F,LL PORT,NS
41518M	CORRAFLEX CORRUGATED TUBING 1.80MT
T7F391M	CIRCUIT WITH FILTER, DUAL LIMB, 1.8M
1680M	CORRAFLEX CURRUGATED 30MT ROLL

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