Pulse Lavage System

Instruction Manual

REF 2500037 EN

General Description

This Pulse Lavage has an easy disposal Alkaline Battery Pack.
It is fully disposable single use, minimizing potential contamination, designed for use in orthopedic surgical procedures, trauma and general wound care.

Orthopaedic Bone Bed Preparation
Pulse lavage is an important equipment used to prepare bone bed before orthopaedic cementing during arthroplasty.
It has been proven to increase cement penetration into the cancellous bone, potentially improving fixation strength leading to longer survival of implants, thus reducing revision rates.

Wound Care
With the adequate steady pressurized, pulsed solution it can achieve a mechanical debridement for wounds, removing dirt, infectious microorganisms, wound exudates, necrotic tissue and debris. Thereby assisting visual examination and optimizing wound healing by stimulating granulation tissue formation.

SUCTION TUBE (OPTIONAL)

Standard package includes:
- 1 Hand piece
- 1 Alkaline battery pack (non rechargeable)
- 1 Small cone tip (ST 01) - Optional
- 1 Coaxial canal tip (CT 02) - Optional
- 1 Suction tube (3M) - Optional

Read the operating instruction carefully before use.
Failure to follow may cause damage or injuries.

Indications

- Arthroplasty bone cleansing
- Wound cleaning

Contraindications

- Should not be performed in the presence of active profuse bleeding (precautionary measures for patients on anticoagulation medication)

⚠ Warnings

Failure to follow may cause harm to person.

- The device must be carried out solely by Orthopedic surgeons and health professionals with wound care knowledge.
- Wear protective equipment, such as fluid-proof gown, gloves, mask or face shield, hair cover to avoid infectious splash backs.
- Use a drape or towel to cover all patient lines, ports, and wounds that aren't being treated.
- When applied on nerves, tendons and bones, should take precautions depending on tissue condition and the judgement of professional personnel.
- Follow standard practices to minimize potential contamination.
- Sterile only if package is unopened and undamaged.
- Single use only! DO NOT Re-sterilize!
- DO NOT operate damaged or malfunctioning instrument.

⚠ Caution

Failure to follow may cause damage to product.

- DO NOT use this product in any other manner other than the purpose or the operation method described in this instruction manual.
- DO NOT submerge instrument in liquid.
- DO NOT drop or hit the instrument.
- Store in a cool dry place.
- No modification of this equipment is allowed
- Do not mix with used or other types of batteries
- Dispose of all single use components, as appropriate, immediately after use.
- Ensure that the device is disposed of safely and in accordance with all current national and international waste disposal directives.
- Some national law restricts this device to sale by or on the order of a physician.

MD STERILE EO CE 0197

Instructions for Use


1. Peel open the package and remove all the items in a sterile way.
2. Attach the required tip to the handpiece and push on until the connection is secured (Fig.A).
 (When using the Intramedullary Tip ensure the **white cap** is removed before use)
3. Spike in the irrigation tip, which is on the battery pack, to the irrigation solution (Fig. B).



4. Attach the suction tube to the handpiece suction adaptor and hospital suction system.

8. Remove the protector from the Toggle Trigger.
 High speed: depress the upper half.
 Low speed: depress the lower half.
 Stop the flow: return to the original position.


NOTE:
- Both irrigation and suction tubes has clamps to prevent leakage of fluids.

In Case of Tip Obstruction: (Optional)

Insert the cleaning rod (HA 250) from the end of the tip and push until releasing the obstruction.

Disposal


- Make sure to remove the battery pack in case of disposal.
- Ensure that the device is disposed of safely and in accordance with all current national and international waste disposal directives.
- Please contact the local authorities responsible for waste disposal if questions arise.

Shelf Life

3 years

Operation

- Temperature: 10°C ~ 40°C
- Relative Humidity: 30%~75%RH
- Atmospheric Pressure: 80kPa ~106kPa
- Altitude limit: 2,000 meter

Transportation and Storage

- Temperature: 5°C ~ 40°C
- Relative Humidity: 15%~75% RH
- Atmospheric Pressure: 80kPa ~106kPa

Specifications

- Adjustable spray speed: High or Low
- Pressure: <15 PSI
- Pulse rate: 1500~1700 cycle/min
- High speed flow rate: 1000ml/min
- Low speed flow rate 700ml/min
- Tube Length 3 meters
- DC 12V / AA x 8 cells Alkaline Battery

(Data may vary according to type of power source, type of tip, distance and length of operation)

Labeling Symbols

	Refer to instruction manual/booklet
	General warning sign / Caution
	Do not use if package is damaged
	Do not re-use
	Keep dry
	Temperature limit
STERILE EO	Sterilized using ethylene oxide
CE 0197	Complies with EU directives
	Manufacturer
IP2X	Solid particle protection:> 12.5mm
LATEX	Not made with latex
	Premature unpacking warning
	Type B applied part
MD	Medical Device
	Direct current
REF	Catalog number
	Use by date
LOT	Batch code
	Date of manufacture
	European Union WEEE Directive Logo
	Do not resterilize
	Humidity limitation
ALKALINE	Recycle alkaline battery
EC REP	Authorized representative in the European Community
CH REP	Swiss Authorised Representative
UK RP	UK Responsible Person

ALKALINE

DO NOT dispose battery in fire
DO NOT recharge

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Medimap Ltd
miCUTec AG
Medline International B.V.
Medline Industries Ltd
Dongguan Kaiser Technology Co., Ltd.