FORMAT: 210x297mm

PMS COLORS

2935C | 2935C 60% | | 431C | | |

Perineal Instant Cold Pack

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Instruction for use

Danyang Rapid Aid Hot and Cold Therapy Products co., ltd,
Glasses Industrial Park, Danfu Road, Danyang, Jiangsu, China 212300

CE 1639

EC REP Lotus NL B.V. Koningin Julianaplein 10, le verd, 2595AA, The Hague, Netherlands

CH REP Casus Switzerland GmbH, Hinterbergstrasse 49, 6312 Steinhausen, Switzerland

CH importer: Medline International Switzerland Sàrl, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

REF CPACKP

11,4x36,2 cm (4.5x14.25 in)

RD23RAL 04-2023

en

Content: Urea blend and water

To activate:

1. Do not remove the pack from the bag prior to activating.
2. Fold the product in half top to bottom, lining up arrows, until the inner water bag ruptures.
3. Shake to thoroughly mix the chemicals and water.
4. Remove from the bag.
5. Remove paper strip to reveal adhesive and apply to undergarment.

Instant cold pack is for used for common aches, pains, swelling and bruising.

⚠ Warnings:

- Do not apply to open wound.
- For external use only.
- Harmful if swallowed.
- If accidently swallowed, drink several glasses of water, contact your local Poison Control Center immediately.
- Do not puncture, strike or step on bag.
- Should contents become exposed to skin or eyes, flush with water immediately and call a physician.
- Keep out of reach of children.
- Do not use on infants.
- Consult a physician before using on sensory-impaired individuals.
- Do not refrigerate or freeze to reuse.
- Check patient and pack regularly during product use.
- Do not apply cold therapy for more than 20 minutes.
- If a pack is leaking, discard immediately.
- Keep away from open flame.
- Protect from freezing, avoid excessive heat.
- Single use only.
- Discard after use.
- Reuse of medical devices intended for single use may result in degraded performance or a loss of functionality.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal Instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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