# INTERTEMP™
# Endocavity Probe Cover: Thermometer Cover
Instructions for Use

## en Endocavity Probe Cover: Thermometer Cover

**Instructions for Use**

**Intended Use / Indications for Use**
A device intended to be used as a physical barrier for protection against the effects of environmental exposure (e.g., body fluids, gels) and/or to maintain the required hygienic level of various diagnostic or surgical procedures which utilize a transducer (probe).
This is a single use device. [Not to be reprocessed.]
The devices are further indicated to be used in or on:
- Natural body orifices (e.g., vagina, rectum, esophagus)

**Intended Patient Population**
Probe covers have no limitation for any definition of the patient group and their estimated state/condition. Probe covers can be used at the discretion of the healthcare professional.

**Contraindications**
There are no known contraindications for these devices.

**Warnings**
- WARNING: CHOKING HAZARD. Don't turn patient (from stomach to back or reverse) during procedure while oral temp measurement is used.
- Do not use if sterile barrier package is damaged. If damage is found, call your Medline-Microtek representative.
- Visually inspect for breaches of packaging integrity (such as tears, cuts, open seals, etc.) prior to use.
- Report serious incidents to legal manufacturer/importer and competent authority.

**Precautions**
- Use care when handling the Probe Cover in close proximity to sharps, cautery devices, and other medical devices or procedures that may damage the integrity of the Probe Cover.

**Product / Device Description**
Probe Covers are composed of natural or synthetic materials and come in a variety of sizes and shapes to allow these devices to properly fit equipment with diverse shapes and uses.
Probe Covers are used in a variety of clinical and surgical settings by healthcare professionals. The devices cover a transducer / probe and, in some instances, the cables connected to the transducers / probes. These coverings offer a physical barrier between the transducer and the patient.

**How Supplied**
- Do not use if package is opened or damaged, store in a cool, dry, dark place.

**Reuse Warning Statement**
For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Handling and Storage**
- Keep away from sunlight
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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©2025 Medline Industries, LP. All rights reserved.
Microtek Medical B.V.