VAPREVENT® ORAL CARE

REF MDS096575EU

REF MDS096575EUB

en

Suction Toothbrushes Individually Wrapped
Intended Use:

Oral Care products are intended to clean oral secretions from the oral cavity, suction swabs and toothbrush also are intended to aspirate oral secretions, during oral care procedures. They can be used alone or in combination with Oral Care solutions to clean the oral cavity. These devices are not intended to be used by lay users.

1. With gloved hands, remove suction toothbrush from packaging. Ensure foam on suction toothbrush is intact.
2. Connect the suction toothbrush to standard suction lines and use the thumb port for suction control. Check for secure connection.
3. Saturate foam with water or desired solution.
4. Gently brush teeth and gums for approximately 1-2 minutes.
5. After use, ensure foam is intact before discarding toothbrush. Inspect oral cavity and remove any particles, if present.

Warning: Do not use for tracheostomy care. Consult physician prior to use to determine appropriateness for the patient. Discontinue use and consult physician if irritation or excessive bleeding occurs. Keep out of the reach of children. Protect from freezing. Avoid excessive heat.

Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established. Disposal Instructions: Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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©2026 Medline Industries, LP
Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.