NebEasy Nebulizer Mask Adult

with mask, nebulizer chamber, tubing
PVC, (DEHP-free)

Rx Only

GaleMed Corporation
87, Li-gong 2nd Road,
Wu-jia, I-Lan 268, Taiwan

Factory:
GaleMed (Xiamen) Co., Ltd.
Xiamen Area of China (Fujian) Pilot Free Trade Zone, 39,
Section 3, Haijing East Road 361026 Fujian Province, China
www.galemed.com Tel: 86-592-6880588
Made in China

CE 2460

EC REP
Obelis s.a.
Bd. Général Wahis 53, 1030
Brussels, Belgium

Distributed by:
Teleflex Medical Europe
IDA Business and Technology Park
Athlone, Co Westmeath, Ireland
www.teleflexmedical.com

P90-00040-02

en

Instructions for Use:
NOTE: These directions are general guidelines intended for use by qualified medical personnel.

1. Place the prescribed solution into the nebulizer chamber.
2. Attach the nebulizer cap to the nebulizer chamber. Attach the cap to the chamber with a clockwise rotation.
3. Fit the outlet of the nebulizer into the delivery device (tee and mouthpiece, face mask, breathing circuit, etc.). A flexible reservoir tube is provided with some models. Discard any items if not used.
4. Connect supply tubing between the nebulizer and the gas source.
5. Set the gas source to the proper flow.
 NebEasy : 8 LPM
6. Check for aerosol mist and gas through the device. If necessary, gently tap the nebulizer until it begins to nebulize.
7. During treatment, periodically tap the nebulizer to minimize residual volume.

WARNING: Be sure all connections are secure.

WARNING: Do not attach the outlet of the nebulizer directly to an endotracheal tube, tracheostomy tube or other device that does not offer an alternate pathway for exhalation.

WARNING:
This product is for single use only. Destroy after use. Do not re-sterilize.
This reduces the risk of cross- infection between patients.

Single Use Products
Reprocesssing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected, or sterilized.

Storage conditions:
Store in clean dry conditions away from heat and light.

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