
CENTURION®
Medical Products

Compass® Thoracentesis

English....2
German....6
French....10
Dutch....14
Italian....18
Spanish....22
Portuguese....26
Swedish....30
Czech....34
Slovak....38
Slovenian....42
Norwegian....46
Polish....50
Finnish....54
Danish....58
Greek....62
Romanian....66
Bulgarian....70
Estonian....74
Hungarian....78
Latvian....82
Lithuanian....86
Croatian....90
Russian....94
Turkish....98
Korean....102
Arabic....106
Hebrew....110
Malay....114
Serbian....118

English

CENTURION®
Medical Products

Compass® Thoracentesis

Intended Use
The Compass® Thoracentesis disposable pressure transducer with integrated digital display is intended for direct measurement of physiological pressure. Compass® devices are intended for use by a licensed healthcare professional within the general population including pediatrics.

Indications for Use
The Compass® Thoracentesis is indicated for use in the measurement of pleural pressure.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Information for Use: Please review prior to use. Sterilized using ethylene oxide. Device is sterile if pouch is unopened and undamaged. This product is designed for single use only. Do not resterilize or reuse. Not for use in oxygen-rich environments.

WARNING: No modification of this equipment is allowed. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

The Compass® device contains a Lithium Manganese button cell battery that cannot be removed from the device housing. The whole Compass® device should be disposed of as medical waste as facilitated by the healthcare facility's protocol and local member state and/or government requirements.

Description

The Compass® is a self-calibrating, single use, sterile device for measuring physiological pressure.

Suggested Use

1. Remove Compass® from its sterile packaging. Remove proximal luer cap if attaching accessories (e.g. stopcock or syringe) to proximal end of Compass®.
2. Activate and calibrate Compass® by depressing power button until display reads "00".

Warning: Turn on Compass® prior to connecting accessories to avoid incorrect baseline reading.

3. Connect Compass® to accessory monitoring equipment (e.g. needle, catheter, stopcocks, syringe, flush devices, extension sets, etc.) per hospital protocol.

NOTE: Attach distal (male luer) side of Compass® to needle or catheter in contact with patient.


Warning: Do not re-infuse blood or other fluid into patient after any such fluid is drawn into the Compass®

Warning: Pressure reading is NOT physiological if either the front or the rear of the device is open to atmospheric pressure.

4. Pressure is displayed on Compass® LCD

NOTE: If Compass® pressure reading does not return to zero after removal from needle or catheter, flush device with saline and then cycle device (turn it off and then back on) to recalibrate.


WARNING: Do not reuse Compass® for more than one procedure in the same patient. Use a new Compass® device for each procedure. Reuse may lead to patient infection or cross contamination due to lack of sterility.

1. Perform procedure per hospital protocol.

Device Power Off

Press button three times in a “PUSH”, “PUSH”, “HOLD” sequence; quickly push button twice and then hold third button push until device turns off.

Compass® Thoracentesis Specifications

Characteristics	Value
Pressure Range	-500 to +500 cm H_2O
Accuracy	±2.7 cm H_2O below 68 cm H_2O ±3% of reading above 68 cm H_2O
Zero Drift	±2 cm H_2O per 4 hours
Operating Life	8 hours
Operating Temperature	15 to 40°C (59 to 104°F)
Storage Temperature	-20 to 49°C (-4 to 120°F)
Light Sensitivity	< 2 cm H_2O
Lithium Battery Voltage	3 Volts

Software Error Codes

Code	Meaning
SF	Software Failure
HF	Hardware Failure
PF	Pressure Sensor Failure
bF	Battery Failure
EE	Environmental Error-Outside operating temperature range
CE	Calibration Error-Outside required accuracy

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©2024 Medline Centurion
Medline International B.V.
BSI Assurance UK Ltd,
Medline Industries Ltd,