VAPREVENT® ORAL CARE 24-HOUR KIT

REF Q400 REF Q800

en Oral Care, Q4/Q8 24-Hour Kit, Instructions for use

Intended Use:

Oral Care products are intended to clean oral secretions from the oral cavity, suction swabs and toothbrush also are intended to aspirate oral secretions, during oral care procedures. These devices are not intended to be used by lay users or with freestanding suction vacuum equipment. They can be used alone or in combination with Oral Care solutions to clean the oral cavity.

1. Open the box as indicated and dispense a tray.
2. With gloved hands, peel back film. Ensure foam on Suction Toothbrush and/or suction swab and DenTips® is intact. Connect Suction Toothbrush/Swab to standard suction lines. Ensure connection is secured.
3. Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
4. Saturate toothbrush/swab with water or desired solution.
5. Gently clean the teeth and gums for approximately 1-2 minutes. Provide suction by using the thumb port.
6. Ensure foam is intact before discarding suction toothbrush/swab. Inspect oral cavity and remove any particles, if present.
7. Open Medline mouth moisturizer and squeeze onto DenTips®. Liberally coat inside of mouth (teeth, lips, tongue and palate).
8. Ensure foam is intact before discarding DenTips®. Inspect oral cavity and remove any particles, if present.

Warning:

Consult physician prior to use to determine appropriateness for the patient. Discontinue use and consult physician if irritation or excessive bleeding occurs. Do not use Dentips® or suction swab to clean a tracheostomy. Keep out of the reach of children.
Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Warning for medline mouth moisturizer: Not for use on children of 12 years old and under, unless directed by a healthcare professional. Must not be swallowed.

Disposal Instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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Packaged in Mexico for:
Medline Industries, LP,
Three Lakes Drive, Northfield, IL 60093 USA.

CE 2797

EC REP Medline International France SAS,
5 rue Charles Lindbergh, 44110
Châteaubriant, France

CH REP **CH Rep and importer:** Medline International
Switzerland SARL, c/o MN & Associés SA, 1 Place de
Longemalle, 1204 Genève, Switzerland

UK CA 0086 **UK representative and importer in UK:**
Medline Industries Ltd, WA4 6HL,UK

Medline International B.V.,
Catharijnesingel 47, 3511 GC
Utrecht, The Netherlands

Medline, Vaprevent and Dentips are registered trademarks of Medline Industries, LP.
©2025 Medline Industries, LP
Pending and U.S. Patent 8, 875, 940

RF25TTC - 06/2025