OR Towel (with x-ray thread)

Instructions for use

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OR Towel (with x-ray thread)

Product composition:

The product is mainly composed of bleached white or dyed gauze and x-ray thread.

Intended use:

It's a kind of invasive dressing, which penetrates inside the body used for absorption of exudates or protecting organs/tissues in surgical operation. X-ray detectable element helps to detect the device left inside patient. This is a professional and single use device. Normally intended for continuous use for not more than 24 hours.

Intended user:

Medical staff

Intended patient population:

Used for all age groups, not recommended for infants to use.

Handling instruction:

- Medical staff should clean and disinfect hands to ensure aseptic operation.
- Check packaging and product are intact before use.
- Open the package take out the product and make it cover the surgical area according to the size and shape of the surgical field.

Warnings and precautions:

- The product is disposable and cannot be reused.
- Stop using if product has any defect.
- Stop using if the product has expired.

Shelf life:

5 years

Contraindications:

Do not use on patients with a known hypersensitivity to the product itself or to its components.

Side effects:

None

Disposal:

Disposal should be handled according to local environment procedures.

Special storage conditions and methods:

The product should be stored in a cool and dry place and avoid direct sunlight. The transportation should prevent the weight, sun and rain.

Notice:

Any serious adverse event related to this product must be report to Allmed Medical Products Co., Ltd or authorized representative or national health authority in your country.

Explanation of figures, symbols and abbreviations used in labels:

	Do not reuse		EO sterilization
	Used-by date		CE symbol
	Batch Code		Catalogue Number
	Single sterile barrier system		Keep away from sunlight
	Do not use if package is damaged		Importer
	Warnings and precautions		Manufacturer name address
	Manufacture date		Consult instructions for use
	European Authorized Representative		Medical device
	Keep dry		Do not resterilize

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CE 0123 EU REP Shanghai International Holding Corp. GmbH (Europe), Eiffestrasse 80, 20537 Hamburg, Germany
CH REP MedNet Swiss GmbH
Medline International B.V.,