HUDSON RCI®

Tracheostomy Masks

REF 1075, 41076

NON STERILE

PHT DEHP

Rx ONLY

Do not reuse. Non riutilizzare. Nicht wiederverwenden. Ne pas réutiliser. Får inte återanvändas. No reutilizar. Não reutilizar. Niet hergebruiken. Ei saa käyttää uudelleen. Må ikke genbruges. Skal ikke gjenbrukes. Μην επαναχρησιμοποιείτε. Nie używać ponownie. Tekrar kullanmayın. Повторно не использовать. Ne uporabite ponovno

CE 1639

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Not made with natural rubber latex.

Directions for Use:

NOTE: These directions are general guidelines intended for use by qualified medical personnel.

Connect aerosol tubing (not supplied) between the mask and the gas source. The mask inlet swivels 360° to position the tubing for supine or upright patients. Set the gas to the proper flow and check for gas flow through the device. Position the elastic strap behind the neck. Gently pull the end of the strap until the mask is secure. When using suction, loosen the mask and drop the mask out of the work area.

REF	Adult Size	Paediatric Size
1075	X	
41076		X

Warning: Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

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© Copyright 2009, 2011, 2013, 2016, 2018 Teleflex
Hudson RCI s a registered trademark of Teleflex Incorporated or its affiliates.