Electrosurgical Cautery Tip Cleaner

de Elektrodenspitzen-Reiniger
fr Grattoir pour pointe de bistouri électrique
nl Elektrochirurgische reinigingspads cauterisatiepunten
it Spugnetta abrasiva pulizia punta elettrobisturi
es Limpiador de punta de electrobisturí
pt Limpador de Pontas de Cautério Eletrocirúrgico
sk Čistič hrotu elektrochirurgického kautera
cs Čistič hrotu Elektrochirurgického kauteru
sv Elektrokirurgisk rengörare för kauterispetsar
fi Sähkökirurgisen kauterikärjen puhdistin
da Rensepads til elektrokirurgiske spidser
no Rengjøringspute for elektrokirurgisk kauterispiss
el Καθαριστικό Άκρων Ηλεκτροχειρουργικού Καυτηριασμού
pl Czyścik do końcówek elektrochirurgicznych
sl Čistilo konice elektrokirurškega kavterja
ro Dispozitiv de curățare pentru vârf electrochirurgical
bg Почистващ препарат за накрайници на електрохирургична каутеризация
hu Elektrokirurgilise kauteriseerimisotsaku puhastaja
lv Elektrokauter hegy tisztító
lt Elektrokautera uzgaļu tīrītāji
hr Elektrochirurginė kauterizacijos antgalių valymo priemonė
et Sredstvo za čišćenje vrha elektrokirurškog kautera
tr Elektrocerrahi Koter Uç Temizleyici
ar منظف طرف الكي الجراحي الكهربائي

REF TIPCLNR

en

Instruction for Use Cautery Tip Cleaners

Note: Please carefully read instruction for use before usage of product

Product name:

Cautery Tip Cleaners

Product Model and specification:

Table 1-1 Model List of Cautery Tip Cleaner

No.	Specification	Reference code
1.	4.5×4.4cm	TIPCLNR

Sterilization:

The product is sterilized by ethylene oxide, for single use only.

Intended use:

As an accessory of electrosurgical products, it can quickly and effectively clean the charred tissue and other adhesions attached to the head of the Electrosurgical Pencils.

Indications:

1. Check that the packaging is not compromised or the content is not expired.
2. Open packaging using aseptic technique.
3. Inspect the device for overall condition and physical integrity.
4. The device can be used to clean charred tissue and other adhesive materials quickly and effectively in electro-coagulation and electro-cutting operation.

Warnings:

1. The device is single use only. Reuse of single devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function, which may lead to injury, illness and/or death of a patient.
2. Do not use if the packaging is damaged.
3. Do not re-sterilize and/or reprocess.
4. If a serious incident has occurred in relation to the device, it should be reported to the supplier and the competent authority of the Member State in which the user and/or patient is established.

Product performance:

1. Working conditions

Ambient temperature:10—40°C

Relative humidity:30%~75%

2. Appearance

The appearance is Cautery Tip Cleaners, the cut edge is neat and smooth, there is no sawtooth dander, no impurity black spots on the back of the adhesive surface, and the adhesive surface of the sponge is tightly bonded, without splitting or sticking.

3. Safety Requirements

a. The cleaning sheet must be used by a professional nurse or doctor.

b. According to clinical requirements, select the appropriate model of the cleaning sheet.

c. Open the box, take out cleaning sheet from packaging, attached to the proper position, then start the operating.

Storage and handling:

The device must be stored in a clean and dark place with no corrosive gases and a good ventilation.

Ambient temperature range:(-20~50)°C

Relative humidity range: ≤80%

Atmospheric pressure range:(86.0~106.0)kPa

Interpretation of medical device symbols:

No.	Symbol	Title of symbol
1		Manufacturer
2		Authorized representative in the European Community
3		Do not use if package is damaged and consult instructions for use
4		Do not re-use
5		Caution
6		Sterilized using ethylene oxide

7		Do not resterilize
8		FRAGILE
9		Keep dry
10		Temperature limit
11		Medical device
12		This product has passed CE certification

Disposal instructions:

Discard after use. Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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Zhejiang Shuyou Surgical Instrument Co., Ltd.
EC REP Luxus Lebenswelt GmbH, Kochstr. 1, 47877, Willich, Germany
CH REP CMC Medical Devices GmbH, Rigistrasse 3, 6300 Zug, Switzerland
UK representative: Kingsmead Service Limited, 19 Mezzanine Floor, 19-21 Crawford Street, London, W1H 1PJ
Importer in UK: Medline Industries Ltd, WA4 6HL, UK
Medline International B.V.,