INSTRUCTIONS FOR USE

Cautery Tip Cleaner

aspen surgical.

Sterile

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone +1 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EC REP
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

CE
2797

For Country of Origin See Product Label

English

Intended Use

Intended to assist in cleaning and removing material from cautery tips during surgical procedures.

Intended User/ Patient Target Groups

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instructions for Use:

- Remove backing and apply to appropriate surface.
- Wipe cautery tip on pad to remove debris.
- Device is intended for use less than or equal to 16 hours.
- Do not use product if product is damaged.

Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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AB'deki Kullanıcılara ve/veya Hastalara Yönelik Bildirim
Emergo Consulting (UK) Limited