REF SENPAD01
REF SENPAD02

en SENSOR PAD - SINGLE PATIENT USE

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CE UK CA

Manufactured by:
Rondish Company Limited, Unit G&H, 4/F, Block 1, Kwai Tak Ind Ctr 15 - 33 Kwai Tak St Kwai Chung, N.T., Hong Kong

Imported and distributed in Europe by:
Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

RH21ROND 2021-08

en

SET-UP INSTRUCTIONS:

1. Battery cover needs to be opened with a screwdriver. Put 3 x AA alkaline batteries into the monitor and replace cover.
2. Loop cord together and thread around strain relief hook near the plug socket. Plug cord into bottom of the monitor.
3. Position user/resident in the bed or chair, so the caregiver will be alerted if they get up unsupervised, or fall out.
4. Test the pad by pressing before positioning the user/resident on the bed or chair. Monitor should beep twice. Active LED should now be blinking. Remove pressure and alarm should beep continuously. Press Reset to silence.

CARE INSTRUCTIONS:

1. Prevent cable pulling and breaking using the strain relief feature at the bottom of the monitor.
2. Monitor will start chirping and Battery Low LED indicates when batteries need to be replaced.
3. Pads should not be immersed in any solutions, especially heated. Cleaning should be using anti-bacterial cleaners and disinfectant “wet” wipes.
4. Replace pad with each new user/resident.
5. Do not fold pads.
6. When using a sensor with a soft seat, always test to make sure it works.

TROUBLESHOOTING:

1. Read the Monitor Instructions to familiarize yourself with the model.
2. Make sure the pad is plugged into the monitor.
3. Make sure your batteries are good, or, if using and AC adapter, ensure it's plugged in.
4. Check the monitor switches.
5. If the above are OK, call your distributor with the monitor and pad set up, to resolve.

STORAGE INSTRUCTIONS:

- Store in a low-humidity warehouse (32-85%) between 10-25C
- Storage period: 2 years
- This device must be disposed of properly according to local ordinances

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