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
## en GENERAL PURPOSE PROBE COVER KIT FEATURING ISOSILK™

### Instructions for Use

**Intended Use / Indications for Use**
A device intended to be used as a physical barrier for protection against the effects of environmental exposure (e.g., body fluids, gels) and/or tomaintain the required hygienic level of various diagnostic or surgical procedures which utilize a transducer (probe).
This is a single use device. [Not to be reprocessed.]
The devices are further indicated to be used in or on:

- Intact or compromised skin
- Natural body orifices (e.g., vagina, rectum, esophagus)
- Intra-operatively (e.g., on internal organs)

**Intended Patient Population**
Probe covers have no limitation for any definition of the patient group and their estimated state/condition. Probe covers can be used at the discretion of the healthcare professional.

**Contraindications**
There are no known contraindications for these devices.

**Warnings**

- Do not use if sterile barrier package is damaged. If damage is found, call your Ecolab-Microtek representative.
- Report serious incidents to legal manufacturer/importer and competent authority.

**Precautions**

- Secure distal tip to probe with rubber bands or ties to prevent gel from migrating away from tip.
- Use care when handling the Probe Cover in close proximity to sharps, cautery devices, and other medical devices or procedures that may damage the integrity of the Probe Cover.

**Product / Device Description**
Probe Covers are composed of natural or synthetic materials and come in a variety of sizes and shapes to allow these devices to properly fit equipment with diverse shapes and uses.
Probe Covers are used in a variety of clinical and surgical settings by healthcare professionals, such as;

- **General Purpose** – For use in biopsy, trauma, emergency room procedures, patients in isolation, catheter and drainage sites
- **Intraoperative** – For use in biopsy, abdominal surgery, cardiac surgery, vascular surgery, and special procedures
- **Endocavity** – For use in transvaginal, transrectal and prostate imaging

The devices cover a transducer / probe and, in some instances, the cables connected to the transducers / probes. These coverings offer a physical barrier between the transducer and the patient.

**How Supplied**

- Do not use if package is opened or damaged, store in a cool, dry, dark place.
- This product may be supplied with ultrasound gel, rubber bands and/or ties.

**Reuse Warning Statement**
For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Handling and Storage**

- Keep away from sunlight
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

**Warranty Statement**
The manufacturer warrants that this product will meet the manufacturer's standard written specifications then in effect and will be fit for the uses and purposes described in manufacturer's literature. As the sole remedy for breach of this warranty, Manufacturer will, in its reasonable discretion, either replace any product that did not meet this warranty or refund the purchase price paid for those products. Handling, storage, cleaning and sterilization of this product as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond manufacturer's control directly affect the product and the results obtained from its use. This warranty does not apply to (i) damage resulting from misuse, neglect, accident or improper use of any product by any person or entity other than manufacturer or (ii) any product altered by any person or entity other than manufacturer. **ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE DISCLAIMED. FURTHER, MANUFACTURER ASSUMES NO LIABILITY WITH RESPECT TO PRODUCTS REUSED, REPROCESSED OR RESTERILIZED AND MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WITH RESPECT TO SUCH PRODUCTS. MANUFACTURER ALSO DISCLAIMS LIABILITY FOR ALL CONSEQUENTIAL, INCIDENTAL, AND/OR SPECIAL DAMAGES CAUSED BY OR RELATED TO THE PRODUCTS.**

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