INSTRUCTIONS FOR USE

Sterion®
Suture Aid Booties

Aspen Surgical.

STERILE EO | MD | R_xONLY

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EC REP
Hillrom SAS
Zone Industrielle de Talhouet
56330 Pluvigner, France

Made in USA

CE 2797

English

Intended use: For one time use to protect, tag and help locate sutures.

SUTURE AID BOOTIE

1. After removing backing from cartridge, apply to any appropriate surface.
2. Insert tip of appropriate instrument into Suture Aid Bootie and remove from cartridge.
3. Dispose or recycle general surgical devices per facility protocol.

EXTRA SUTURE AID BOOTIE

1. Apply Extra Aid Booties over tip of appropriate instrument.
2. Dispose or recycle general surgical devices per facility protocol.

Precaution: Device is single use only. Reuse of device could result in infection/contamination and/ or device faliure which could lead to patient harm. Device is intended for limited contact duration, less than or equal to 24 hours.

Intended User/ Patient Target Groups: Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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