Aerosol Mask

Instructions for use

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DEHP free, DEHP-frei, Sans DEHP, DEHP vrij, Senza DEHP, Sin DEHP, Sem DEHP, Bez DEHP, Neobsahuje DEHP, DEHP-fri, DEHP-vapaa, DÉHP-fri, DEHP-fri, Χωρίς DEHP, Nie zawiera DEHP, Brez DEHP, Fără DEHP, не съдържа DEHP, DEHP mentes, DEHP içermez.

Ningbo Shengyurui Medical Appliances Co., Ltd.
No. 138, Binhaisi Road, Hangzhou Bay New Zone. 315336, Ningbo, Zhejiang Province, P.R CHINA

CE 1639

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

CH REP MedNet Swiss GmbH
D4 Platz 4, 6039 Root D4, Switzerland

CH importer: Medline International Switzerland Sàrl,
c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK representative: Etheria Medical Ltd, The Old brush Factory Unit 2D Whickham Industrial Estate Swalwell Newcastle Upon Tyne NE16 3DA United Kingdom
Importer in UK: Medline Industries Ltd, WA4 6HL, UK

Medline International B.V.,
Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Distributed in Australia by: Medline International
Two Australia Pty. Ltd., Sydney, Australia

RL24SWY 2024-11

REF 1083S

REF 41085S

en Aerosol Mask

Intended use:
Aerosol Masks are intended to deliver aerosol therapy to the patient for inhalation of nebulized medications or for humidification of the respiratory tract.

These devices are intended for use by medical personnel, lay persons, or patients that have been prescribed aerosol therapy.

1083S: Intended for use on adults only.

41085S: Intended for use on pediatrics only.

Intended to be used for a maximum of 30 days.

Instructions for use:
1. Connect aerosol tubing or a nebulizer to the mask inlet.
2. Check for secure connections by confirming flow of aerosol through device.
3. Place mask on face with the elastic strap below the ears and around the neck.
4. Gently pull the ends of the strap until the mask is secure.
5. Mold the metal strip on the mask to fit across the bridge of the nose.

Disposal Instructions:
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warnings:
Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Non-Sterile. Do not sterilize prior to use.

Do not use if device is soiled or damaged.

If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Caution should be exercised when delivering oxygen via an aerosol mask to patients with chronic obstructive pulmonary disease or facial trauma.

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