Fig.1 Fig.2 Fig.3

GB

SwabFlush®
Manufactured by: Excelsior Medical
0.9% Sodium Chloride Injection, USP with Integrated SwabCap®
ZR® Flush Syringe

Description: SwabFlush® is a saline I.V. flush syringe with an integrated SwabCap® A SwabCap® is a luer access valve cap that contains 70% isopropyl alcohol as a disinfectant. They are integrated together for the convenience of the user.

Available Sizes: 3 ml, fill in 10 ml, syringe, 5 ml, fill in 10ml, syringe, 10 ml, fill in 10ml, syringe

Flush Syringe

Sterile fluid and fluid pathway. Each syringe, unopened and undamaged, contains a sterile, non-pyrogenic isotonic solution of 0.9% Sodium Chloride Injection USP, 9mg NaCl per ml, with an osmolarity of 0.31mOsm/mL, pH 4.5 – 7.0. The syringe is a luer lock design that is luer lock fitting which is intended to connect to standard female luer lock fittings, including luer access valves.

Pharmacology: 0.9% Sodium Chloride Injection USP is a sterile aqueous injection having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues.

Indications for use: 0.9% Sodium Chloride Injection, USP prefilled syringes are intended for flushing of IV catheters and IV tubing only.

Contraindications: Hypernatremia and fluid retention, when the administration of sodium can produce edema and/or congestive heart failure.

Precautions: If compatibility is in doubt, consult the appropriate specialized literature. **Do not use if solution is discolored or contains a precipitate.** The unused SwabCap® in the syringe should be reserved for use to cover the accessed valve once it has been flushed using the SwabFlush® syringe.Use by expiration date on package. Prior to use, check with the drug, IV tubing and/or catheters' manufacturer instructions for use to ensure compatibility with 0.9% sodium chloride solution. If 0.9% sodium chloride is not compatible, follow the drug, IV tubing and/or catheters' manufacturer instructions for flushing practices.

Adverse reactions: None known in these quantities.

Overdosage: None known in these quantities.

Dosage and Administration:: This product should be **inspected visually for particulate matter and discoloration prior to use.** Maintenance of patency of intravenous access devices: Solution is injected as a single dose via the injection site of the vascular access device. The volume injected is equivalent to that of the vascular access device. After each use of the vascular access device for injection or infusion of medication, the solution should be injected to restore the patency of the access device. Aspirate before administering any solution via the access device lock in order to confirm patency and location of the catheter tip.

SwabCap®

Each syringe contains a sterile SwabCap® containing 70% isopropyl alcohol.

Indications for use: SwabCap® is intended for use on swab-able luer access valves as a cover to protect the luer access valves from potential contamination. The SwabCap® acts as a physical barrier to contamination between the luer access valve once it has been flushed using the SwabFlush® syringe.Use by expiration date on package. Prior to use, check with the drug, IV tubing and/or catheters' manufacturer instructions for flushing practices.

Adverse reactions: None known in these quantities.

Overdosage: None known in these quantities.

Dosage and Administration:: This product should be **inspected visually for particulate matter and discoloration prior to use.** Maintenance of patency of intravenous access devices: Solution is injected as a single dose via the injection site of the vascular access device. The volume injected is equivalent to that of the vascular access device. After each use of the vascular access device for injection or infusion of medication, the solution should be injected to restore the patency of the access device. Aspirate before administering any solution via the access device lock in order to confirm patency and location of the catheter tip.

SwabCap®

Each SwabCap® contains a sterile, single use, unopened and undamaged, sterile, non-pyrogenic isotonic solution of 0.9% Sodium Chloride Injection USP, 9mg NaCl per ml, with an osmolarity of 0.31mOsm/mL, pH 4.5 – 7.0. The syringe luer lock fitting which is intended to connect to standard female luer lock fittings, including luer access valves.

Pharmacology: 0.9% Sodium Chloride Injection USP is a sterile aqueous injection having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues.

Indications for use: 0.9% Sodium Chloride Injection, USP prefilled syringes are intended for flushing of IV catheters and IV tubing only.

Contraindications: Hypernatremia and fluid retention, when the administration of sodium can produce edema and/or congestive heart failure.

Precautions: If compatibility is in doubt, consult the appropriate specialized literature. **Do not use if solution is discolored or contains a precipitate.** The unused SwabCap® in the syringe should be reserved for use to cover the accessed

Fig.1 Fig.2 Fig.3

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