Sterilization Instructions for Kit Packers

DO NOT COPY STERILIZATION INSTRUCTIONS FOR END USERS

en

DO NOT COPY STERILIZATION INSTRUCTIONS FOR END USERS
Sterilization Instructions for Kit Packers
This drape is offered NON-STERILE, in bulk packaging intended for sale to kit packers for inclusion in surgical packs / kits which must be sterilized prior to distribution as a finished product. The kit packer is responsible to validate and verify their chosen packaging and sterilization conditions to ensure the finished product is fully sterilized and meets the functional performance and stability requirements of all local, national and international standards and regulations. The kit packer is responsible to label the finished product with the sterilization method used and the Use By Date provided by Medline-Microtek on the polybag label. Medline-Microtek has packaged items in Poly/tyvek peel pouches within a single wall corrugated cardboard shipping container prior to Ethylene Oxide (EO) sterilization under the following conditions: Gas Dwell Temperature: 125 - 140°F, Relative Humidity: 55-85%, EO Gas Dwell: 4 to 4.75 Hours, EO Gas Concentration 400 – 600 mg/L. When sterilized in this way, items can be EO sterilized up to two (2) times. The kit packer should: 1) inform the end user not to re-sterilize or reuse the drapes 2) ensure the final kit expiration date does not exceed the date on this product carton label, and 3) provide a finished sterile product and the instructions for use (above) with appropriate translations as needed. The end-user instructions for these products are provided in the section titled "Instructions for Use." It is the responsibility of the kit packer to provide end-user instructions and translations as appropriate.

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©2025 Medline Industries, LP. All rights reserved.
Microtek Medical LLC
Northfield, IL 60093 U.S.A.
Microtek Medical B.V.
Medline Industries Ltd